

Method 508Date Extracted 10/27/01

Page _____

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count
	LB	300ml		1ml			2ml		7	
	6FCH				1ml					
	MDL CH				1ml					
	PCL 1232				1ml					
	25599									6
	25600									7
	25601									8
	25602									9
	25603									10
	25604									11
Remarks										
Reviewed by, <u>CB</u> 1/13/02										
drip rate must be between 5-15 ml/min for liq-liq extractors										
Extraction type - <u>Son Fun</u> Solv. lot# <u>V29E68</u> NaCl lot# <u>V11699</u> Na2SO4 lot# <u>V32603</u>										
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	<u>Chloroform 1993</u>									
MDL/MDRF										
O. CHECK										
INT. STD.										
SURROGATE	<u>2064L</u>									
LPC										
ENDRIN/DDT										

Analyst C.H.
revision 4/12/01

Verified by _____

Reviewed by DW
Date 12/3/11