

Method 508 mod.

Date Extracted 1/29/02

Page 52

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
LB		500mL		1mL			1mL	2B	7			
REF					5mL			2A				
1960								3B				
1961								3B				
1962								3A				
1963								1A				
1964								4A				
1965								1A				
1966								1B				
1967								2A				
1968								4A				
1969								5A				
1970								1B				
Remarks: many failures, will be re-ex. entire batch												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type -		Solv. lot# <u>V41470</u>		NaCl lot# <u>V1599</u>		Na2SO4 lot# <u>V39H00</u>						
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE	<u>91579</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.	<u>022202</u>											
SURROGATE	<u>GCMS 011602</u>											
LPC												
ENDRIN/DDT												

Analyst JA

Reviewed by DW
Date 3/22/02

Sequence Name: C:\HPCHEM\1\SEQUENCE\FEB2502.S

Comment:

Operator:

Data Path: C:\HPCHEM\1\DATA\feb2502\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

(X) Full Method

(X) Inject Anyway

() Reprocessing Only

() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
------	------	------	----------	--------	-------------

1	Sample	1	DFTPP	GCMS0208	
2	Sample	2	CAL25	GCMS0208	
3	Sample	3	LB	GCMS0208	gc/ms scan 1/25
4	Sample	4	REF	GCMS0208	gc/ms scan 1/25
5	Sample	5	1643	GCMS0208	gc/ms scan 1/25
6	Sample	6	1650	GCMS0208	gc/ms scan 1/25
7	Sample	7	1651	GCMS0208	gc/ms scan 1/25
8	Sample	8	1687	GCMS0208	gc/ms scan 1/25
9	Sample	9	1688	GCMS0208	gc/ms scan 1/25
10	Sample	10	1689	GCMS0208	gc/ms scan 1/25
11	Sample	11	1690	GCMS0208	gc/ms scan 1/25
12	Sample	12	1691	GCMS0208	gc/ms scan 1/25
13	Sample	13	1692	GCMS0208	gc/ms scan 1/25
14	Sample	14	1693	GCMS0208	gc/ms scan 1/25
15	Sample	15	1694	GCMS0208	gc/ms scan 1/25
16	Sample	16	1732	GCMS0208	gc/ms scan 1/25
17	Sample	17	LBA	GCMS0208	gc/ms scan 1/25
18	Sample	18	REFA	GCMS0208	gc/ms scan 1/25
19	Sample	19	1733	GCMS0208	gc/ms scan 1/25
20	Sample	20	1734	GCMS0208	gc/ms scan 1/25
21	Sample	21	1735	GCMS0208	gc/ms scan 1/25 fb
22	Sample	22	1796	GCMS0208	gc/ms scan 1/25
23	Sample	23	1797	GCMS0208	gc/ms scan 1/25
24	Sample	24	1798	GCMS0208	gc/ms scan 1/25
25	Sample	25	1799	GCMS0208	gc/ms scan 1/25
26	Sample	1	DFTPP1	GCMS0208	
27	Sample	2	CAL25A	GCMS0208	
28	Sample	26	1800	GCMS0208	gc/ms scan 1/25
29	Sample	27	1801	GCMS0208	gc/ms scan 1/25
30	Sample	28	1802	GCMS0208	gc/ms scan 1/25
31	Sample	29	1803	GCMS0208	gc/ms scan 1/25 fb
32	Sample	30	LBB	GCMS0208	gc/ms scan 1/29
33	Sample	31	REFB	GCMS0208	gc/ms scan 1/29
34	Sample	32	1960	GCMS0208	gc/ms scan 1/29
35	Sample	33	1961	GCMS0208	gc/ms scan 1/29
36	Sample	34	1962	GCMS0208	gc/ms scan 1/29
37	Sample	35	1963	GCMS0208	gc/ms scan 1/29
38	Sample	36	1964	GCMS0208	gc/ms scan 1/29
39	Sample	37	1965	GCMS0208	gc/ms scan 1/29
40	Sample	38	1966	GCMS0208	gc/ms scan 1/29
41	Sample	39	1967	GCMS0208	gc/ms scan 1/29
42	Sample	40	1968	GCMS0208	gc/ms scan 1/29

Last Modified: Tue Feb 26 12:24:47 2002

Page: 1

43 Sample 41 1969 GCMS0208 gc/ms scan 1/29 |

Sequence Name: C:\HPCHEM\1\SEQUENCE\FEB2502.S

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	42 1970	GCMS0208	gc/ms	scan 1/29