

Method 508 mod

Date Extracted 12/24/01

Page 87

ref

GCMS-16612

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count																				
LB		500ml		1ml			1ml		7																					
Ref					0.5ml																									
31341	see below		1/28/02																											
31250																														
31277																														
31299																														
31276																														
31395																														
31342																														
31275																														
31342																														
31346																														
31340																														
Remarks Reviewed by, <u>LB 5/15/02</u> drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type <u>sp fun</u> Solv. lot# <u>V39406</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V39100</u> <table border="1"> <tr> <td>STANDARD</td> <td>INT STD LOG#</td> </tr> <tr> <td>CALIBRATION</td> <td><u>HPD/3092</u></td> </tr> <tr> <td>REF/SPIKE</td> <td><u>91/592</u></td> </tr> <tr> <td>MDL/MDRF</td> <td></td> </tr> <tr> <td>Q. CHECK</td> <td></td> </tr> <tr> <td>INT. STD.</td> <td><u>150/3002</u></td> </tr> <tr> <td>SURROGATE</td> <td><u>12/101</u></td> </tr> <tr> <td>LPC</td> <td></td> </tr> <tr> <td>ENDRIN/DDT</td> <td></td> </tr> </table>													STANDARD	INT STD LOG#	CALIBRATION	<u>HPD/3092</u>	REF/SPIKE	<u>91/592</u>	MDL/MDRF		Q. CHECK		INT. STD.	<u>150/3002</u>	SURROGATE	<u>12/101</u>	LPC		ENDRIN/DDT	
STANDARD	INT STD LOG#																													
CALIBRATION	<u>HPD/3092</u>																													
REF/SPIKE	<u>91/592</u>																													
MDL/MDRF																														
Q. CHECK																														
INT. STD.	<u>150/3002</u>																													
SURROGATE	<u>12/101</u>																													
LPC																														
ENDRIN/DDT																														

Analyst LF
revision 4/12/01

Verified by nglan

Reviewed by DW
Date 2/12/02

Sequence Name: C:\HPCHEM\1\SEQUENCE\JAN1002.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\jan0902\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line Type	Vial	DataFile	Method	Sample Name
1 Sample	1	DFTPP	GCMS0103	
2 Sample	2	CAL40	GCMS0103	calib mix
3 Sample	21	31216	GCMS0103	gcms scan 12/21
4 Sample	22	31218	GCMS0103	gcms scan 12/21
5 Sample	23	31133	GCMS0103	gcms scan 12/21
6 Sample	24	31217	GCMS0103	gcms scan 12/21
7 Sample	25	31135	GCMS0103	gcms scan 12/21
8 Sample	26	31129A	GCMS0103	gcms scan 12/21
9 Sample	27	31130A	GCMS0103	gcms scan 12/21
10 Sample	28	31219	GCMS0103	gcms scan 12/21
11 Sample	4	LB	GCMS0103	gcms scan 12/24
12 Sample	5	REF	GCMS0103	gcms scan 12/24
13 Sample	6	31250-	GCMS0103	gcms scan 12/24
14 Sample	7	31340-	GCMS0103	gcms scan 12/24
15 Sample	8	31342-	GCMS0103	gcms scan 12/24
16 Sample	9	31275-	GCMS0103	gcms scan 12/24
17 Sample	10	31343	GCMS0103	gcms scan 12/24
18 Sample	11	31345-	GCMS0103	gcms scan 12/24
19 Sample	12	31346-	GCMS0103	gcms scan 12/24
20 Sample	13	31277-	GCMS0103	gcms scan 12/24
21 Sample	14	31276-	GCMS0103	gcms scan 12/24
22 Sample	15	31249-	GCMS0103	gcms scan 12/24
23 Sample	16	31278	GCMS0103	gcms scan 12/24

DFTPP

Data File : C:\HPCHEM\1\DATA\JAN0902\DFTPP.D

Acq On : 10 Jan 2002 8:39 am

Sample :

Misc :

MS Integration Params: RTEINT.P

Method : C:\HPCHEM\1\METHODS\GCMS1126.M (RTE Integrator)


Title : 209 625/TCLP calibration table

Vial: 1

Operator: 209 GALOS

Inst : GC/MS 209

Multiplr: 1.00



Spectrum Information: Average of 26.048 to 26.085 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	43.3	136530	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	53.2	167563	PASS
70	69	0.00	2	0.1	176	PASS
127	198	40	60	49.0	154363	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	315014	PASS
199	198	5	9	6.6	20669	PASS
275	198	10	30	20.5	64451	PASS
365	198	1	100	2.3	7289	PASS
441	443	0.01	100	89.5	31048	PASS
442	198	40	110	57.0	179646	PASS
443	442	17	23	19.3	34686	PASS