

Method 508 Modified Date Extracted 12/19/01Page 77

REF

#GCMs - 16503

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
LB ✓		500ml		1ml			1ml	9A	7			
REF ✓					0.5ml			9B				
* 26532	30983 ✓							6B				
* 26533	30984 ✓							9A				
* 26534	30985 ✓							10B				
* 26536	30986 ✓							10A				
30694	FB ✓							10B				
30695								10A				
30696								9B				
30711								8A				
30712								8B				
30713								7B				
30714								7A				
30605	(re-extracted) ✓							6A				

Remarks

* = New Numbers for FB. eg. (30983).

30605 re-extracted from sample dated 12/18/01.

drip rate must be between 5-15 ml/min for liq-liq extractors

Extraction type -

Solv. lot#

NaCl lot#

Na2SO4 lot#

STANDARD INT STD LOG#

CALIBRATION

REF/SPIKE 91759

MDL/MDRF

Q. CHECK

INT. STD.

SURROGATE 10/901

LPC

ENDRIN/DDT

Reviewed by, C/B 5/15/02

Analyst J. Heath
revision 4/12/01Verified by myReviewed by DR 1/25/02

DFTPP

Data File : C:\HPCHEM\1\DATA\JAN0802\DFTPP.D

Acq On : 8 Jan 2002 10:26 am

Sample :

Misc :

MS Integration Params: GOOFY.P

Method : C:\HPCHEM\1\METHODS\GCMS1126.M (RTE Integrator)


Title : 209 625/TCLP calibration table

Vial: 1

Operator: 209 GALOS

Inst : GC/MS 209

Multiplr: 1.00



Spectrum Information: Average of 26.052 to 26.107 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	47.0	251832	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	56.1	300536	PASS
70	69	0.00	2	0.1	298	PASS
127	198	40	60	52.9	283283	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	535854	PASS
199	198	5	9	6.6	35536	PASS
275	198	10	30	19.7	105482	PASS
365	198	1	100	2.5	13226	PASS
441	443	0.01	100	92.7	41667	PASS
442	198	40	110	43.0	230189	PASS
443	442	17	23	19.5	44937	PASS

Sequence Name: C:\HPCHEM\1\SEQUENCE\jan0802.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\jan0802\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0103	
2	Sample	2	CAL40	GCMS0103	calib mix
3	Sample	3	28355	GCMS0103	
4	Sample	5	28327	GCMS0103	508
5	Sample	6	29466	GCMS0103	
6	Sample	4	28355A	GCMS0103	
7	Sample	7	LB	GCMS0103	gcms scan 12/20
8	Sample	8	REF	GCMS0103	gcms scan 12/20
9	Sample	9	31076	GCMS0103	gcms scan 12/20
10	Sample	10	31073	GCMS0103	gcms scan 12/20
11	Sample	11	31075	GCMS0103	gcms scan 12/20
12	Sample	12	31074	GCMS0103	gcms scan 12/20
13	Sample	13	31017	GCMS0103	gcms scan 12/20
14	Sample	14	31018	GCMS0103	gcms scan 12/20
15	Sample	15	31020	GCMS0103	gcms scan 12/20
16	Sample	16	31024	GCMS0103	gcms scan 12/20
17	Sample	17	30605	GCMS0103	gcms scan 12/20
18	Sample	18	LBB	GCMS0103	gcms scan 12/20a
19	Sample	19	REFB	GCMS0103	gcms scan 12/20a
20	Sample	20	30902	GCMS0103	gcms scan 12/20a
21	Sample	1	DFTPP1	GCMS0103	
22	Sample	21	31022	GCMS0103	gcms scan 12/20a
23	Sample	22	31021	GCMS0103	gcms scan 12/20a
24	Sample	23	30903	GCMS0103	gcms scan 12/20a
25	Sample	24	31019	GCMS0103	gcms scan 12/20a
26	Sample	25	31023	GCMS0103	gcms scan 12/20a
27	Sample	2	CAL40A	GCMS0103	calib mix
28	Sample	26	27288	GCMS0103	gcms scan 12/20a
29	Sample	4	LBC	GCMS0103	gcms scan 12/24
30	Sample	5	REFC	GCMS0103	gcms scan 12/24
31	Sample	6	31278	GCMS0103	gcms scan 12/24
32	Sample	7	31129	GCMS0103	gcms scan 12/24
33	Sample	8	31132	GCMS0103	gcms scan 12/24
34	Sample	9	31255	GCMS0103	gcms scan 12/24
35	Sample	10	31253	GCMS0103	gcms scan 12/24
36	Sample	11	31344	GCMS0103	gcms scan 12/24
37	Sample	12	31254	GCMS0103	gcms scan 12/24
38	Sample	13	31248	GCMS0103	gcms scan 12/24
39	Sample	14	31252	GCMS0103	gcms scan 12/24
40	Sample	15	31130	GCMS0103	gcms scan 12/24
41	Sample	16	31251	GCMS0103	gcms scan 12/24
42	Sample	17	LBD	GCMS0103	gcms scan 12/21

Last Modified: Thu Jan 10 08:16:55 2002

Page: 1