

GC/MS unknown by 508

Method 508Date Extracted 20 Nov 01

Page _____

Alternate Method

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH = 7	Count		
✓	2B	500mL		1mL			1mL					
✓	Ref				1/2mL							
✓	28143									1		
✓	28145									2		
✓	28147									3		
✓	28153									4		
✓	28149									5		
✓	28151									6		
✓	28155									7		
✓	28157									8		
✓	28159									9		
✓	28161									10		
✓	28141									11		
✓	28163									12		
✓	28265									13		
✓	28325									14		
✓	28326									15		
✓	28327									16		
✓	28334									17		
✓	28335									18		
✓	28230									19		
✓	28229									20		
Remarks 1/2mL of Ref 91759 added to "Ref" 28163 reextracted												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type <u>Sep funnel</u> Solv. lot# <u>V 13 E07</u> NaCl lot# <u>V 11599</u> Na2SO4 lot# <u>V 32603</u>												
STANDARD	INT STD LOG#											
CALIBRATION	<u>HY 112601</u>											
REF/SPIKE	<u>91759</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.	<u>15112601</u>											
SURROGATE	<u>24000</u>											
LPC												
ENDRIN/DDT												

Reviewed by

5/11/02

Manny NAW-DEP

Analyst LF/RW
revision 4/12/01

Verified by _____

Reviewed by

Date

12/6/01

Sequence Name: C:\HPCHEM\1\SEQUENCE\NOV2601.S
 Comment: INST 209; 30Mx0.25umX0.50
 Operator: 209 GALOS
 Data Path: C:\HPCHEM\1\data\nov2601\
 Pre-Seq Cmd:
 Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch
 (X) Full Method (X) Inject Anyway
 () Reprocessing Only () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	NP103001	
2	Sample	2	CAL5	NP103001	calib mix
3	Sample	3	CAL10	NP103001	calib mix
4	Sample	4	CAL25	NP103001	calib mix
5	Sample	5	CAL40	NP103001	calib mix
6	Sample	6	CAL80	NP103001	calib mix
7	Sample	7	CAL100	NP103001	calib mix
8	Sample	8	LB	NP103001	10/22
9	Sample	9	LBA	NP103001	11/20 lab blk scan
10	Sample	10	28145	NP103001	11/20 scan
11	Sample	11	28155	NP103001	11/20 scan
12	Sample	12	28141	NP103001	11/20 scan
13	Sample	13	28147	NP103001	11/20 scan
14	Sample	14	28143	NP103001	11/20 scan
15	Sample	15	28149	NP103001	11/20 scan
16	Sample	16	28157	NP103001	11/20 scan
17	Sample	17	28153	NP103001	11/20 scan
18	Sample	18	28151	NP103001	11/20 scan
19	Sample	19	28161	NP103001	11/20 scan
20	Sample	20	28265	NP103001	11/20 scan
21	Sample	21	28159	NP103001	11/20 scan
22	Sample	22	REF	NP103001	11/20 scan
23	Sample	23	LBB	NP103001	11/21 scan
24	Sample	1	DFTPP1	NP103001	
25	Sample	2	24DDD	NP103001	gc/ms scan surrogate 5 ug/ml
26	Sample	24	28269	NP103001	11/21 scan
27	Sample	25	28267	NP103001	11/21 scan
28	Sample	26	28271	NP103001	11/21 scan
29	Sample	27	28273	NP103001	11/21 scan
30	Sample	28	28277	NP103001	11/21 scan
31	Sample	29	28473	NP103001	11/21 scan
32	Sample	30	REFA	NP103001	11/21 scan
33	Sample	2	CAL40A	NP103001	calib
34	Sample	3	28141A	NP103001	11/20 re injected
35	Sample	4	LBC	NP103001	11/23 scan
36	Sample	5	28473A	NP103001	11/23 scan
37	Sample	6	28477	NP103001	11/23 scan
38	Sample	7	28479	NP103001	11/23 scan
39	Sample	8	REFB	NP103001	11/23 scan
40	Sample	9	LBD	NP103001	11/24 scan
41	Sample	10	28682	NP103001	11/24 scan
42	Sample	11	28680	NP103001	11/24 scan
43	Sample	12	REFC	NP103001	11/24 scan

Health Dept