

Method 508 modDate Extracted 11/28/01Page 1

Ref

\$GCMS 15453

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count
LB ✓		500ml		1.0ml			1 ml		7	
Ref ✓					0.5ml					
28860 ✓										
28861 ✓										
28862 ✓										
28897 ✓										
28898 ✓										
28899 ✓										
28900 ✓										
28901 ✓										
28902 ✓										
28650 ✓	re-extract									
Remarks										
28650 was reextracted										
Reviewed by C/B 5/11/02										
drip rate must be between 5-15 ml/min for liq-liq extractors										
drip rate -										
Extraction type - <u>sep fin</u> Solv. lot# <u>V43E07</u> NaCl lot# <u>V3206</u> Na2SO4 lot# <u>V11599</u>										
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	91759									
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE	GCMS 112601									
LPC										
ENDRIN/DDT										

Analyst LF
Revision 4/12/01

Verified by _____

Reviewed by DW
Date 12/16/01