

Method 508 MOD.

Date Extracted 11/30/01

Page _____

Ref

\$60 M-515971

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Sol add	IS add	Final vol.	T.V. pos.	pH	Count
✓	LB	✓ 500ml		1ml			1ml		7	
✓	REF				1ml					
✓	29023									1
	29087									2
✓	29022									3
✓	29141									4
✓	29139									5
✓	29134									6
✓	29190									7
✓	29188									8
Remarks										
Reviewed by, <u>G/B 3/11/02</u>										
drip rate must be between 5-15 ml/min for liq-liq extractors										
drip rate -										
Extraction type <u>SPFume</u>		Solv. lot# <u>V39466</u>		NaCl lot# <u>V11599</u>		Na2SO4 lot# <u>V39400</u>				
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	<u>91759</u>									
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE	<u>24000</u>									
LPC										
ENDRIN/DDT										

Analyst M
revision 4/12/01

Verified by _____

Reviewed by DN
Date 12/17/11

Sequence Name: C:\HPCHEM\1\SEQUENCE\DCE1001.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\dec1001\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

☒ Full Method☒ Inject Anyway☐ Reprocessing Only☐ Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP2	GCMS1126	
2	Sample	2	CAL40	GCMS1126	
3	Sample	3	LB	GCMS1126	gc/ms scan 11/30
4	Sample	4	REF	GCMS1126	gc/ms scan 11/30
5	Sample	1	DFTPP1	GCMS1126	
6	Sample	5	29188	GCMS1126	gc/ms scan 11/30
7	Sample	6	29187	GCMS1126	gc/ms scan 11/30
8	Sample	7	29023	GCMS1126	gc/ms scan 11/30
9	Sample	8	29141	GCMS1126	gc/ms scan 11/30
10	Sample	9	29022	GCMS1126	gc/ms scan 11/30
11	Sample	10	29134	GCMS1126	gc/ms scan 11/30
12	Sample	11	29139	GCMS1126	gc/ms scan 11/30
13	Sample	12	29190	GCMS1126	gc/ms scan 11/30
14	Sample	13	LBB	GCMS1126	gc/ms scan 12/3
15	Sample	14	REFB	GCMS1126	gc/ms scan 12/3
16	Sample	15	29233	GCMS1126	gc/ms scan 12/3
17	Sample	16	29232	GCMS1126	gc/ms scan 12/3
18	Sample	17	29287	GCMS1126	gc/ms scan 12/3
19	Sample	18	29230	GCMS1126	gc/ms scan 12/3
20	Sample	19	29136	GCMS1126	gc/ms scan 12/3
21	Sample	20	29189	GCMS1126	gc/ms scan 12/3
22	Sample	21	29596	GCMS1126	gc/ms scan 12/3
23	Sample	22	29135	GCMS1126	gc/ms scan 12/3
24	Sample	23	29239	GCMS1126	gc/ms scan 12/3
25	Sample	1	DFTPP3	GCMS1126	
26	Sample	2	CAL40A	GCMS1126	
27	Sample	24	29231	GCMS1126	gc/ms scan 12/3
28	Sample	25	LBC	GCMS1126	gc/ms scan 11/30a
29	Sample	26	REFA	GCMS1126	gc/ms scan 11/30a
30	Sample	27	29088	GCMS1126	gc/ms scan 11/30a
31	Sample	28	29027	GCMS1126	gc/ms scan 11/30a
32	Sample	29	29025	GCMS1126	gc/ms scan 11/30a
33	Sample	30	29026	GCMS1126	gc/ms scan 11/30a
34	Sample	31	29024	GCMS1126	gc/ms scan 11/30a
35	Sample	32	29028	GCMS1126	gc/ms scan 11/30a
36	Sample	33	29137	GCMS1126	gc/ms scan 11/30a
37	Sample	34	29029	GCMS1126	gc/ms scan 11/30a
38	Sample	35	29086	GCMS1126	gc/ms scan 11/30a
39	Sample	36	29089	GCMS1126	gc/ms scan 11/30a
40	Sample	37	29140	GCMS1126	gc/ms scan 11/30a
41	Sample	38	29138	GCMS1126	gc/ms scan 11/30a
42	Sample	3	LBD	GCMS1126	gc/ms scan 12/4

Last Modified: Wed Dec 12 10:37:46 2001

Page: 1