

Method 508 ModifiedDate Extracted 12/4/01Page 38

REF

FGCNS14131

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
✓ LB ✓	—	500ML		1 M	0.5 M		1 M		7			
✓ REF ✓	—											
✓ 29291 ✓	—											
✓ 29229 ✓	—											
✓ 29235 ✓	—											
✓ 29295 ✓	—											
✓ 29293 ✓	—											
✓ 29294 ✓	—											
✓ 29296 ✓	—											
✓ 29285 ✓	—											
✓ 29455 ✓	—											
✓ 29456 ✓	—											
Remarks 29285 Very Acidic. 2 ml 6M NaOH was added to bring pH 7.												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type <u>Sep Fun</u>		Solv. lot# <u>V39466</u>	NaCl lot# <u>V1159</u>		Na2SO4 lot# <u>V39400</u>							
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE	<u>91759</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	<u>GMCS 112601</u>											
LPC												
ENDRIN/DDT												

Analyst E. Smith
revision 4/12/01

Verified by _____

Reviewed by DNDate 12/28/11