

Method 508 modDate Extracted 12/7/01Page 51


Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count	rec'd
LB ✓		500 mL		1 mL			1 mL	10A	7		
REP ✓					0.5 mL			10B			
29563 ✓								1B			
29564 ✓								2B			
29666 ✓								1A			
29522 ✓								9B			✓
29670 ✓								9B			
29669 ✓								1B			
29523 ✓								9A			
29668 ✓								1A			
29667 ✓								10A			
29671 ✓								9A			
29734 ✓											
Remarks											
drip rate must be between 5-15 ml/min for liq-liq extractors											
drip rate -											
Extraction type - <u>sep fun</u>		Solv. lot# <u>V39966</u>	NaCl lot# <u>V11599</u>		Na2SO4 lot# <u>V39H00</u>						
STANDARD	INT STD LOG#										
CALIBRATION											
REF/SPIKE	<u>91759</u>										
MDL/MDRF											
Q. CHECK											
INT. STD.											
SURROGATE	<u>GCMS 112601</u>										
LPC											
ENDRIN/DDT											

Analyst LF
revision 4/12/01Verified by phReviewed by DN
Date 1/28/2

DFTPP

Data File : C:\HPCHEM\1\DATA\DEC2701\DFTPP2.D
 Acq On : 29 Dec 2001 1:08 am
 Sample :
 Misc :
 MS Integration Params: GOOFY.P
 Method : C:\HPCHEM\1\METHODS\GCMS1126.M (RTE Integrator)
 Title : 209 625/TCLP calibration table

Vial: 1
 Operator: 209 GALOS
 Inst : GC/MS 209
 Multiplr: 1.00



Spectrum Information: Average of 26.214 to 26.288 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	46.8	343895	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	57.5	422513	PASS
70	69	0.00	2	0.0	28	PASS
127	198	40	60	53.0	389202	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	734309	PASS
199	198	5	9	6.7	49165	PASS
275	198	10	30	20.5	150823	PASS
365	198	1	100	2.7	19826	PASS
441	443	0.01	100	93.9	65199	PASS
442	198	40	110	49.0	359580	PASS
443	442	17	23	19.3	69441	PASS

Sequence Name: C:\HPCHEM\1\SEQUENCE\Dec2701.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\dec2701\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

(X) Full Method

(X) Inject Anyway

() Reprocessing Only

() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS1126	
2	Sample	2	CAL100	GCMS1126	
3	Sample	3	LB	GCMS1126	GC/MS SCAN 12/12
4	Sample	4	LBA	GCMS1126	GC/MS SCAN 12/13
5	Sample	4	REF	GCMS1126	GC/MS SCAN 12/13
6	Sample	5	30425	GCMS1126	GC/MS SCAN 12/13
7	Sample	6	30354	GCMS1126	GC/MS SCAN 12/13
8	Sample	35	SURR	GCMS1126	GC/MS SCAN 12/13
9	Sample	7	30355	GCMS1126	GC/MS SCAN 12/13
10	Sample	8	30251	GCMS1126	GC/MS SCAN 12/13
11	Sample	9	30023	GCMS1126	GC/MS SCAN 12/13
12	Sample	10	26477	GCMS1126	GC/MS SCAN 12/13
13	Sample	11	26827	GCMS1126	GC/MS SCAN 12/13
14	Sample	12	30009	GCMS1126	GC/MS SCAN 12/12A
15	Sample	13	30010	GCMS1126	GC/MS SCAN 12/12A
16	Sample	14	30110	GCMS1126	GC/MS SCAN 12/12A
17	Sample	15	30111	GCMS1126	GC/MS SCAN 12/12A
18	Sample	16	30107	GCMS1126	GC/MS SCAN 12/12A
19	Sample	17	30353	GCMS1126	GC/MS SCAN 12/13
20	Sample	18	30357	GCMS1126	GC/MS SCAN 12/13
21	Sample	19	30352	GCMS1126	GC/MS SCAN 12/13
22	Sample	20	30351	GCMS1126	GC/MS SCAN 12/13
23	Sample	21	30356	GCMS1126	GC/MS SCAN 12/13
24	Sample	1	DFTPP1	GCMS1126	
25	Sample	22	30252	GCMS1126	GC/MS SCAN 12/13
26	Sample	2	CAL100A	GCMS1126	
27	Sample	23	30426	GCMS1126	GC/MS SCAN 12/13
28	Sample	24	LBB	GCMS1126	GC/MS SCAN 12/17
29	Sample	25	30631	GCMS1126	GC/MS SCAN 12/17
30	Sample	26	30632	GCMS1126	GC/MS SCAN 12/17
31	Sample	27	30633	GCMS1126	GC/MS SCAN 12/17
32	Sample	28	30634	GCMS1126	GC/MS SCAN 12/17
33	Sample	29	30601	GCMS1126	GC/MS SCAN 12/17
34	Sample	30	30602	GCMS1126	GC/MS SCAN 12/17
35	Sample	31	30603	GCMS1126	GC/MS SCAN 12/17
36	Sample	32	30608	GCMS1126	GC/MS SCAN 12/17
37	Sample	33	30607	GCMS1126	GC/MS SCAN 12/17
38	Sample	34	30604	GCMS1126	GC/MS SCAN 12/17
39	Sample	36	29563	GCMS1126	GC/MS SCAN 12/7
40	Sample	37	29523	GCMS1126	GC/MS SCAN 12/7
41	Sample	1	DFTPP2	GCMS1126	
42	Sample	4	REFA	GCMS1126	GC/MS SCAN 12/13

Last Modified: Fri Dec 28 16:44:24 2001

Page: 1

Sequence Name: C:\HPCHEM\1\SEQUENCE\DEC3101.S

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	42	30906	GCMS1126	gcms scan 12/19a
45	Sample	43	30907	GCMS1126	gcms scan 12/19a
46	Sample	44	30904	GCMS1126	gcms scan 12/19a
47	Sample	45	30905	GCMS1126	gcms scan 12/19a
48	Sample	1	DFTPP2	GCMS1126	
49	Sample	46	30705	GCMS1126	gcms scan 12/19a
50	Sample	47	30908	GCMS1126	gcms scan 12/19a
51	Sample	48	30935	GCMS1126	gcms scan 12/19a
52	Sample	2	CAL40B	GCMS1126	
53	Sample	49	30936	GCMS1126	gcms scan 12/19a
54	Sample	50	30933	GCMS1126	gcms scan 12/19a
55	Sample	51	30934	GCMS1126	gcms scan 12/19a
56	Sample	52	30909	GCMS1126	gcms scan 12/19a
57	Sample	53	REFG	GCMS1126	gcms scan 12/17
58	Sample	54	REFH	GCMS1126	gcms scan 12/17
59	Sample	55	REFI	GCMS1126	gcms scan 12/13
60	Sample	56	29667	GCMS1126	gcms scan 12/7
61	Sample	57	30257	GCMS1126	gcms scan 12/12
62	Sample	58	30116	GCMS1126	gcms scan 12/12
63	Sample	59	LBF	GCMS1126	gcms scan 12/7
64	Sample	60	REFF	GCMS1126	gcms scan 12/7
65	Sample	61	29669	GCMS1126	gcms scan 12/7
66	Sample	62	29564	GCMS1126	gcms scan 12/7
67	Sample	63	29522	GCMS1126	gcms scan 12/7
68	Sample	64	29666	GCMS1126	gcms scan 12/7
69	Sample	65	29668	GCMS1126	gcms scan 12/7

Sequence Name: C:\HPCHEM\1\SEQUENCE\jan0302.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\jan0302\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

☒ Full Method☐ Reprocessing Only

On A Barcode Mismatch

☒ Inject Anyway☐ Don't Inject

Line Type	Vial	DataFile	Method	Sample Name
1 Sample	1	DFTPP	GCMS1126	
2 Sample	2	CAL5	GCMS1126	mix
3 Sample	3	CAL10	GCMS1126	mix
4 Sample	4	CAL25	GCMS1126	mix
5 Sample	5	CAL40	GCMS1126	mix
6 Sample	6	CAL80	GCMS1126	mix
7 Sample	7	CAL100	GCMS1126	mix
8 Sample	8	CALCHK	GCMS1126	20
9 Sample	9	SURR	GCMS1126	
10 Sample	10	30117	GCMS1126	gcms scan 12/12
11 Sample	11	30118	GCMS1126	gcms scan 12/12
12 Sample	12	29791	GCMS1126	gcms scan 12/10
13 Sample	13	29671	GCMS1126	gcms scan 12/7
14 Sample	14	29734	GCMS1126	gcms scan 12/7
15 Sample	15	29670	GCMS1126	gcms scan 12/7
16 Sample	16	29520	GCMS1126	gcms scan
17 Sample	17	29517	GCMS1126	gcms scan
18 Sample	18	LB	GCMS1126	gcms scan 1/2
19 Sample	19	REF	GCMS1126	gcms scan 1/2
20 Sample	20	29228	GCMS1126	gcms scan 1/2
21 Sample	21	29562	GCMS1126	gcms scan 1/2
22 Sample	22	29519	GCMS1126	gcms scan 1/2
23 Sample	23	29292	GCMS1126	gcms scan 1/2
24 Sample	24	29284	GCMS1126	gcms scan 1/2
25 Sample	25	29343	GCMS1126	gcms scan 1/2
26 Sample	26	29344	GCMS1126	gcms scan 1/2
27 Sample	27	30632	GCMS1126	gcms scan 1/2
28 Sample	28	29465	GCMS1126	gcms scan 1/2
29 Sample	29	29286	GCMS1126	gcms scan 1/2
30 Sample	30	30631	GCMS1126	gcms scan 1/2
31 Sample	31	LBB	GCMS1126	gcms scan 1/2a
32 Sample	32	REFB	GCMS1126	gcms scan 1/2a
33 Sample	33	31642	GCMS1126	gcms scan 1/2a
34 Sample	34	30111	GCMS1126	gcms scan 1/2a
35 Sample	35	30253	GCMS1126	gcms scan 1/2a
36 Sample	36	31640	GCMS1126	gcms scan 1/2a
37 Sample	37	30023	GCMS1126	gcms scan 1/2a
38 Sample	38	31586	GCMS1126	gcms scan 1/2a
39 Sample	39	31641	GCMS1126	gcms scan 1/2a
40 Sample	40	31633	GCMS1126	gcms scan 1/2a
41 Sample	41	30010	GCMS1126	gcms scan 1/2a
42 Sample	42	31632	GCMS1126	gcms scan 1/2a

Last Modified: Thu Jan 03 16:34:24 2002

Page: 1