

Method 508 mod

Date Extracted 1/02/02

Page 2

ref

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
LB		500ml		1 mL			1 mL	4B	7			
Ref					0.5 mL			1A				
31549								2B				
31587								3B				
31499								3A				
31498								4B				
31551								2A				
31591								2B				
31589								1B				
31590								1B				
31588								1A				
31592								2A				
31630								2B				
Remarks												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type - <u>sep fun</u> Solv. lot# <u>V39466</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V39H00</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE	<u>91759</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	<u>119101</u>											
LPC												
ENDRIN/DDT												

Analyst LF
revision 4/12/01

Verified by _____

Reviewed by _____
Date _____

Sequence Name: C:\HPCHEM\1\SEQUENCE\JAN2402.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\jan2402\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0103	
2	Sample	2	CAL25	GCMS0103	
3	Sample	16	31134	GCMS0103	REINJECT
4	Sample	17	31218	GCMS0103	REINJECT
5	Sample	3	LB	GCMS0103	GC/MS SCAN 1/3
6	Sample	31	SURR	GCMS0103	GC/MS SCAN 1/3
7	Sample	4	REF	GCMS0103	GC/MS SCAN 1/3
8	Sample	5	31663	GCMS0103	GC/MS SCAN 1/3
9	Sample	6	7	GCMS0103	GC/MS SCAN 1/3
10	Sample	7	31662	GCMS0103	GC/MS SCAN 1/3
11	Sample	8	6	GCMS0103	GC/MS SCAN 1/3
12	Sample	9	102	GCMS0103	GC/MS SCAN 1/3
13	Sample	10	8	GCMS0103	GC/MS SCAN 1/3
14	Sample	11	9	GCMS0103	GC/MS SCAN 1/3
15	Sample	12	5	GCMS0103	GC/MS SCAN 1/3
16	Sample	13	10	GCMS0103	GC/MS SCAN 1/3
17	Sample	14	107	GCMS0103	GC/MS SCAN 1/3
18	Sample	15	104	GCMS0103	GC/MS SCAN 1/3
19	Sample	18	LBA	GCMS0103	GC/MS SCAN 1/2
20	Sample	19	REFA	GCMS0103	GC/MS SCAN 1/2
21	Sample	20	31549	GCMS0103	GC/MS SCAN 1/2
22	Sample	21	31587	GCMS0103	GC/MS SCAN 1/2
23	Sample	22	31499	GCMS0103	GC/MS SCAN 1/2
24	Sample	23	31498	GCMS0103	GC/MS SCAN 1/2
25	Sample	24	31551	GCMS0103	GC/MS SCAN 1/2
26	Sample	25	31591	GCMS0103	GC/MS SCAN 1/2
27	Sample	26	31589	GCMS0103	GC/MS SCAN 1/2
28	Sample	27	31590	GCMS0103	GC/MS SCAN 1/2
29	Sample	28	31588	GCMS0103	GC/MS SCAN 1/2
30	Sample	29	31592	GCMS0103	GC/MS SCAN 1/2
31	Sample	30	31630	GCMS0103	GC/MS SCAN 1/2

Sequence Name: C:\HPCHEM\1\SEQUENCE\JAN2502.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\jan2502\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

☒ Full Method☐ Reprocessing Only

On A Barcode Mismatch

☒ Inject Anyway☐ Don't Inject

Line Type	Vial	DataFile	Method	Sample Name
1 Sample	1	DFTPP	GCMS0103	
2 Sample	2	CAL25	GCMS0103	
3 Sample	3	30936A	GCMS0103	REINJECT
4 Sample	23	31498	GCMS0103	GC/MS SCAN 1/2
5 Sample	24	31551	GCMS0103	GC/MS SCAN 1/2
6 Sample	25	31591	GCMS0103	GC/MS SCAN 1/2
7 Sample	26	31589	GCMS0103	GC/MS SCAN 1/2
8 Sample	27	31590	GCMS0103	GC/MS SCAN 1/2
9 Sample	28	31588	GCMS0103	GC/MS SCAN 1/2
10 Sample	29	31592	GCMS0103	GC/MS SCAN 1/2
11 Sample	30	31630	GCMS0103	GC/MS SCAN 1/2
12 Sample	32	REF	GCMS0103	REINJECT
13 Sample	33	30642	GCMS0103	REINJECT
14 Sample	34	30649	GCMS0103	REINJECT
15 Sample	35	30651	GCMS0103	REINJECT
16 Sample	36	30704	GCMS0103	REINJECT

Sequence Name: C:\HPCHEM\1\SEQUENCE\FEB2102.S

Comment:

Operator:

Data Path: C:\HPCHEM\1\DATA\feb2102\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0208	
2	Sample	2	CAL40	GCMS0208	
3	Sample	29	SURR	GCMS0208	
4	Sample	3	LB	GCMS0208	gc/ms scan 1/22
5	Sample	4	REF	GCMS0208	gc/ms scan 1/22
6	Sample	5	1329	GCMS0208	gc/ms scan 1/22
7	Sample	6	472	GCMS0208	gc/ms scan 1/22
8	Sample	7	1435	GCMS0208	gc/ms scan 1/22
9	Sample	8	1436	GCMS0208	gc/ms scan 1/22
10	Sample	9	1437	GCMS0208	gc/ms scan 1/22
11	Sample	10	1438	GCMS0208	gc/ms scan 1/22
12	Sample	11	1439	GCMS0208	gc/ms scan 1/22
13	Sample	12	1440	GCMS0208	gc/ms scan 1/22
14	Sample	13	1441	GCMS0208	gc/ms scan 1/22
15	Sample	14	1442	GCMS0208	gc/ms scan 1/22
16	Sample	15	1462	GCMS0208	gc/ms scan 1/22
17	Sample	1	DFTPP2	GCMS0208	gc/ms scan 1/22
18	Sample	2	CAL40A	GCMS0208	gc/ms scan 1/22
19	Sample	16	LBA	GCMS0208	gc/ms scan 1/22
20	Sample	17	REFA	GCMS0208	gc/ms scan 1/22
21	Sample	18	1463	GCMS0208	gc/ms scan 1/22
22	Sample	19	1464	GCMS0208	gc/ms scan 1/22
23	Sample	20	1465	GCMS0208	gc/ms scan 1/22
24	Sample	21	1498	GCMS0208	gc/ms scan 1/22
25	Sample	22	1499	GCMS0208	gc/ms scan 1/22
26	Sample	23	1500	GCMS0208	gc/ms scan 1/22
27	Sample	24	1501	GCMS0208	gc/ms scan 1/22
28	Sample	25	1502	GCMS0208	gc/ms scan 1/22
29	Sample	26	1503	GCMS0208	gc/ms scan 1/22
30	Sample	27	1504	GCMS0208	gc/ms scan 1/22
31	Sample	28	1505	GCMS0208	gc/ms scan 1/22
32	Sample	30	LBB	GCMS0208	gc/ms scan 2/21
33	Sample	31	REFB	GCMS0208	gc/ms scan 2/21
34	Sample	32	REFB2	GCMS0208	gc/ms scan 2/21
35	Sample	33	31549	GCMS0208	gc/ms scan 2/21
36	Sample	34	31587	GCMS0208	gc/ms scan 2/21
37	Sample	35	31499	GCMS0208	gc/ms scan 2/21
38	Sample	36	31498	GCMS0208	gc/ms scan 2/21
39	Sample	37	31551	GCMS0208	gc/ms scan 2/21
40	Sample	38	31591	GCMS0208	gc/ms scan 2/21
41	Sample	39	31589	GCMS0208	gc/ms scan 2/21
42	Sample	40	31590	GCMS0208	gc/ms scan 2/21

Last Modified: Fri Feb 22 17:04:26 2002

Page: 1

* 43 Sample 41 31588 GCMS0208 gc/ms scan 2/21 |

Sequence Name: C:\HPCHEM\1\SEQUENCE\FEB2102.S

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	42	31592	GCMS0208	gc/ms scan 2/21
45	Sample	43	31630	GCMS0208	gc/ms scan 2/21 fb
46	Sample	1	DFTPP3	GCMS0208	
47	Sample	2	CAL40B	GCMS0208	
48	Sample	44	LBC	GCMS0208	gc/ms scan 1/14
49	Sample	45	REFC	GCMS0208	gc/ms scan 1/14
50	Sample	46	REFC2	GCMS0208	gc/ms scan 1/14
51	Sample	47	3385	GCMS0208	gc/ms scan 1/14
52	Sample	48	3386	GCMS0208	gc/ms scan 1/14
53	Sample	49	3387	GCMS0208	gc/ms scan 1/14
54	Sample	50	3388	GCMS0208	gc/ms scan 1/14
55	Sample	51	3389	GCMS0208	gc/ms scan 1/14
56	Sample	52	3390	GCMS0208	gc/ms scan 1/14
57	Sample	53	3391	GCMS0208	gc/ms scan 1/14
58	Sample	54	3434	GCMS0208	gc/ms scan 1/14
59	Sample	55	3435	GCMS0208	gc/ms scan 1/14
60	Sample	56	3436	GCMS0208	gc/ms scan 1/14
61	Sample	57	LBD	GCMS0208	gc/ms scan 2/21
62	Sample	58	REFD	GCMS0208	gc/ms scan 2/21
63	Sample	59	REFD2	GCMS0208	gc/ms scan 2/21
64	Sample	60	938	GCMS0208	gc/ms scan 2/21
65	Sample	61	939	GCMS0208	gc/ms scan 2/21
66	Sample	62	940	GCMS0208	gc/ms scan 2/21
67	Sample	63	941	GCMS0208	gc/ms scan 2/21
68	Sample	64	942	GCMS0208	gc/ms scan 2/21
69	Sample	65	1143	GCMS0208	gc/ms scan 2/21
70	Sample	66	1153	GCMS0208	gc/ms scan 2/21
71	Sample	67	1155	GCMS0208	gc/ms scan 2/21
72	Sample	68	1154	GCMS0208	gc/ms scan 2/21
73	Sample	69	1137	GCMS0208	gc/ms scan 2/21
74	Sample	1	DFTPP4	GCMS0208	
75	Sample	2	CAL40C	GCMS0208	
76	Sample	70	LBE	GCMS0208	gc/ms scan 2/22
77	Sample	71	REFE	GCMS0208	gc/ms scan 2/22
78	Sample	72	REFE2	GCMS0208	gc/ms scan 2/22
79	Sample	73	3833	GCMS0208	gc/ms scan 2/22
80	Sample	74	3834	GCMS0208	gc/ms scan 2/22
81	Sample	75	3835	GCMS0208	gc/ms scan 2/22 fb
82	Sample	76	3914	GCMS0208	gc/ms scan 2/22
83	Sample	77	3915	GCMS0208	gc/ms scan 2/22
84	Sample	78	3916	GCMS0208	gc/ms scan 2/22 fb
85	Sample	79	4024	GCMS0208	gc/ms scan 2/22
86	Sample	80	4025	GCMS0208	gc/ms scan 2/22
87	Sample	81	4026	GCMS0208	gc/ms scan 2/22
88	Sample	82	4027	GCMS0208	gc/ms scan 2/22

DFTPP

Data File : C:\HPCHEM\1\DATA\FEB2102\DFTPP2.D

Acq On : 22 Feb 2002 4:52 am

Sample : gc/ms scan 1/22

Misc :

MS Integration Params: GOOFY.P

Method : C:\HPCHEM\1\METHODS\GCMS0208.M (RTE Integrator)


Title : 209 625/TCLP calibration table

Vial: 1

Operator:

Inst : GC/MS 209

Multiplr: 1.00



Spectrum Information: Average of 26.912 to 26.977 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	40.4	112716	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	52.1	145189	PASS
70	69	0.00	2	0.1	95	PASS
127	198	40	60	45.8	127743	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	278805	PASS
199	198	5	9	6.7	18558	PASS
275	198	10	30	22.5	62635	PASS
365	198	1	100	2.6	7382	PASS
441	443	0.01	100	94.0	27216	PASS
442	198	40	110	54.0	150626	PASS
443	442	17	23	19.2	28962	PASS