

Method 508 MOD.

Date Extracted 2-26-02

Page 108

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count
L B		500 ml		1 ml			1 ml	6A	7	
Rd 1					1 ml			1A		
Rd 2					1 ml			1B		
1335								2B		
1548								2A		
3438								3A		
3439								7A		
3440								7B		
3441	FB							8A		
3454	FB, Q							3B		
3471								8B		
3472		450 ml						4A		
917		500 ml						4B		
Remarks - NaOH WAS ADDED TO 3454 TO BRING PH TO SEVEN - - 3472 WAS ONLY 450 mL Reviewed by, <u>CB 6/5/02</u>										
drip rate must be between 5-15 ml/min for liq-liq extractors										
Extraction type - <u>SEP. Fub</u> Solv. lot# <u>NaCl lot# 39618</u> Na2SO4 lot# <u>46623</u>										
STANDARD	INT STD LOG#									
CALIBRATION										
REF SPIKE	2193									
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE	21602									
LPC										
ENDRIN/DDT										

Analyst A.D.
revision 4/12/01

Verified by _____

Reviewed by 3/5/12
Date