

Method 508 MOD.Date Extracted 3-18-02Page 11

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
L.B.		500ml		1ml			1ml	8B	7			
Ref 1					1ml			7A				
Ref 2					1ml			8A				
6303								10A				
6304	F.B.							4A				
6305								4B				
6306								10B				
6307								9B				
6308	F.B.							7B				
6333								9A				
6334								6A				
6335								10B				
6336								4B				
Remarks												
Reviewed by, <u>JB 3/6/02</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type - <u>Sep. Fun</u> Solv. lot# <u>M 17262</u> NaCl lot# <u>V39618</u> Na2SO4 lot# <u>V46623</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	2193											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	2160 A											
LPC												
ENDRIN/DDT												

Analyst A.D.
Revision 4/12/01

Verified by _____

Reviewed by DN
Date 3/18/02