

westchester county labs

Folder: RUNS2

Page 497

Protocol: hg-cv

POST-RUN REPORT

Line	Conc.	Units	SD/RSD	1	2	3	4	5
*** Standard: 1 Rep: 1								
				Seq: 749		15:21:01	05 Apr 02	HG
Hg	.000	ppb	-11297					
*** Standard: 2 Rep: 1								
				Seq: 750		15:23:03	05 Apr 02	HG
Hg	.400	ppb	15939					
*** Standard: 3 Rep: 1								
				Seq: 751		15:24:42	05 Apr 02	HG
Hg	1.00	ppb	57207					
*** Standard: 4 Rep: 1								
				Seq: 752		15:26:29	05 Apr 02	HG
Hg	3.00	ppb	189519					
*** Standard: 5 Rep: 1								
				Seq: 753		15:28:09	05 Apr 02	HG
Hg	5.00	ppb	324829					
*** Standard: 6 Rep: 1								
				Seq: 754		15:29:52	05 Apr 02	HG
Hg	10.0	ppb	641848					
*** Check Standard: 2 Ck21.0 qc								
Line	Flag	%Rcv.	Found	True	Units	SD/RSD		
Hg		93.6	.936	1.00	ppb	.000		
*** Sample ID: 7765 pot								
				Seq: 756		15:33:43	05 Apr 02	HG
Hg	-.021	ppb	.000	-.021				
*** Sample ID: 7765 spk								
				Seq: 757		15:35:21	05 Apr 02	HG
Hg	5.12	ppb	.000	5.12				
*** Sample ID: 7765d								
				Seq: 758		15:37:10	05 Apr 02	HG
Hg	-.056	ppb	.000	-.056				
*** Sample ID: 7765d spk								
				Seq: 759		15:38:49	05 Apr 02	HG
Hg	5.23	ppb	.000	5.23				
*** Sample ID: 7766								
				Seq: 760		15:40:38	05 Apr 02	HG
Hg	-.030	ppb	.000	-.030				
*** Sample ID: 7767								
				Seq: 761		15:42:19	05 Apr 02	HG
Hg	-.030	ppb	.000	-.030				
*** Sample ID: 7768								
				Seq: 762		15:43:55	05 Apr 02	HG
Hg	-.028	ppb	.000	-.028				

Reviewed by,

westchester county labs

Folder: RUNS2

Page 499

Protocol: hg-cv

POST-RUN REPORT

Line	Conc.	Units	SD/RSD	1	2	3	4	5
*** Sample ID: 7978d spl ¹								
Hg	5.41	ppb	.000	5.41				
Seq: 776 16:10:11 05 Apr 02 HG								
*** Check Standard: 2 Ck21.0 qc								
Line	Flag	%Rcv.	Found	True	Units	SD/RSD		
Hg		110.7	1.10	1.00	ppb	.000		
Seq: 777 16:11:58 05 Apr 02 HG								
*** Sample ID: 8012								
Hg	-.060	ppb	.000	-.060				
Seq: 778 16:13:49 05 Apr 02 HG								
*** Sample ID: 8013								
Hg	-.027	ppb	.000	-.027				
Seq: 779 16:15:47 05 Apr 02 HG								
*** Sample ID: 8014								
Hg	-.044	ppb	.000	-.044				
Seq: 780 16:17:29 05 Apr 02 HG								
*** Sample ID: 8016								
Hg	-.042	ppb	.000	-.042				
Seq: 781 16:19:15 05 Apr 02 HG								
*** Sample ID: 8017								
Hg	-.030	ppb	.000	-.030				
Seq: 782 16:20:51 05 Apr 02 HG								
*** Sample ID: 8018								
Hg	-.028	ppb	.000	-.028				
Seq: 783 16:22:27 05 Apr 02 HG								
*** Sample ID: 8063								
Hg	-.023	ppb	.000	-.023				
Seq: 784 16:24:45 05 Apr 02 HG								
*** Sample ID: 8064								
Hg	-.025	ppb	.000	-.025				
Seq: 785 16:26:25 05 Apr 02 HG								
*** Sample ID: 8065								
Hg	-.029	ppb	.000	-.029				
Seq: 786 16:28:02 05 Apr 02 HG								
*** Sample ID: 8066								
Hg	-.040	ppb	.000	-.040				
Seq: 787 16:29:39 05 Apr 02 HG								
*** Check Standard: 2 Ck21.0 qc								
Line	Flag	%Rcv.	Found	True	Units	SD/RSD		
Hg		112.7	1.12	1.00	ppb	.000		
Seq: 788 16:31:47 05 Apr 02 HG								

westchester county labs

Folder: RUNS2

Page 500

Protocol: hg-cv

POST-RUN REPORT

Line	Conc.	Units	SD/RSD	1	2	3	4	5
*** Sample ID: 8067								
				Seq: 789	16:33:25	05 Apr 02	HG	
Hg	.017	ppb	.000	.017				
*** Sample ID: 8068								
				Seq: 790	16:35:15	05 Apr 02	HG	
Hg	-.053	ppb	.000	-.053				
*** Sample ID: 8069								
				Seq: 791	16:37:38	05 Apr 02	HG	
Hg	-.079	ppb	.000	-.079				
*** Check Standard: 2 Ck21.0 qc								
				Seq: 799	16:52:35	05 Apr 02	HG	
Line	Flag	%Rcv.	Found	True	Units	SD/RSD		
Hg		109. ✓	1.09	1.00	ppb	.000		

WinHg Database 1.3

FileUtilityHelp

RN<RN>?

Protocolhg-cvDataset/ProtoRUNS2/hg-cv

ProtocolLine infoCal CurveReportCtrl ChartViewer

Reset

Calib Coeffs

New Cal

Update Coeffs

Spike Coeffs

A

B1.52895e-5

C1.29448e-1

Rho.999891

TypeLinear

IncludeS6Rep 1☒

Rel. Abs.641848

Accepted

☒ Calibrated

☒ Accepted

Accept

☐ New

S	Conc.	Calc.	Dev.	Mean	SD or %RSD	Rep 1	Rep 2	Rep 3
01	.00000	-.043	-.043	-11298	1	-11297		
02	.40000	.373	-.027	15939	0%	15939		
03	1.0000	1.00	.004	57208	0%	57207		
04	3.0000	3.03	.027	189519	0%	189519		
05	5.0000	5.10	.096	324830	0%	324829		
06	10.000	9.94	-.057	641848	0%	641848		

ReadyNUM