

508 MOD.

Date Extracted

4/4/02

Page

44

ref

\$6CMS-18024

Analyst A.D.
revision 4/12/01

Verified by

Reviewed by _____
Date _____

Sequence Name: C:\HPCHEM\1\SEQUENCE\APR1102.S

Comment:

Operator:

Data Path: C:\HPCHEM\1\DATA\apr1102\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

(X) Full Method

() Reprocessing Only

On A Barcode Mismatch

(X) Inject Anyway


() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0408	
2	Sample	2	CAL5	GCMS0408	
3	Sample	3	CAL10	GCMS0408	
4	Sample	4	CAL25	GCMS0408	
5	Sample	5	CAL40	GCMS0408	
6	Sample	6	CAL80	GCMS0408	
7	Sample	7	CAL100	GCMS0408	
8	Sample	7	CAL100A	GCMS0408	really
9	Sample	26	CAL25DM	GCMS0408	4/1
10	Sample	27	CAL40DM	GCMS0408	4/1
11	Sample	8	CAL40DM2	GCMS0408	4/10
12	Sample	29	CAL40MAR	GCMS0408	4/10
13	Sample	28	CAL80DM	GCMS0408	4/1
14	Sample	9	LB	GCMS0408	gc/ms scan 4/3
15	Sample	10	REF1	GCMS0408	gc/ms scan 4/3
16	Sample	11	REF2	GCMS0408	gc/ms scan 4/3
17	Sample	12	7636	GCMS0408	gc/ms scan 4/3
18	Sample	13	7637	GCMS0408	gc/ms scan 4/3
19	Sample	14	7765	GCMS0408	gc/ms scan 4/3
20	Sample	15	LBA	GCMS0408	gc/ms scan 4/4
21	Sample	16	REFA1	GCMS0408	gc/ms scan 4/4
22	Sample	17	REFA2	GCMS0408	gc/ms scan 4/4
23	Sample	18	8012	GCMS0408	gc/ms scan 4/4
24	Sample	19	8013	GCMS0408	gc/ms scan 4/4
25	Sample	20	8014	GCMS0408	gc/ms scan 4/4
26	Sample	21	8015	GCMS0408	gc/ms scan 4/4
27	Sample	22	8016	GCMS0408	gc/ms scan 4/4
28	Sample	23	8017	GCMS0408	gc/ms scan 4/4
29	Sample	24	8018	GCMS0408	gc/ms scan 4/4
30	Sample	25	8019	GCMS0408	gc/ms scan 4/4

DFTPP

Data File : C:\HPCHEM\1\DATA\APR1102\DFTPP.D
 Acq On : 11 Apr 2002 9:29 am
 Sample :
 Misc :
 MS Integration Params: GOOFY.P
 Method : C:\HPCHEM\1\METHODS\GCMS0412.M (RTE Integrator)
 Title : 209 625/TCLP calibration table

Vial: 1
 Operator:
 Inst : 209
 Multiplr: 1.00



Spectrum Information: Scan 2370

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	34.3	132032	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	46.1	177728	PASS
70	69	0.00	2	0.0	0	PASS
127	198	40	60	43.0	165568	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	385280	PASS
199	198	5	9	6.5	25000	PASS
275	198	10	30	23.1	88984	PASS
365	198	1	100	2.4	9228	PASS
441	443	0.01	100	97.1	50280	PASS
442	198	40	110	71.8	276800	PASS
443	442	17	23	18.7	51800	PASS