

Method 508 mod

Date Extracted 4-10-02

Page 54

REF

\$GEMS-18296

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Sol add	IS add	Final vol.	T.V. pos.	pH	Count		
LB		500ml		1ml			1ml	1B	7			
REF ₁					1ml			4B				
REF ₂					↓			1A				
8482 ✓								7A				
8483 ✓								4A				
8484 ✓								2A				
8485								8A				
8486 ✓								7B				
8487								6A				
8488 ✓ FB								6B				
8489 ✓ FB								2B				
Remarks												
Reviewed by, <u>JP 9/16/02</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type <u>SEP Fun</u>		Solv. lot# <u>X05E07</u>		NaCl lot# <u>V39614</u>		Na2SO4 lot# <u>V46648</u>						
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE	<u>2193</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	<u>New Surrogate</u>											
LPC												
ENDRIN/DDT												

Analyst E.H.
Revision 4/12/01

Verified by _____

Reviewed by DN
Date 5/7/02

Sequence Name: C:\MSDchem\1\sequence\020419.S
Comment: GCMS scan
Operator: Nefliu
Data Path: C:\MSDchem\1\DATA\020419\
Pre-Seq Cmd:
Post-Seq Cmd:

Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line Type	Vial	DataFile	Method	Sample Name
1 Sample	1	DFTPP	508FAST	dftpp
2 Sample	2	0419C3	508SCAN1	ccv
3 Sample	3	0409LB	508SCAN1	04/09 lab blk
4 Sample	4	0410LB	508SCAN1	04/10 lab blk
5 Sample	5	0409REF1	508SCAN1	04/09 ref
6 Sample	6	0410REF1	508SCAN1	04/10 ref
7 Sample	7	0409REF2	508SCAN1	04/09 ref
8 Sample	8	0410REF2	508SCAN1	04/10 ref
9 Sample	9	08271	508SCAN1	sample
10 Sample	10	08323	508SCAN1	sample
11 Sample	11	08324	508SCAN1	sample
12 Sample	12	08425	508SCAN1	sample
13 Sample	13	08426	508SCAN1	sample
14 Sample	14	08427	508SCAN1	sample
15 Sample	15	08428	508SCAN1	sample
16 Sample	16	08429	508SCAN1	sample
17 Sample	17	08430	508SCAN1	sample
18 Sample	18	08431	508SCAN1	sample
19 Sample	19	08482	508SCAN1	sample
20 Sample	20	08483	508SCAN1	sample
21 Sample	21	08484	508SCAN1	sample
22 Sample	22	08485	508SCAN1	sample
23 Sample	23	08486	508SCAN1	sample
24 Sample	24	08487	508SCAN1	sample
25 Sample	25	08488	508SCAN1	sample
26 Sample	26	08489	508SCAN1	sample

Method 508 GC MS Scan
Instrument 214

File Directory 020419

Storage Directory _____

Lab #	Filename	Description	Vial #		Inj. ul	Standards used - Lot #			Int. Std.	Remarks	Initials
			front	rear		Calibration	Reference	Surrogate			
	A77PA	d7pp	1		1						MN
	041903	cal	2			2206/22010/2206					
	0409LB	cal	3								
	0410LB		4								
	0409Ry1	cal	5				293				
	0410Ry1		6								
	0409Ry2		7								
	0410Ry2		8								
	08271	sample	9								
	08323		10								
	08324		11								
	08425		12								
	08426		13								
	08427		14								
	08428		15								
	08429		16								
	08430		17								
	08431		18								
	08432		19								
	08433		20								
	08434		21								
	08435		22								
	08436		23								
	08437		24								
	08438		25								
	08439		26								

Date 04/19/02

Extraction date 04/10, 09/02

Page 4