

Method 508 mod

Date Extracted 1/15/02

Page

ref

#GCMS-16990

see 2/22
On
previous
batch
check
W/B.

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
LB		500 mL		1 mL			1 mL					
Ref					0.5 mL							
364.	re-ext	↓										
945.	re-ext	↓										
830.		↓										
831.		490 mL										
832.		500 mL										
833.		↓										
1036.		↓										
1037.		↓										
1038.		↓										
1039.		↓										
1040.		↓										
1042.		↓										
Remarks												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type -		Solv. lot#		NaCl lot#		Na2SO4 lot#						
STANDARD	INT STD LOG#											
CALIBRATION	HP013002											
REF/SPIKE												
MDL/MDRF												
Q. CHECK												
INT. STD.	B1110320 1/30											
SURROGATE	011602											
LPC												
ENDRIN/DDT												

Reviewed by,

C/B 5/18/02

Analyst LF
revision 4/12/01

Verified by my

Reviewed by DN
Date 2/26/02

Comment: INST 209; 30MX0.25umX0.50
 Operator: 209 GALOS
 Data Path: C:\HPCHEM\1\data\feb1102\
 Pre-Seq Cmd:
 Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch
 (X) Full Method (X) Inject Anyway
 () Reprocessing Only () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0208	
2	Sample	2	CAL80	GCMS0208	
3	Sample	3	LB	GCMS0208	GC/MS SCAN 2/8
4	Sample	4	REF	GCMS0208	GC/MS SCAN 2/8
5	Sample	5	REF1	GCMS0208	GC/MS SCAN 2/8
6	Sample	6	31432	GCMS0208	GC/MS SCAN 2/8
7	Sample	7	31433	GCMS0208	GC/MS SCAN 2/8
8	Sample	8	31500	GCMS0208	GC/MS SCAN 2/8
9	Sample	9	31501	GCMS0208	GC/MS SCAN 2/8
10	Sample	10	31502	GCMS0208	GC/MS SCAN 2/8
11	Sample	11	31503	GCMS0208	GC/MS SCAN 2/8
12	Sample	12	31504	GCMS0208	GC/MS SCAN 2/8
13	Sample	13	31505	GCMS0208	GC/MS SCAN 2/8 FB
14	Sample	14	31550	GCMS0208	GC/MS SCAN 2/8
15	Sample	15	31552	GCMS0208	GC/MS SCAN 2/8 FB
16	Sample	16	LBA	GCMS0208	GC/MS SCAN 1/15
17	Sample	17	REFA	GCMS0208	GC/MS SCAN 1/15
18	Sample	18	364	GCMS0208	GC/MS SCAN 1/15
19	Sample	19	945	GCMS0208	GC/MS SCAN 1/15
20	Sample	20	830	GCMS0208	GC/MS SCAN 1/15
21	Sample	21	831	GCMS0208	GC/MS SCAN 1/15
22	Sample	22	832	GCMS0208	GC/MS SCAN 1/15
23	Sample	23	833	GCMS0208	GC/MS SCAN 1/15
24	Sample	24	1036	GCMS0208	GC/MS SCAN 1/15
25	Sample	25	1037	GCMS0208	GC/MS SCAN 1/15
26	Sample	26	1038	GCMS0208	GC/MS SCAN 1/15
27	Sample	27	1039	GCMS0208	GC/MS SCAN 1/15
28	Sample	28	1040	GCMS0208	GC/MS SCAN 1/15
29	Sample	29	1042	GCMS0208	GC/MS SCAN 1/15
30	Sample	1	DFTPP2	GCMS0208	
31	Sample	2	CAL80A	GCMS0208	
32	Sample	30	LBB	GCMS0208	GC/MS SCAN 1/14
33	Sample	31	REFB	GCMS0208	GC/MS SCAN 1/14
34	Sample	32	898	GCMS0208	GC/MS SCAN 1/14
35	Sample	33	895	GCMS0208	GC/MS SCAN 1/14
36	Sample	34	899	GCMS0208	GC/MS SCAN 1/14
37	Sample	35	897	GCMS0208	GC/MS SCAN 1/14
38	Sample	36	896	GCMS0208	GC/MS SCAN 1/14
39	Sample	37	892	GCMS0208	GC/MS SCAN 1/14
40	Sample	38	915	GCMS0208	GC/MS SCAN 1/14 FB
41	Sample	39	918	GCMS0208	GC/MS SCAN 1/14
42	Sample	40	916	GCMS0208	GC/MS SCAN 1/14


Last Modified: Wed Feb 13 12:48:02 2002

Page: 1

DFTPP

Data File : C:\HPCHEM\1\DATA\FEB1102\DFTPP.D
Acq On : 11 Feb 2002 10:14 am
Sample :
Misc :
MS Integration Params: GOOFY.P
Method : C:\HPCHEM\1\METHODS\GCMS0208.M (RTE Integrator)
Title : 209 625/TCLP calibration table

Vial: 1
Operator: 209 GALOS
Inst : GC/MS 209
Multiplr: 1.00



Spectrum Information: Average of 26.921 to 26.976 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	39.6	104003	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	51.1	134419	PASS
70	69	0.00	2	0.1	181	PASS
127	198	40	60	47.5	124768	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	262811	PASS
199	198	5	9	6.6	17348	PASS
275	198	10	30	22.3	58681	PASS
365	198	1	100	2.6	6765	PASS
441	443	0.01	100	88.6	23524	PASS
442	198	40	110	52.7	138495	PASS
443	442	17	23	19.2	26557	PASS