

Method 508

Date Extracted 20 Nov 01

Page _____

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count
	LB	500ml		1ml			2ml			
	Ref ChL									
	Ref TOX									
	Ref Pest									
27955	S ChL				1ml					
27955	SD ChL									
27951	S TOX									
27951	SD TOX									
27957	S Pest									
27957	SD Pest									
	PRL 1248									
	27872									1
	27873									2
	27874									3
	27951									4
	27952									5
	27953									6
	27954									7
	27955									8
	27956									9
	27957									10
Remarks										
drip rate must be between 5-15 ml/min for liq-liq extractors										
drip rate -										
Extraction type - <u>Sep funnel</u> Solv. lot# <u>V43 E07</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V32603</u>										
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	<u>TOX: 2109</u>	<u>ChL</u>	<u>2108</u>	<u>Pest: 2062</u>						
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE	<u>18396</u>									
LPC										
ENDRIN/DDT										

Analyst AD/EA
revision 4/12/01

Verified by _____

Reviewed by DN
Date 11/7/2