

Method 508 MOD.

Date Extracted 3-11-08

Page 145

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
<u>L.B.</u>		<u>500ml</u>		<u>1ml</u>			<u>1ml</u>	<u>1B</u>	<u>7</u>			
<u>Ref 1</u>					<u>1ml</u>			<u>10B</u>				
<u>Ref 2</u>					<u>1ml</u>			<u>10A</u>				
<u>5303</u>								<u>9A</u>				
<u>5304</u>								<u>6A</u>				
<u>5305</u>								<u>8B</u>				
<u>5306</u>								<u>7B</u>				
<u>5307</u>								<u>2A</u>				
<u>5308</u>	<u>F.B.</u>							<u>4A</u>				
<u>5309</u>	<u>F.B.</u>							<u>8A</u>				
<u>5512</u>								<u>2B</u>				
<u>5513</u>	<u>Q</u>							<u>1A</u>				
<u>5514</u>								<u>7A</u>				
<u>5302</u>								<u>4B</u>				
Remarks <u>Na OH WAS ADDED TO THE Q-SAMPLE (5513) TO BRING pH TO 7.</u>												
Reviewed by, <u>C/B 3/11/08</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type <u>-SEP. Fun</u> Solv. lot# <u>V41470</u> NaCl lot# <u>V39618</u> Na2SO4 lot# <u>V46623</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	<u>2193</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	<u>21600</u>											
LPC												
ENDRIN/DDT												

Analyst A.D.
 Vision 4/12/01

Verified by

Reviewed by DW 3/29/08

Sequence Name: C:\MSDchem\1\sequence\031702.s
Comment:
Operator: McLean
Data Path: C:\MSDCHEM\1\DATA\031702\
Pre-Seq Cmd:
Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch
(X) Full Method (X) Inject Anyway
() Reprocessing Only () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFT0317	508SCAN	DFTPP
2	Sample	2	CAL40XA	508SCAN	40x
3	Sample	3	311REF1	508SCAN	3/11/02
4	Sample	4	311REF2	508SCAN	3/11/02
5	Sample	5	311LB	508SCAN	3/11/02
6	Sample	6	5303	508SCAN	3/11/02
7	Sample	7	5304	508SCAN	3/11/02
8	Sample	8	5305	508SCAN	3/11/02
9	Sample	9	5306	508SCAN	3/11/02
10	Sample	10	5307	508SCAN	3/11/02
11	Sample	11	5308FB	508SCAN	3/11/02
12	Sample	12	5309FB	508SCAN	3/11/02
13	Sample	13	5512	508SCAN	3/11/02
14	Sample	14	5513Q	508SCAN	3/11/02
15	Sample	15	5514	508SCAN	3/11/02
16	Sample	16	5302	508SCAN	3/11/02
17	Sample	1	DFTPP	508SCAN	
18	Sample	2	40XCAL	508SCAN	
19	Sample	17	314REF1	508SCAN	
20	Sample	18	314REF2	508SCAN	3/14/02
21	Sample	19	314LB	508SCAN	3/14/02
22	Sample	20	5910	508SCAN	3/14/02
23	Sample	21	5909	508SCAN	3/14/02
24	Sample	22	3914RE	508SCAN	3/14/02RE EXT
25	Sample	23	4027FBR	508SCAN	3/14/02RE EXT
26	Sample	24	3915RE	508SCAN	3/14/02RE EXT
27	Sample	25	5902FB	508SCAN	3/14/02
28	Sample	26	5911	508SCAN	3/14/02
29	Sample	27	5908	508SCAN	3/14/02
30	Sample	28	5912	508SCAN	3/14/02
31	Sample	29	5913	508SCAN	3/14/02

still to
process
✓ 3/25

DFTPP

Data File : C:\MSDCHEM\1\DATA\031702\DFT0317.D
 Acq On : 17 Mar 2002 8:48 am
 Sample : DFTPP
 Misc :
 MS Integration Params: BENZO.P

Vial: 1
 Operator: McLean
 Inst : Instrumen
 Multiplr: 1.00

Method : C:\MSDCHEM\1\5973N\5080302.M (RTE Integrator)
 Title : Quantitation For Method 508gcscan

Spectrum Information: Scan 2288

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	10	80	56.8	461120	PASS
68	69	0.00	2	1.7	6535	PASS
69	198	0.00	100	46.4	376640	PASS
70	69	0.00	2	0.5	1909	PASS
127	198	10	80	50.1	406784	PASS
197	198	0.00	2	0.0	0	PASS
198	198	50	100	100.0	811200	PASS
199	198	5	9	6.5	52832	PASS
275	198	10	60	22.3	181120	PASS
365	198	1	100	2.7	21952	PASS
441	443	0.01	100	78.3	95904	PASS
442	198	50	100	79.4	644352	PASS
443	442	15	24	19.0	122440	PASS

DFT0317.D 5080302.M Thu Mar 21 14:54:56 2002

File : C:\MSDCHEM\1\DATA\031702\DFT0317.D
 Operator : McLean
 Acquired : 17 Mar 2002 8:48 am using AcqMethod 508SCAN
 Instrument : Instrumen
 Sample Name: DFTPP
 Misc Info :
 Vial Number: 1

