

NYC-WTC_000060562

508SCAW

Method 525.2Instrument 213File Directory 040502 Storage Directory _____

Lab #	Filename	Description	Vial #		Inj. ul	Standards used - Lot #			Int. Std.	Remarks	Initials
			front	rear		Calibration	Reference	Surrogate			
	DFAA405	DFTPP	1		1	2110b					TS
	CSAA405	C5	2			2121d/2286d					
	C25A405	C25	52			1					
	SWRA405	Surrogate	3			2249					
	LBAD405	LB	4							ext 3/15/02	
	REFAD405	Ref 1	5								
	REFBO405	Ref 2	6								
	MCA405	MC	7								
	6049	Sample	8							ext 3/15/02	
	6050		9								
	6051		10								
	6052		11								
	6053		12								
	6054		13								
	6055		14								
	6056		15								
	6057		16								
	6089		17								
	DFTBO405	DFTPP	18			2110b					
	C10A405	C10	19			2121d/2286d					
	C10AA405	1	19			1					
	SWRBO405	Surrogate	20			2249					
	LBBO405	LB	21							ext 3/20/02	
	REFCO405	Ref 1	22								
	REFDO405	Ref 2	23								
	MCBO405	MC	24								
	3613RE	Sample	25							ext 3/20/02	
	6458	1	26								

Date 4/5/02 - 4/7/02Extraction date 3/15/02
3/20/02
3/22/02

Page _____

5085CAN

Method 625-2

Instrument 213

File Directory 040502 Storage Directory

Lab #	Filename	Description	Vial #		Inj.	Standards used - Lot #			Int. Std.	Remarks	Initials
			front	rear		Calibration	Reference	Surrogate			
	6459	sample	27		1					ext 3/20/02	TS
	6529		28								
	6530		29								
	6531		30								
	6581		31								
	6582		32								
	6583		33								
	6584		34								
	DFT0405	DFTPP	35			21106					
	C40A0405	C40	36			21212	22860				
	C40A0405	1	36			1					
	648ED405	surrogate	37			2249					
	LBC0405	LB	38							ext 3/22/02	
	REFED405	Ref 1	39								
	REFED405	Ref 2	40								
	MCC0405	MC	41								
	6800	sample	42							ext 3/24/02	
	6850		43								
	6851		44								
	6852		45								
	6853		46								
	6854		47								
	6855		48								
	6856		49								
	5908RE		50								
	3800RE		51								

Date 4/5-4/7/02

Extraction date 3/20/02
3/22/02

Page

Sequence Name: C:\MSDchem\1\sequence\040502.s

Comment:

Operator: Bohlk

Data Path: C:\MSDCHEM\1\DATA\040502\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch

(X) Full Method (X) Inject Anyway

() Reprocessing Only () Don't Inject

Line Type	Vial	DataFile	Method	Sample Name	
1 Sample	1	DFAAA0405	508SCAN	DFTPP	
2 Sample	2	C5AA0405	508SCAN	C5	
3 Sample	52	C25A0405	508SCAN	C25	
4 Sample	3	SURA0405	508SCAN	SURROGATE	
5 Sample	4	LBA0405	508SCAN	LB 3/15/02	
6 Sample	5	REFA0405	508SCAN	REF1 3/15/02	
7 Sample	6	REFB0405	508SCAN	REF2 3/15/02	
8 Sample	7	MCA0405	508SCAN	MC	
9 Sample	8	6049	508SCAN	SAMPLE 6049 3/15/02	
10 Sample	9	6050	508SCAN	SAMPLE 6050 3/15/02	
11 Sample	10	6051	508SCAN	SAMPLE 6051 3/15/02	
12 Sample	11	6052	508SCAN	SAMPLE 6052 3/15/02	
13 Sample	12	6053	508SCAN	SAMPLE 6053 3/15/02	
14 Sample	13	6054	508SCAN	SAMPLE 6054 3/15/02	
15 Sample	14	6055	508SCAN	SAMPLE 6055 3/15/02	
16 Sample	15	6056	508SCAN	SAMPLE 6056 3/15/02	
17 Sample	16	6057	508SCAN	SAMPLE 6057 3/15/02	
18 Sample	17	6089	508SCAN	SAMPLE 6089 3/15/02	
19 Sample	18	DFTB0405	508SCAN	DFTPP	
20 Sample	19	C10A0405	508SCAN	C10	
21 Sample	19	C10AA0405	508SCAN	C10	
22 Sample	20	SURB0405	508SCAN	SURROGATE	
23 Sample	21	LBB0405	508SCAN	LB 3/20/02	
24 Sample	22	REFC0405	508SCAN	REF1 3/20/02	
25 Sample	23	REFD0405	508SCAN	REF2 3/20/02	
26 Sample	24	MCB0405	508SCAN	MC	
27 Sample	25	3613RE	508SCAN	SAMPLE 3613 3/20/02	RE-EXT
28 Sample	26	6458	508SCAN	SAMPLE 6458 3/20/02	
29 Sample	27	6459	508SCAN	SAMPLE 6459 3/20/02	
30 Sample	28	6529	508SCAN	SAMPLE 6529 3/20/02	
31 Sample	29	6530	508SCAN	SAMPLE 6530 3/20/02	
32 Sample	30	6531	508SCAN	SAMPLE 6531 3/20/02	
33 Sample	31	6581	508SCAN	SAMPLE 6581 3/20/02	
34 Sample	32	6582	508SCAN	SAMPLE 6582 3/20/02	
35 Sample	33	6583	508SCAN	SAMPLE 6583 3/20/02	
36 Sample	34	6584	508SCAN	SAMPLE 6584 3/20/02	
37 Sample	35	DFTC0405	508SCAN	DFTPP	
38 Sample	36	C40A0405	508SCAN	C40	
39 Sample	36	C40AA0405	508SCAN	C40	
40 Sample	37	SURC0405	508SCAN	SURROGATE	
41 Sample	38	LBC0405	508SCAN	LB 3/22/02	
42 Sample	39	REFE0405	508SCAN	REF1 3/22/02	
43 Sample	40	REFF0405	508SCAN	REF2 3/22/02	

Sequence Name: C:\MSDchem\1\sequence\040502.s

Line	Type	Vial	DataFile	Method	Sample Name	
44	Sample	41	MCC0405	508SCAN	MC	
45	Sample	42	6800	508SCAN	SAMPLE 6800	3/22/02
46	Sample	43	6850	508SCAN	SAMPLE 6850	3/22/02
47	Sample	44	6851	508SCAN	SAMPLE 6851	3/22/02
48	Sample	45	6852	508SCAN	SAMPLE 6852	3/22/02
49	Sample	46	6853	508SCAN	SAMPLE 6853	3/22/02
50	Sample	47	6854	508SCAN	SAMPLE 6854	3/22/02
51	Sample	48	6855	508SCAN	SAMPLE 6855	3/22/02
52	Sample	49	6856	508SCAN	SAMPLE 6856	3/22/02
53	Sample	50	5908RE	508SCAN	SAMPLE 5908	3/22/02 RE-EXT
54	Sample	51	3800RE	508SCAN	SAMPLE 3800	3/22/02 RE-EXT

DFTPP

Data File : C:\MSDCHEM\1\DATA\040502\DFTC0405.D

Acq On : 6 Apr 2002 9:57 pm

Sample : DFTPP

Misc :

MS Integration Params: BENZO.P

Vial: 35

Operator: Bohlk

Inst : Instrumen

Multiplr: 1.00

Method : C:\MSDCHEM\1\5973N\5080402.M (RTE Integrator)

Title : Quantitation For Method 508gcscan

Spectrum Information: Scan 3518

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	44.1	320576	PASS
68	69	0.00	2	1.8	6164	PASS
69	198	0.00	100	47.7	347200	PASS
70	69	0.00	2	0.5	1563	PASS
127	198	40	60	51.0	370816	PASS
197	198	0.00	1	0.6	4347	PASS
198	198	100	100	100.0	727232	PASS
199	198	5	9	6.5	47176	PASS
275	198	10	30	21.2	154304	PASS
365	198	1	100	2.4	17632	PASS
441	443	0.01	100	82.1	65520	PASS
442	198	40	110	58.1	422592	PASS
443	442	17	23	18.9	79760	PASS

DFTC0405.D 5080402.M

Mon Apr 08 05:07:18 2002

File : C:\MSDCHEM\1\DATA\040502\DFTC0405.D
 Operator : Bohlk
 Acquired : 6 Apr 2002 9:57 pm using AcqMethod 508SCAN
 Instrument : Instrumen
 Sample Name: DFTPP
 Misc Info :
 Vial Number: 35

