

Method 508 MOD.

Date Extracted 3-25-02

Page 26

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
L.B.		500 mL		1 mL			1 mL	1A	7			
Ref 1					1 mL			4A				
Ref 2					1 mL			2B				
6580								9A				
6857	FB.							4A				
6957								2A				
6958								4B				
6959								9B				
6960								4A				
6961								10B				
6962								9B				
6963								9A				
6964	FB.							10A				
Remarks: <u>extract for FB</u> Sample 6857 contained water - <u>eluted</u> at three sodium sulfate. Reviewed by,												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type - <u>SEP. FUM.</u>		Solv. lot# <u>X01E19</u>		NaCl lot# <u>V39618</u>		Na2SO4 lot# <u>V46623</u>						
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	<u>2193</u>											
MDL/MORF												
STD CHECK												
INT STD.	<u>2286 m</u>											
SURROGATE	<u>4249</u>											
PLC												
ENDRIN/ODT												

Analyst A.D.
Date 4/12/01

Verified by _____

Reviewed by DN
Date 4/23/02

Method 525.2

Instrument 213

File Directory 041102

Storage Directory

Lab #	Filename	Description	Vial #		Inj. ul	Standards used - Lot #			Int. Std.	Remarks	Initials
			front	rear		Calibration	Reference	Surrogate			
	DFAA0417	DFA0P	1		1	21106					TB
	C10A0417	C10	2			020412				based on	
	C80A0417	C80	3			6				100%	
	SURAD0417	Surrogate	4			2249				4/12/02	
	LBAD0417	LB	5						2286m		
	6857	FB	6								
	REFAD0417	Rep 1	7								
	REFBD0417	Rep 2	8								
	MCA0417	MC	9								
	6580	Sample	10						2286m		
	6957		11								
	6958		12								
	6959		13								
	6960		14								
	6961		15								
	6962		16								
	6963		17								
	6964		18								

Date

4/17/02

Extraction date

3/25/02

Page

Sequence Name: C:\MSDCHEM\1\SEQUENCE\041702.S

Comment:

Operator: Bohlk

Data Path: C:\MSDCHEM\1\DATA\041702\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch

☒ (X) Full Method ☒ (X) Inject Anyway☐ () Reprocessing Only ☐ () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTAA417	508SCAN	DFTPP
2	Sample	2	C10A0417	508SCAN	C10
3	Sample	3	C80A0417	508SCAN	C80
4	Sample	4	SURA0417	508SCAN	SURROGATE, DDD
5	Sample	5	LBA0417	508SCAN	LB 3/25/02
6	Sample	6	6857	508SCAN	SAMPLE 6857, FB 3/25/02
7	Sample	7	REFA0417	508SCAN	REF1 3/25/02
8	Sample	8	REFB0417	508SCAN	REF2 3/25/02
9	Sample	9	MCA0417	508SCAN	MC
10	Sample	10	6580	508SCAN	SAMPLE 6580 3/25/02
11	Sample	11	6957	508SCAN	SAMPLE 6957 3/25/02
12	Sample	12	6958	508SCAN	SAMPLE 6958 3/25/02
13	Sample	13	6959	508SCAN	SAMPLE 6959 3/25/02
14	Sample	14	6960	508SCAN	SAMPLE 6960 3/25/02
15	Sample	15	6961	508SCAN	SAMPLE 6961 3/25/02
16	Sample	16	6962	508SCAN	SAMPLE 6962 3/25/02
17	Sample	17	6963	508SCAN	SAMPLE 6963 3/25/02
18	Sample	18	6964	508SCAN	SAMPLE 6964, FB 3/25/02