

Sequence Name: M:\INSTRU~1\104\SEQUENCE\515\010915.S

Comment: 515 SAMPLES

Operator: BOHLK

Data Path: C:\HPCHEM\1\DATA\515F\010904H\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

(X) Full Method

(X) Inject Anyway

() Reprocessing Only

() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	2nd Info	1	-----	515AF	
2	Sample	2	HEXANE1	515AF	hex
3	2nd Info	3	-----	515AF	
4	Sample	4	0915LPC	515AF	LPC
5	2nd Info	5	-----	515AF	
6	Sample	6	0915C1B	515AF	C1B
7	2nd Info	7	-----	515AF	
8	Sample	8	0915C1	515AF	C1
9	2nd Info	9	-----	515AF	
10	Sample	10	0915C2B	515AF	C2B
11	2nd Info	11	-----	515AF	
12	Sample	12	0915C2	515AF	C2
13	2nd Info	13	-----	515AF	
14	Sample	14	0915C3	515AF	C3
15	2nd Info	15	-----	515AF	
16	Sample	16	0915UNR	515AF	UNEXTRACTED 515 REF
17	2nd Info	17	-----	515AF	
18	Sample	18	0915UNX	515AF	UNEXTRACTED DCPA REF
19	2nd Info	19	-----	515AF	
20	Sample	20	0910LB	515AF	09/10 lab blk
21	2nd Info	21	-----	515AF	
22	Sample	22	0910MDL	515AF	09/10 MDL
23	2nd Info	23	-----	515AF	
24	Sample	24	0910REF	515AF	09/10 REF STD, EXTD
25	2nd Info	1	-----	515AF	
26	Sample	2	HEXANE2	515AF	HEXANE
27	2nd Info	25	-----	515AF	
28	Sample	26	21394	515AF	SAMPLE 221394
29	2nd Info	27	-----	515AF	
30	Sample	28	21395	515AF	SAMPLE 21395
31	2nd Info	29	-----	515AF	
32	Sample	30	21395S	515AF	SAMPLE 21395 spike
33	2nd Info	31	-----	515AF	
34	Sample	32	21395SD	515AF	SAMPLE 21395 spike dupe
35	2nd Info	1	-----	515AF	
36	Sample	2	HEXANE3	515AF	HEXANE
37	2nd Info	33	-----	515AF	
38	Sample	34	0912LB	515AF	09/12 lab blk
39	2nd Info	35	-----	515AF	
40	Sample	36	0912REF	515AF	09/12 ref
41	2nd Info	1	-----	515AF	
42	Sample	2	HEXANE4	515AF	HEXANE

Last Modified: Sat Sep 15 14:15:54 2001

Page: 1

43 2nd Info 37 ----- 515AF

Sequence Name: M:\INSTRU~1\104\SEQUENCE\515\010915.S

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	38	22166	515AF	sample 22166
45	2nd Info	39	-----	515AF	
46	Sample	40	22167	515AF	sample 22167
47	2nd Info	41	-----	515AF	
48	Sample	42	22168	515AF	sample 22168
49	2nd Info	43	-----	515AF	
50	Sample	44	22169	515AF	sample 22169
51	2nd Info	45	-----	515AF	
52	Sample	46	22170	515AF	sample 22170
53	2nd Info	47	-----	515AF	
54	Sample	48	22171	515AF	sample 22166
55	2nd Info	1	-----	515AF	
56	Sample	2	HEXANE5	515AF	hexane
57	2nd Info	49	-----	515AF	
58	Sample	50	C2	515AF	c2