

Method 508

Date Extracted 9/18/01

Page 86

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH		
	Lab Blank	500		✓			2ml		7		next
70X	MDL				✓						
70X	REF				✓						
1248	PRL Check				✓						
	22548									10	
	22549									11	
	22550									12	
	22551									13	
	22552									14	
	22553									15	
	22558									16	✓
	22577	✓		✓			✓		✓	17	
QC Reviewed by: <u>11/13/01</u>											
Remarks											
drip rate must be between 5-15 ml/min for liq-liq extractors Extraction type - <u>Sept.</u> Solv. lot# <u>9262-03</u> NaCl lot# <u>T72617</u> Na2SO4 lot# <u>U18617</u> HEX <u>92614-03</u>											
STANDARD	INT STD LOG#										
CALIBRATION											
REF/SPIKE											
MDL/MDRF											
Q. CHECK											
INT. STD.											
SURROGATE											
LPC											
ENDRIN/DDT											

Analyst MMF/AD
revision 10/15/99

Verified by _____

Supv. approved DN
9/11

