

Method 508

Date Extracted 9/21/01

Page

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH = 7	Sp count
	Lab 11h	500ml	—	✓		—	2ml			
	Ref Post									
	MDL Post									
	REL 1016									
22976	spike post				✓					
22976	spike dup post				✓					
22977	spike post				✓					
22977	spike dup post				✓					
	22976									1
	22977									2
	22978									3
	22979									4
	23079									5
	23080									6
	23081									7
	23082									8
	23083									9
	23084									10
	23085									11
	23123									12
Remarks										
QC Reviewed by <u>[Signature]</u> 9/21/01										
drip rate must be between 5-15 ml/min for liq-liq extractors			drip rate -							
Extraction type <u>SEP</u>		Solv. lot#	NaCl lot#		Na2SO4 lot#					
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	<u>post 2061W</u>	<u>1021078W</u>								
MDL/MDRF	<u>2061W</u>	<u>1678W</u>								
Q. CHECK										
INT. STD.										
SURROGATE	<u>20640</u>									
LPC										
ENDRIN/DDT										

Analyst my AD/CM.
revision 10/15/99

Verified by

Supv. approved [Signature] 10/15/01