

Method


508

Date Extracted

9/27/01

Page

6

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
	leukthly	Stone		✓			2ml					
	pest chloroform											
	MDL chloro											
	PRC 1248											
23473	spike amp pest				✓							
23473	spike pest				✓							
23497	spike amp pest				✓							
23497	spike pest				✓							
	23498									6		
	23499									7		
	23500									8		
	23501									9		
	23502									10		
	23503									11		
	23605									12		
	23606									13		
	23607									14		
	23650	✓	✓				✓		✓	15		
	23473	✓	✓				✓		✓	16		
	23497	✓	✓				✓		✓	17		
Remarks												
QC Reviewed by,  12/13/01												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type - <u>Sep t.</u> Solv. lot# NaCl lot# <u>T12617</u> Na2SO4 lot# <u>V18617</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE												
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE												
LPC												
ENDRIN/DDT												

Analyst


CM/AD

revision 4/12/01

Verified by

Reviewed by

Date


 10/17/01