

Sequence-Name: C:\HPCHEM\1\SEQUENCE\508F\011012.S

Comment: samples 508

Operator: Nefliu

Data Path: C:\HPCHEM\1\DATA\508F\011012\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch
 (X) Full Method (X) Inject Anyway
 () Reprocessing Only () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	2nd Info	1	-----	508F	
2	Sample	2	HEXANE1	508F	hexane
3	2nd Info	3	-----	508F	
4	Sample	4	1012END	508F	endrin decomposition
5	2nd Info	5	-----	508F	
6	Sample	6	1012LPC	508F	lpc
7	2nd Info	7	-----	508F	
8	Sample	8	1248C2	508F	1248 c2
9	2nd Info	9	-----	508F	
10	Sample	10	CHLC1	508F	chlordan ccv1
11	2nd Info	11	-----	508F	
12	Sample	12	TOXC1	508F	toxaphene ic1
13	2nd Info	13	-----	508F	
14	Sample	14	TOXC2	508F	toxaphene ic2
15	2nd Info	15	-----	508F	
16	Sample	16	TOXC3	508F	toxaphene ic3
17	2nd Info	17	-----	508F	
18	Sample	18	TOXUR	508F	toxaphene unextracted ref
19	2nd Info	19	-----	508F	
20	Sample	20	PESTC1	508F	pest ic1
21	2nd Info	21	-----	508F	
22	Sample	22	PESTC2	508F	pest ic2
23	2nd Info	23	-----	508F	
24	Sample	24	PESTC3	508F	pest ic3
25	2nd Info	25	-----	508F	
26	Sample	26	PESTUR	508F	pest unextracted ref
27	2nd Info	85	-----	508F	
28	Sample	86	2END	508F	endrin decomposition
29	2nd Info	27	-----	508F	
30	Sample	28	0927LB	508F	09/27 lab blk
31	2nd Info	29	-----	508F	
32	Sample	30	0927MDL	508F	09/27 mdl chlordan
33	2nd Info	31	-----	508F	
34	Sample	32	0927REF	508F	09/27 ref chlordan
35	2nd Info	33	-----	508F	
36	Sample	34	PRL1248	508F	prl 1248
37	2nd Info	35	-----	508F	
38	Sample	36	23473	508F	sample
39	2nd Info	37	-----	508F	
40	Sample	38	23497	508F	sample
41	2nd Info	39	-----	508F	
42	Sample	40	23498	508F	sample
43	2nd Info	41	-----	508F	

Sequence Name: C:\HPCHEM\1\SEQUENCE\508F\011012:S

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	42	23499	508F	sample
45	2nd Info	43	-----	508F	
46	Sample	44	23500	508F	sample
47	2nd Info	45	-----	508F	
48	Sample	46	23501	508F	sample
49	2nd Info	47	-----	508F	
50	Sample	48	23502	508F	sample
51	2nd Info	49	-----	508F	
52	Sample	50	23503	508F	sample
53	2nd Info	51	-----	508F	
54	Sample	52	23605	508F	sample
55	2nd Info	53	-----	508F	
56	Sample	54	23606	508F	sample
57	2nd Info	55	-----	508F	
58	Sample	56	23607	508F	sample
59	2nd Info	57	-----	508F	
60	Sample	58	23650	508F	sample
61	2nd Info	59	-----	508F	
62	Sample	60	CHLC3	508F	chlordan cccv2
63	2nd Info	61	-----	508F	
64	Sample	62	23497ST	508F	23497 spike tox
65	2nd Info	63	-----	508F	
66	Sample	64	23497STD	508F	23497 spike tox dupe
67	2nd Info	65	-----	508F	
68	Sample	66	2TOXC2	508F	tox ccv
69	2nd Info	67	-----	508F	
70	Sample	68	23473SP	508F	23473 spike pest
71	2nd Info	69	-----	508F	
72	Sample	70	23473SPD	508F	23473 spike pest dupe
73	2nd Info	71	-----	508F	
74	Sample	72	2PESTC2	508F	pest ccv
75	2nd Info	73	-----	508F	
76	Sample	74	1016C2	508F	1016 c2
77	2nd Info	75	-----	508F	
78	Sample	76	1221C2	508F	1221 c2
79	2nd Info	77	-----	508F	
80	Sample	78	1232C2	508F	1232 c2
81	2nd Info	79	-----	508F	
82	Sample	80	1242C2	508F	1242 c2
83	2nd Info	81	-----	508F	
84	Sample	82	1254C2	508F	1254 c2
85	2nd Info	83	-----	508F	
86	Sample	84	1260C2	508F	1260 c2