

batches 16390-91-16402

Method 200.8 - Summary Report**Sample ID: Blank**

Sample Date/Time: Friday, January 11, 2002 12:31:45

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth


Dataset File: d:\elandata\Dataset\daily performance\Blank.2016

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

**Concentration Results**

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.333	18.368				ug/L
>	Sc-1	45	634943.467	4.280				ug/L
[	Mn	55	626.028	5.124				ug/L
>	In	115	1318775.684	1.001				ug/L
	Sb	121	35.334	10.715				ug/L
[	Sb	123	26.516	9.471				ug/L
>	Ho	165	1796034.933	2.155				ug/L
[	Tl	205	216.671	5.545				ug/L

Reviewed by,

 3/26/02

Sample ID: Blank

Report Date/Time: Friday, January 11, 2002 12:32:18

Page 1

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45					
[	Mn	55					
[>	In	115					
[	Sb	121					
[	Sb	123					
[>	Ho	165					
[	Tl	205					

Sample ID: Blank

Report Date/Time: Friday, January 11, 2002 12:32:18

Page 2

Method 200.8 - Summary Report

Sample ID: Standard 1

Sample Date/Time: Friday, January 11, 2002 12:34:44

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\Standard 1.2017

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14805.195	1.384	0.024	30.00000	3.97	ug/L
>	Sc-1	45	629139.254	4.226	629139.254			ug/L
[	Mn	55	538181.772	1.552	0.855	30.00000	3.30	ug/L
[>	In	115	1310300.195	1.370	1310300.195			ug/L
	Sb	121	254292.433	1.369	0.194	30.00000	1.96	ug/L
[	Sb	123	194033.326	1.557	0.148	30.00000	0.74	ug/L
[>	Ho	165	1797239.228	2.146	1797239.228			ug/L
[	Tl	205	845580.927	2.652	0.470	30.00000	1.48	ug/L

Sample ID: Standard 1

Report Date/Time: Friday, January 11, 2002 12:35:17

Page 1

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
] >	Sc-1	45					
[	Mn	55					
] >	In	115					
[	Sb	121					
[	Sb	123					
] >	Ho	165					
[	Tl	205					

Sample ID: Standard 1

Report Date/Time: Friday, January 11, 2002 12:35:17

Page 2

Method 200.8 - Summary Report

Sample ID: Standard 2

Sample Date/Time: Friday, January 11, 2002 12:36:49

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\Standard 2.2018

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	29638.868	0.753	0.048	60.11729	2.44	ug/L
[>	Sc-1	45	623391.727	2.152	623391.727			ug/L
[	Mn	55	1082436.790	1.457	1.736	60.17630	2.17	ug/L
[>	In	115	1325048.942	0.803	1325048.942			ug/L
[	Sb	121	511784.959	1.466	0.386	59.94161	2.28	ug/L
[	Sb	123	393353.657	1.447	0.297	60.02907	1.39	ug/L
[>	Ho	165	1778879.442	3.503	1778879.442			ug/L
[	Tl	205	1645173.875	2.308	0.925	59.79728	2.01	ug/L

Sample ID: Standard 2

Report Date/Time: Friday, January 11, 2002 12:37:21

Page 1

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45					
[	Mn	55					
[>	In	115					
	Sb	121					
[	Sb	123					
[>	Ho	165					
[	Tl	205					

Sample ID: Standard 2

Report Date/Time: Friday, January 11, 2002 12:37:21

Page 2

Method 200.8 - Summary Report

Sample ID: ICV CAL BK

Sample Date/Time: Friday, January 11, 2002 12:40:24

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\ICV CAL BK.2019

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	13.333	22.913	0.000	0.00447	140.93	ug/L
[>	Sc-1	45	624229.826	1.694	624229.826			ug/L
[	Mn	55	664.698	4.108	0.000	0.00274	59.11	ug/L
[>	In	115	1304267.467	2.297	1304267.467			ug/L
[	Sb	121	309.008	7.053	0.000	0.03262	8.39	ug/L
[	Sb	123	232.521	1.881	0.000	0.03201	4.68	ug/L
[>	Ho	165	1762499.547	0.959	1762499.547			ug/L
[	Tl	205	356.677	3.506	0.000	0.00528	8.06	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		98.313			
[	Mn	55					
[>	In	115		98.900			
[	Sb	121					
[	Sb	123					
[>	Ho	165		98.133			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: ICV CAL BK

Report Date/Time: Friday, January 11, 2002 12:40:57

Page 1

Method 200.8 - Summary Report**Sample ID: ICV STD1 30 PPB**

Sample Date/Time: Friday, January 11, 2002 12:42:30

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\ICV STD1 30 PPB.2020

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14597.859	1.793	0.023	29.37907	1.72	ug/L
>	Sc-1	45	628003.738	3.010	628003.738			ug/L
[	Mn	55	536322.932	1.185	0.853	29.58427	2.04	ug/L
>	In	115	1292264.235	1.204	1292264.235			ug/L
	Sb	121	250358.476	1.904	0.194	30.06790	3.04	ug/L
[	Sb	123	193040.943	0.529	0.149	30.20873	1.73	ug/L
>	Ho	165	1754612.620	0.937	1754612.620			ug/L
[	Tl	205	845628.678	2.907	0.482	31.14386	2.32	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		98.907			
[	Mn	55					
>	In	115		97.990			
	Sb	121					
[	Sb	123					
>	Ho	165		97.694			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: ICV STD1 30 PPB

Report Date/Time: Friday, January 11, 2002 12:43:02

Page 1

Method 200.8 - Summary Report

Sample ID: ICV STD2 60 PPB

Sample Date/Time: Friday, January 11, 2002 12:44:34

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\ICV STD2 60 PPB.2021

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	29011.574	4.643	0.046	58.29368	3.90	ug/L
[>	Sc-1	45	629087.683	2.602	629087.683			ug/L
[	Mn	55	1058967.550	2.266	1.683	58.33427	1.94	ug/L
[>	In	115	1309371.244	1.029	1309371.244			ug/L
[	Sb	121	500269.411	1.078	0.382	59.29144	1.50	ug/L
[	Sb	123	388564.474	1.277	0.297	60.01432	2.04	ug/L
[>	Ho	165	1770388.046	1.027	1770388.046			ug/L
[	Tl	205	1632140.770	3.756	0.922	59.58820	3.74	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		99.078			
[	Mn	55					
[>	In	115		99.287			
[	Sb	121					
[	Sb	123					
[>	Ho	165		98.572			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: ICV STD2 60 PPB

Report Date/Time: Friday, January 11, 2002 12:45:06

Page 1

Method 200.8 - Summary Report

Sample ID: QCS 50 PPB

Sample Date/Time: Friday, January 11, 2002 12:46:42

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\QCS 50 PPB.2022

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	24993.562	0.956	0.039	49.84527	1.78	ug/L
>	Sc-1	45	633931.052	2.427	633931.052			ug/L
[	Mn	55	912005.018	2.410	1.438	49.86739	3.96	ug/L
>	In	115	1322601.056	0.379	1322601.056			ug/L
	Sb	121	437663.670	0.856	0.331	51.34776	0.49	ug/L
[	Sb	123	335725.092	2.369	0.254	51.32553	2.03	ug/L
>	Ho	165	1795447.295	3.374	1795447.295			ug/L
[	Tl	205	1408723.964	4.261	0.784	50.71003	2.47	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		99.841			
[	Mn	55					
>	In	115		100.290			
	Sb	121					
[	Sb	123					
>	Ho	165		99.967			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: QCS 50 PPB

Report Date/Time: Friday, January 11, 2002 12:47:15

Page 1

Method 200.8 - Summary Report

Sample ID: LFB 50 PPB

Sample Date/Time: Friday, January 11, 2002 12:48:45

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\LFB 50 PPB.2023

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	24149.954	2.363	0.038	48.33262	3.03	ug/L
[>	Sc-1	45	631658.292	1.682	631658.292			ug/L
[	Mn	55	887138.513	2.066	1.404	48.65866	1.87	ug/L
[>	In	115	1291137.928	0.843	1291137.928			ug/L
[	Sb	121	417638.140	1.011	0.323	50.19610	1.51	ug/L
[	Sb	123	325186.052	0.814	0.252	50.92871	0.59	ug/L
[>	Ho	165	1759072.604	1.840	1759072.604			ug/L
[	Tl	205	1356864.518	1.839	0.771	49.85534	0.86	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		99.483			
[	Mn	55					
[>	In	115		97.904			
[	Sb	121					
[	Sb	123					
[>	Ho	165		97.942			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: LFB 50 PPB

Report Date/Time: Friday, January 11, 2002 12:49:18

Page 1

Method 200.8 - Summary Report

Sample ID: LRB

Sample Date/Time: Friday, January 11, 2002 12:51:12

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\LRB.2024

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.000	9.091	-0.000	-0.00007	2652.54	ug/L
[>	Sc-1	45	618246.301	2.131	618246.301			ug/L
[	Mn	55	851.717	23.655	0.000	0.01366	85.82	ug/L
[>	In	115	1285385.390	1.338	1285385.390			ug/L
[	Sb	121	247.339	5.808	0.000	0.02571	7.34	ug/L
[	Sb	123	200.314	7.996	0.000	0.02746	10.05	ug/L
[>	Ho	165	1730347.275	1.874	1730347.275			ug/L
[	Tl	205	515.020	10.904	0.000	0.01145	18.60	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		97.370			
[	Mn	55					
[>	In	115		97.468			
[	Sb	121					
[	Sb	123					
[>	Ho	165		96.343			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: LRB

Report Date/Time: Friday, January 11, 2002 12:51:45

Page 1

Method 200.8 - Summary Report

Sample ID: AD31456

Sample Date/Time: Friday, January 11, 2002 12:53:22

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31456.2025

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.000	30.000	-0.000	-0.00285	193.31	ug/L
>	Sc-1	45	639338.937	2.933	639338.937			ug/L
[	Mn	55	1581530.953	1.582	2.473	85.74657	1.41	ug/L
[>	In	115	1323670.557	1.115	1323670.557			ug/L
	Sb	121	712.368	3.278	0.001	0.07935	2.40	ug/L
[	Sb	123	551.575	6.599	0.000	0.08018	6.20	ug/L
[>	Ho	165	1811548.525	1.555	1811548.525			ug/L
[	Tl	205	356.677	0.648	0.000	0.00493	3.90	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		100.692			
[	Mn	55					
[>	In	115		100.371			
	Sb	121					
[	Sb	123					
[>	Ho	165		100.864			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31456

Report Date/Time: Friday, January 11, 2002 12:53:54

Page 1

Method 200.8 - Summary Report

Sample ID: AD31456

Sample Date/Time: Friday, January 11, 2002 12:55:49

Sample Type: QC Duplicate

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31456.2026

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL): 5e+001

Diluted To Volume (mL): 1e+002

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	13.000	7.692	0.000	0.00610	44.31	ug/L
[>	Sc-1	45	641069.108	2.537	641069.108			ug/L
[	Mn	55	867469.024	3.467	1.352	93.76459	2.98	ug/L
[>	In	115	1326488.544	1.654	1326488.544			ug/L
[	Sb	121	390.012	3.353	0.000	0.08294	3.51	ug/L
[	Sb	123	310.973	3.685	0.000	0.08666	2.26	ug/L
[>	Ho	165	1816750.613	1.696	1816750.613			ug/L
[	Tl	205	278.674	5.229	0.000	0.00424	26.89	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					551.674
[>	Sc-1	45		100.965			
[	Mn	55					✓ 8.933
[>	In	115		100.585			
[	Sb	121					4.433
[	Sb	123					7.771
[>	Ho	165		101.153			
[	Tl	205					15.071

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31456

Report Date/Time: Friday, January 11, 2002 12:56:21

Page 1

Method 200.8 - Summary Report

Sample ID: AD31456

Sample Date/Time: Friday, January 11, 2002 12:58:15

Sample Type: Spike - 2

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31456.2027

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	24583.791	2.747	0.039	49.56342	2.44	ug/L
[>	Sc-1	45	627069.761	3.011	627069.761			ug/L
[	Mn	55	2149294.564	1.617	3.429	118.86127	3.39	ug/L
[>	In	115	1317500.528	1.913	1317500.528			ug/L
[	Sb	121	432671.863	1.039	0.328	50.97222	2.25	ug/L
[	Sb	123	331602.005	1.240	0.252	50.90131	1.54	ug/L
[>	Ho	165	1788221.263	1.780	1788221.263			ug/L
[	Tl	205	1282097.666	1.840	0.717	46.33863	0.51	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9			99.133		
[>	Sc-1	45		98.760			
[	Mn	55			66.229	70-130	
[>	In	115		99.903			
[	Sb	121			101.786		
[	Sb	123			101.642		
[>	Ho	165		99.565			
[	Tl	205			92.667		

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31456

Report Date/Time: Friday, January 11, 2002 12:58:48

Page 1

Method 200.8 - Summary Report

Sample ID: AD31456

Sample Date/Time: Friday, January 11, 2002 13:01:05

Sample Type: Spike - 2

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31456.2028

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	23978.497	2.094	0.038	48.03863	4.72	ug/L
[>	Sc-1	45	631502.406	3.496	631502.406			ug/L
[	Mn	55	2143030.380	0.953	3.395	117.68436	2.63	ug/L
[>	In	115	1283355.998	1.307	1283355.998			ug/L
[	Sb	121	426202.587	0.962	0.332	51.54004	1.86	ug/L
[	Sb	123	331909.008	0.688	0.259	52.30318	1.66	ug/L
[>	Ho	165	1767965.648	0.521	1767965.648			ug/L
[	Tl	205	1284073.024	3.191	0.726	46.93812	2.75	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9			96.083		
[>	Sc-1	45		99.458			
[	Mn	55			63.876		
[>	In	115		97.314			
[	Sb	121			102.921		
[	Sb	123			104.446		
[>	Ho	165		98.437			
[	Tl	205			93.866		

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31456

Report Date/Time: Friday, January 11, 2002 13:01:38

Page 1

Method 200.8 - Summary Report

Sample ID: AD31457

Sample Date/Time: Friday, January 11, 2002 13:03:51

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31457.2029

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	12.667	44.891	0.000	0.00326	372.26	ug/L
[>	Sc-1	45	622954.194	4.636	622954.194			ug/L
[	Mn	55	556469.844	3.538	0.893	30.95786	3.89	ug/L
[>	In	115	1318514.786	0.409	1318514.786			ug/L
[	Sb	121	826.379	2.844	0.001	0.09311	3.31	ug/L
[	Sb	123	633.878	7.103	0.000	0.09316	7.43	ug/L
[>	Ho	165	1793513.106	2.332	1793513.106			ug/L
[	Tl	205	419.014	0.955	0.000	0.00731	6.07	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		98.112			
[	Mn	55					
[>	In	115		99.980			
[	Sb	121					
[	Sb	123					
[>	Ho	165		99.860			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31457

Report Date/Time: Friday, January 11, 2002 13:04:24

Page 1

Method 200.8 - Summary Report

Sample ID: AD31458

Sample Date/Time: Friday, January 11, 2002 13:05:55

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31458.2030

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	12.667	44.891	0.000	0.00325	373.20	ug/L
[>	Sc-1	45	623371.579	3.352	623371.579			ug/L
[	Mn	55	571800.160	2.926	0.917	31.77653	2.81	ug/L
[>	In	115	1321710.386	1.613	1321710.386			ug/L
[	Sb	121	638.029	4.934	0.000	0.07075	4.83	ug/L
[	Sb	123	483.946	1.500	0.000	0.07000	3.05	ug/L
[>	Ho	165	1777002.096	2.935	1777002.096			ug/L
[	Tl	205	358.344	8.756	0.000	0.00525	25.51	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		98.177			
[	Mn	55					
[>	In	115		100.223			
[	Sb	121					
[	Sb	123					
[>	Ho	165		98.940			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31458

Report Date/Time: Friday, January 11, 2002 13:06:27

Page 1

Method 200.8 - Summary Report

Sample ID: AD31460

Sample Date/Time: Friday, January 11, 2002 13:08:37

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31460.2031

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.333	51.701	0.000	0.00061	2090.51	ug/L
[>	Sc-1	45	624799.477	2.979	624799.477			ug/L
[	Mn	55	1623164.640	0.099	2.598	90.08186	2.94	ug/L
[>	In	115	1274541.806	0.340	1274541.806			ug/L
[	Sb	121	682.366	3.162	0.001	0.07892	3.00	ug/L
[	Sb	123	523.876	4.041	0.000	0.07905	4.00	ug/L
[>	Ho	165	1747210.007	0.628	1747210.007			ug/L
[	Tl	205	346.343	10.045	0.000	0.00502	27.12	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		98.402			
[	Mn	55					
[>	In	115		96.646			
[	Sb	121					
[	Sb	123					
[>	Ho	165		97.282			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31460

Report Date/Time: Friday, January 11, 2002 13:09:10

Page 1

Method 200.8 - Summary Report

Sample ID: AD31461

Sample Date/Time: Friday, January 11, 2002 13:11:01

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31461.2032

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	15.333	9.962	0.000	0.00909	24.58	ug/L
[>	Sc-1	45	611621.307	2.867	611621.307			ug/L
[	Mn	55	555921.791	2.537	0.908	31.48699	2.87	ug/L
[>	In	115	1265884.688	1.668	1265884.688			ug/L
[	Sb	121	474.350	5.499	0.000	0.05397	4.12	ug/L
[	Sb	123	374.437	3.894	0.000	0.05577	4.92	ug/L
[>	Ho	165	1732258.570	2.817	1732258.570			ug/L
[	Tl	205	354.677	5.649	0.000	0.00544	12.88	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		96.327			
[	Mn	55					
[>	In	115		95.989			
[	Sb	121					
[	Sb	123					
[>	Ho	165		96.449			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31461

Report Date/Time: Friday, January 11, 2002 13:11:34

Page 1

Method 200.8 - Summary Report

Sample ID: AD31462

Sample Date/Time: Friday, January 11, 2002 13:13:05

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31462.2033

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	20.145	-0.000	-0.00150	279.04	ug/L
[>	Sc-1	45	620153.458	1.773	620153.458			ug/L
[	Mn	55	562478.364	0.260	0.906	31.41517	1.52	ug/L
[>	In	115	1295023.483	1.610	1295023.483			ug/L
[	Sb	121	470.684	3.862	0.000	0.05223	2.59	ug/L
[	Sb	123	364.733	2.900	0.000	0.05289	2.78	ug/L
[>	Ho	165	1770990.687	0.426	1770990.687			ug/L
[	Tl	205	314.009	3.978	0.000	0.00366	13.05	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		97.671			
[	Mn	55					
[>	In	115		98.199			
[	Sb	121					
[	Sb	123					
[>	Ho	165		98.606			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31462

Report Date/Time: Friday, January 11, 2002 13:13:37

Page 1

Method 200.8 - Summary Report

Sample ID: AD31464

Sample Date/Time: Friday, January 11, 2002 13:15:18

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31464.2034

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	13.667	15.232	0.000	0.00525	76.88	ug/L
>	Sc-1	45	620802.390	1.872	620802.390			ug/L
[	Mn	55	1549138.993	2.234	2.495	86.47995	1.73	ug/L
>	In	115	1286724.147	1.758	1286724.147			ug/L
	Sb	121	660.697	2.667	0.000	0.07554	3.38	ug/L
[	Sb	123	511.324	4.074	0.000	0.07627	2.59	ug/L
>	Ho	165	1756930.999	0.726	1756930.999			ug/L
[	Tl	205	318.009	6.897	0.000	0.00391	22.90	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		97.773			
[	Mn	55					
>	In	115		97.570			
	Sb	121					
[	Sb	123					
>	Ho	165		97.823			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31464

Report Date/Time: Friday, January 11, 2002 13:15:51

Page 1

Method 200.8 - Summary Report

Sample ID: AD31465

Sample Date/Time: Friday, January 11, 2002 13:18:12

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31465.2035

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.333	18.368	0.000	0.00166	327.17	ug/L
[>	Sc-1	45	595675.652	4.902	595675.652			ug/L
[	Mn	55	550689.152	1.787	0.925	32.06312	4.90	ug/L
[>	In	115	1240941.770	0.300	1240941.770			ug/L
[	Sb	121	439.015	7.772	0.000	0.05075	8.59	ug/L
[	Sb	123	343.811	3.100	0.000	0.05196	3.17	ug/L
[>	Ho	165	1690168.737	1.788	1690168.737			ug/L
[	Tl	205	284.341	1.332	0.000	0.00308	10.80	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		93.816			
[	Mn	55					
[>	In	115		94.098			
[	Sb	121					
[	Sb	123					
[>	Ho	165		94.106			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31465

Report Date/Time: Friday, January 11, 2002 13:18:44

Page 1

Method 200.8 - Summary Report

Sample ID: AD31466

Sample Date/Time: Friday, January 11, 2002 13:20:22

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31466.2036

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.667	21.571	0.000	0.00150	362.87	ug/L
[>	Sc-1	45	613787.940	4.085	613787.940			ug/L
[	Mn	55	584822.558	2.382	0.953	33.03505	4.84	ug/L
[>	In	115	1274639.732	0.714	1274639.732			ug/L
[	Sb	121	437.015	5.839	0.000	0.04906	7.04	ug/L
[	Sb	123	350.102	3.122	0.000	0.05147	2.73	ug/L
[>	Ho	165	1744774.553	2.276	1744774.553			ug/L
[	Tl	205	266.673	6.670	0.000	0.00208	27.22	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		96.668			
[	Mn	55					
[>	In	115		96.653			
[	Sb	121					
[	Sb	123					
[>	Ho	165		97.146			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31466

Report Date/Time: Friday, January 11, 2002 13:20:54

Page 1

Method 200.8 - Summary Report

Sample ID: AD31549

Sample Date/Time: Friday, January 11, 2002 13:22:35

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31549.2037

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.667	68.750	0.000	0.00188	891.01	ug/L
[>	Sc-1	45	604952.530	1.874	604952.530			ug/L
[	Mn	55	543170.776	2.087	0.897	31.09290	0.75	ug/L
[>	In	115	1256259.067	1.859	1256259.067			ug/L
[	Sb	121	448.016	1.464	0.000	0.05121	3.38	ug/L
[	Sb	123	357.648	5.195	0.000	0.05349	4.50	ug/L
[>	Ho	165	1717615.729	0.728	1717615.729			ug/L
[	Tl	205	261.673	11.442	0.000	0.00204	51.59	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		95.277			
[	Mn	55					
[>	In	115		95.259			
[	Sb	121					
[	Sb	123					
[>	Ho	165		95.634			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31549

Report Date/Time: Friday, January 11, 2002 13:23:07

Page 1

Method 200.8 - Summary Report

Sample ID: CCV CAL BK

Sample Date/Time: Friday, January 11, 2002 13:25:36

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV CAL BK.2038

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	9.000	40.062	-0.000	-0.00324	251.74	ug/L
[>	Sc-1	45	591573.710	2.343	591573.710			ug/L
[	Mn	55	638.029	5.285	0.000	0.00323	72.41	ug/L
[>	In	115	1246301.858	1.787	1246301.858			ug/L
[	Sb	121	65.334	19.015	0.000	0.00399	40.97	ug/L
[	Sb	123	50.432	17.426	0.000	0.00414	38.31	ug/L
[>	Ho	165	1680632.082	1.949	1680632.082			ug/L
[	Tl	205	249.006	2.811	0.000	0.00178	5.68	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		93.170			
[	Mn	55					
[>	In	115		94.504			
[	Sb	121					
[	Sb	123					
[>	Ho	165		93.575			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV CAL BK

Report Date/Time: Friday, January 11, 2002 13:26:09

Page 1

Method 200.8 - Summary Report

Sample ID: CCV STD1 30 PPB

Sample Date/Time: Friday, January 11, 2002 13:27:39

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV STD1 30 PPB.2039

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14530.418	1.918	0.024	30.28752	2.05	ug/L
[>	Sc-1	45	606299.017	2.217	606299.017			ug/L
[	Mn	55	519996.541	2.387	0.857	29.70913	3.30	ug/L
[>	In	115	1267695.534	1.447	1267695.534			ug/L
[	Sb	121	239714.720	0.194	0.189	29.34464	1.46	ug/L
[	Sb	123	184922.097	0.777	0.146	29.50090	2.16	ug/L
[>	Ho	165	1716049.251	1.592	1716049.251			ug/L
[	Tl	205	800386.965	3.105	0.466	30.14204	2.63	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		95.489			
[	Mn	55					
[>	In	115		96.127			
[	Sb	121					
[	Sb	123					
[>	Ho	165		95.547			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV STD1 30 PPB

Report Date/Time: Friday, January 11, 2002 13:28:12

Page 1

Method 200.8 - Summary Report

Sample ID: CCV STD2 60 PPB

Sample Date/Time: Friday, January 11, 2002 13:29:43

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV STD2 60 PPB.2040

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	28543.051	3.044	0.046	58.73098	3.10	ug/L
[>	Sc-1	45	614331.198	0.975	614331.198			ug/L
[	Mn	55	1023326.710	0.580	1.665	57.71608	0.66	ug/L
[>	In	115	1268802.467	2.129	1268802.467			ug/L
[	Sb	121	482539.933	0.818	0.380	59.03756	2.81	ug/L
[	Sb	123	376970.200	0.706	0.297	60.10107	2.73	ug/L
[>	Ho	165	1729949.127	0.983	1729949.127			ug/L
[	Tl	205	1592831.126	4.616	0.920	59.49670	3.64	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		96.754			
[	Mn	55					
[>	In	115		96.211			
[	Sb	121					
[	Sb	123					
[>	Ho	165		96.320			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV STD2 60 PPB

Report Date/Time: Friday, January 11, 2002 13:30:16

Page 1

Method 200.8 - Summary Report

Sample ID: AD31550

Sample Date/Time: Friday, January 11, 2002 13:32:18

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31550.2041

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	9.667	21.534	-0.000	-0.00298	135.75	ug/L
[>	Sc-1	45	623434.522	1.518	623434.522			ug/L
[	Mn	55	556842.039	2.208	0.892	30.92879	0.90	ug/L
[>	In	115	1303383.301	1.192	1303383.301			ug/L
[	Sb	121	759.373	5.855	0.001	0.08623	5.11	ug/L
[	Sb	123	582.789	3.879	0.000	0.08637	4.36	ug/L
[>	Ho	165	1756507.796	0.856	1756507.796			ug/L
[	Tl	205	414.014	3.321	0.000	0.00744	5.53	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		98.187			
[	Mn	55					
[>	In	115		98.833			
[	Sb	121					
[	Sb	123					
[>	Ho	165		97.799			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31550

Report Date/Time: Friday, January 11, 2002 13:32:51

Page 1

Method 200.8 - Summary Report

Sample ID: AD31550

Sample Date/Time: Friday, January 11, 2002 13:34:21

Sample Type: QC Duplicate

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31550.2042

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL): 5e+001

Diluted To Volume (mL): 1e+002

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	12.000	22.048	0.000	0.00526	226.25	ug/L
[>	Sc-1	45	604289.618	4.805	604289.618			ug/L
[	Mn	55	289529.771	1.857	0.479	33.20330	5.13	ug/L
[>	In	115	1287098.498	1.147	1287098.498			ug/L
[	Sb	121	367.344	3.011	0.000	0.08027	3.61	ug/L
[	Sb	123	293.560	3.405	0.000	0.08414	4.53	ug/L
[>	Ho	165	1749485.888	3.022	1749485.888			ug/L
[	Tl	205	277.340	9.294	0.000	0.00488	31.33	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					722.645
[>	Sc-1	45		95.172			
[	Mn	55					7.093
[>	In	115		97.598			
[	Sb	121					7.152
[	Sb	123					2.617
[>	Ho	165		97.408			
[	Tl	205					41.415

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31550

Report Date/Time: Friday, January 11, 2002 13:34:54

Page 1

Method 200.8 - Summary Report

Sample ID: AD31550

Sample Date/Time: Friday, January 11, 2002 13:36:25

Sample Type: Spike - 2

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31550.2043

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	25628.011	2.539	0.042	53.29290	2.38	ug/L
[>	Sc-1	45	608065.632	3.564	608065.632			ug/L
[	Mn	55	1368824.148	2.323	2.252	78.07170	4.23	ug/L
[>	In	115	1265669.867	1.538	1265669.867			ug/L
[	Sb	121	426656.085	0.552	0.337	52.31847	1.89	ug/L
[	Sb	123	327508.002	0.757	0.259	52.32872	0.99	ug/L
[>	Ho	165	1728238.334	1.450	1728238.334			ug/L
[	Tl	205	1312741.780	1.958	0.759	49.09252	1.11	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9			106.592		
[>	Sc-1	45		95.767			
[	Mn	55			94.286		
[>	In	115		95.973			
[	Sb	121			104.464		
[	Sb	123			104.485		
[>	Ho	165		96.225			
[	Tl	205			98.170		

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31550

Report Date/Time: Friday, January 11, 2002 13:36:58

Page 1

Method 200.8 - Summary Report

Sample ID: AD31551

Sample Date/Time: Friday, January 11, 2002 13:38:38

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31551.2044

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	13.000	27.735	0.000	0.00456	174.46	ug/L
[>	Sc-1	45	607599.681	4.234	607599.681			ug/L
[	Mn	55	1595267.020	1.940	2.626	91.04597	2.30	ug/L
[>	In	115	1252902.750	0.472	1252902.750			ug/L
[	Sb	121	1084.076	5.353	0.001	0.13012	5.54	ug/L
[	Sb	123	820.731	2.352	0.001	0.12842	2.91	ug/L
[>	Ho	165	1735985.857	1.636	1735985.857			ug/L
[	Tl	205	434.348	10.766	0.000	0.00837	19.56	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		95.694			
[	Mn	55					
[>	In	115		95.005			
[	Sb	121					
[	Sb	123					
[>	Ho	165		96.657			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31551

Report Date/Time: Friday, January 11, 2002 13:39:10

Page 1

Method 200.8 - Summary Report

Sample ID: AD31631

Sample Date/Time: Friday, January 11, 2002 13:40:41

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31631.2045

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14.667	27.555	0.000	0.00825	107.33	ug/L
>	Sc-1	45	602858.993	2.006	602858.993			ug/L
[	Mn	55	1846394.029	1.329	3.063	106.17341	2.45	ug/L
>	In	115	1240528.218	1.460	1240528.218			ug/L
	Sb	121	776.374	4.582	0.001	0.09296	4.46	ug/L
[	Sb	123	578.619	2.693	0.000	0.09026	2.18	ug/L
>	Ho	165	1683919.229	1.086	1683919.229			ug/L
[	Tl	205	322.676	3.414	0.000	0.00459	11.45	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		94.947			
[	Mn	55					
>	In	115		94.067			
	Sb	121					
[	Sb	123					
>	Ho	165		93.758			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31631

Report Date/Time: Friday, January 11, 2002 13:41:14

Page 1

Method 200.8 - Summary Report

Sample ID: AD31632

Sample Date/Time: Friday, January 11, 2002 13:42:45

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31632.2046

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.667	23.593	0.000	0.00004	16325.26	ug/L
[>	Sc-1	45	599934.589	4.453	599934.589			ug/L
[	Mn	55	564397.033	0.328	0.941	32.62399	4.55	ug/L
[>	In	115	1237704.732	2.181	1237704.732			ug/L
[	Sb	121	515.020	2.027	0.000	0.06046	4.54	ug/L
[	Sb	123	404.857	6.303	0.000	0.06206	5.21	ug/L
[>	Ho	165	1692199.238	1.032	1692199.238			ug/L
[	Tl	205	280.007	6.888	0.000	0.00290	28.30	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rei. % Diff
[	Be	9					
[>	Sc-1	45		94.486			
[	Mn	55					
[>	In	115		93.853			
[	Sb	121					
[	Sb	123					
[>	Ho	165		94.219			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31632

Report Date/Time: Friday, January 11, 2002 13:43:18

Page 1

Method 200.8 - Summary Report

Sample ID: AD31633

Sample Date/Time: Friday, January 11, 2002 13:46:02

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31633.2047

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	12.000	33.333	0.000	0.00239	351.37	ug/L
[>	Sc-1	45	608730.599	1.231	608730.599			ug/L
[	Mn	55	597033.695	2.975	0.980	33.96487	2.22	ug/L
[>	In	115	1245539.049	1.357	1245539.049			ug/L
[	Sb	121	485.018	3.407	0.000	0.05628	4.01	ug/L
[	Sb	123	374.861	8.041	0.000	0.05685	10.11	ug/L
[>	Ho	165	1700040.067	1.597	1700040.067			ug/L
[	Tl	205	237.672	3.999	0.000	0.00124	33.05	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		95.872			
[	Mn	55					
[>	In	115		94.447			
[	Sb	121					
[	Sb	123					
[>	Ho	165		94.655			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31633

Report Date/Time: Friday, January 11, 2002 13:46:34

Page 1

Method 200.8 - Summary Report

Sample ID: AD31670

Sample Date/Time: Friday, January 11, 2002 14:32:46

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31670.2048

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14.000	28.571	0.000	0.00596	137.34	ug/L
[>	Sc-1	45	620464.421	1.478	620464.421			ug/L
[	Mn	55	663858.528	2.251	1.069	37.05550	0.78	ug/L
[>	In	115	1311679.635	2.272	1311679.635			ug/L
[	Sb	121	427.681	11.497	0.000	0.04638	10.33	ug/L
[	Sb	123	328.895	1.704	0.000	0.04665	2.44	ug/L
[>	Ho	165	1798452.695	0.280	1798452.695			ug/L
[	Tl	205	193.337	5.995	-0.000	-0.00085	48.26	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		97.720			
[	Mn	55					
[>	In	115		99.462			
[	Sb	121					
[	Sb	123					
[>	Ho	165		100.135			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31670

Report Date/Time: Friday, January 11, 2002 14:33:19

Page 1

Method 200.8 - Summary Report

Sample ID: AD31671

Sample Date/Time: Friday, January 11, 2002 14:34:50

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31671.2049

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	48.709	-0.000	-0.00077	1448.25	ug/L
>	Sc-1	45	603362.754	4.705	603362.754			ug/L
[	Mn	55	635766.270	0.462	1.054	36.55182	4.92	ug/L
>	In	115	1261517.429	0.396	1261517.429			ug/L
	Sb	121	431.681	10.444	0.000	0.04896	11.74	ug/L
[	Sb	123	326.435	4.967	0.000	0.04826	5.16	ug/L
>	Ho	165	1715429.977	1.524	1715429.977			ug/L
[	Tl	205	259.006	6.380	0.000	0.00197	39.49	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		95.026			
[	Mn	55					
>	In	115		95.658			
	Sb	121					
[	Sb	123					
>	Ho	165		95.512			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31671

Report Date/Time: Friday, January 11, 2002 14:35:22

Page 1

Method 200.8 - Summary Report

Sample ID: AD31672

Sample Date/Time: Friday, January 11, 2002 14:36:54

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31672.2050

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.667	27.553	0.000	0.00214	353.40	ug/L
[>	Sc-1	45	601728.559	3.770	601728.559			ug/L
[	Mn	55	1973139.184	1.768	3.280	113.71453	2.67	ug/L
[>	In	115	1263095.415	1.098	1263095.415			ug/L
[	Sb	121	598.692	6.607	0.000	0.06941	7.38	ug/L
[	Sb	123	440.156	6.210	0.000	0.06642	7.06	ug/L
[>	Ho	165	1764627.242	2.420	1764627.242			ug/L
[	Tl	205	270.673	2.621	0.000	0.00212	1.01	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		94.769			
[	Mn	55					
[>	In	115		95.778			
[	Sb	121					
[	Sb	123					
[>	Ho	165		98.251			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31672

Report Date/Time: Friday, January 11, 2002 14:37:26

Page 1

Method 200.8 - Summary Report

Sample ID: AD31674

Sample Date/Time: Friday, January 11, 2002 14:39:21

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31674.2051

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	65.873	-0.000	-0.00099	1482.70	ug/L
[>	Sc-1	45	611137.154	2.439	611137.154			ug/L
[	Mn	55	629329.997	1.868	1.029	35.67634	2.55	ug/L
[>	In	115	1284473.376	1.497	1284473.376			ug/L
[	Sb	121	419.014	4.932	0.000	0.04645	3.95	ug/L
[	Sb	123	335.483	2.698	0.000	0.04877	3.99	ug/L
[>	Ho	165	1721809.156	1.286	1721809.156			ug/L
[	Tl	205	233.339	7.671	0.000	0.00096	70.07	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		96.251			
[	Mn	55					
[>	In	115		97.399			
[	Sb	121					
[	Sb	123					
[>	Ho	165		95.867			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31674

Report Date/Time: Friday, January 11, 2002 14:39:54

Page 1

Method 200.8 - Summary Report

Sample ID: AD31675

Sample Date/Time: Friday, January 11, 2002 14:41:25

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31675.2052

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	12.667	12.059	0.000	0.00435	99.80	ug/L
[>	Sc-1	45	597741.768	4.292	597741.768			ug/L
[	Mn	55	631058.983	4.489	1.057	36.63943	7.76	ug/L
[>	In	115	1270172.843	0.092	1270172.843			ug/L
[	Sb	121	414.680	1.371	0.000	0.04651	1.58	ug/L
[	Sb	123	317.531	4.297	0.000	0.04649	4.57	ug/L
[>	Ho	165	1738778.128	1.924	1738778.128			ug/L
[	Tl	205	226.672	4.003	0.000	0.00063	79.28	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rei. % Diff
[	Be	9					
[>	Sc-1	45		94.141			
[	Mn	55					
[>	In	115		96.315			
[	Sb	121					
[	Sb	123					
[>	Ho	165		96.812			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31675

Report Date/Time: Friday, January 11, 2002 14:41:57

Page 1

Method 200.8 - Summary Report**Sample ID: CCV CAL BK**

Sample Date/Time: Friday, January 11, 2002 14:45:11

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV CAL BK.2053

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	5.587	-0.000	-0.00063	219.64	ug/L
[>	Sc-1	45	595900.955	2.238	595900.955			ug/L
[	Mn	55	594.692	5.149	0.000	0.00040	262.48	ug/L
[>	In	115	1228192.331	0.734	1228192.331			ug/L
[	Sb	121	41.667	6.040	0.000	0.00111	29.09	ug/L
[	Sb	123	28.813	25.212	0.000	0.00068	175.05	ug/L
[>	Ho	165	1678240.017	1.887	1678240.017			ug/L
[	Tl	205	224.005	11.391	0.000	0.00084	134.19	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		93.851			
[	Mn	55					
[>	In	115		93.131			
[	Sb	121					
[	Sb	123					
[>	Ho	165		93.441			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV CAL BK

Report Date/Time: Friday, January 11, 2002 14:45:44

Page 1

Method 200.8 - Summary Report

Sample ID: CCV STD1 30 PPB

Sample Date/Time: Friday, January 11, 2002 14:47:40

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV STD1 30 PPB.2054

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14571.486	2.488	0.024	30.16189	3.13	ug/L
>	Sc-1	45	610570.746	1.842	610570.746			ug/L
[	Mn	55	515866.105	2.877	0.844	29.26006	2.95	ug/L
>	In	115	1250930.307	1.164	1250930.307			ug/L
	Sb	121	240288.062	0.911	0.192	29.80908	1.90	ug/L
[	Sb	123	186575.146	0.871	0.149	30.15966	1.27	ug/L
>	Ho	165	1707289.635	2.036	1707289.635			ug/L
[	Tl	205	798910.931	1.262	0.468	30.24462	1.05	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		96.161			
[	Mn	55					
>	In	115		94.855			
	Sb	121					
[	Sb	123					
>	Ho	165		95.059			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV STD1 30 PPB

Report Date/Time: Friday, January 11, 2002 14:48:12

Page 1

Method 200.8 - Summary Report

Sample ID: CCV STD2 60 PPB

Sample Date/Time: Friday, January 11, 2002 14:49:44

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV STD2 60 PPB.2055

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	28471.467	0.634	0.046	58.59843	1.91	ug/L
>	Sc-1	45	614339.940	2.480	614339.940			ug/L
[	Mn	55	1022493.227	1.785	1.664	57.68816	2.93	ug/L
>	In	115	1251073.326	0.792	1251073.326			ug/L
	Sb	121	480390.858	1.087	0.384	59.58585	1.16	ug/L
[	Sb	123	374396.171	1.121	0.299	60.51730	1.60	ug/L
>	Ho	165	1705348.059	2.267	1705348.059			ug/L
[	Tl	205	1570377.700	0.653	0.921	59.53357	1.75	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
>	Sc-1	45		96.755			
[	Mn	55					
>	In	115		94.866			
	Sb	121					
[	Sb	123					
>	Ho	165		94.951			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
>	In	115	Linear Thru Zero	
	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV STD2 60 PPB

Report Date/Time: Friday, January 11, 2002 14:50:16

Page 1

Method 200.8 - Summary Report

Sample ID: LFB 50 PPB

Sample Date/Time: Friday, January 11, 2002 14:51:47

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\LFB 50 PPB.2056

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	24107.165	1.115	0.041	51.76273	1.03	ug/L
[>	Sc-1	45	588717.080	1.918	588717.080			ug/L
[	Mn	55	849376.603	2.593	1.442	49.99860	3.48	ug/L
[>	In	115	1244025.503	0.486	1244025.503			ug/L
[	Sb	121	404708.221	0.239	0.325	50.48208	0.72	ug/L
[	Sb	123	312858.482	1.132	0.251	50.85195	0.66	ug/L
[>	Ho	165	1692831.455	2.277	1692831.455			ug/L
[	Tl	205	1303607.485	3.105	0.770	49.77272	2.27	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		92.720			
[	Mn	55					
[>	In	115		94.332			
[	Sb	121					
[	Sb	123					
[>	Ho	165		94.254			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: LFB 50 PPB

Report Date/Time: Friday, January 11, 2002 14:52:20

Page 1

Method 200.8 - Summary Report

Sample ID: AD31676

Sample Date/Time: Friday, January 11, 2002 14:53:52

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31676.2057

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	13.333	8.660	0.000	0.00619	50.75	ug/L
[>	Sc-1	45	587216.298	2.644	587216.298			ug/L
[	Mn	55	2190415.690	1.797	3.730	129.30314	0.98	ug/L
[>	In	115	1237395.628	1.616	1237395.628			ug/L
[	Sb	121	822.045	5.772	0.001	0.09891	4.79	ug/L
[	Sb	123	612.840	2.093	0.000	0.09608	0.56	ug/L
[>	Ho	165	1682055.527	2.264	1682055.527			ug/L
[	Tl	205	520.354	8.579	0.000	0.01218	10.44	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		92.483			
[	Mn	55					
[>	In	115		93.829			
[	Sb	121					
[	Sb	123					
[>	Ho	165		93.654			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31676

Report Date/Time: Friday, January 11, 2002 14:54:25

Page 1

Method 200.8 - Summary Report

Sample ID: AD31676

Sample Date/Time: Friday, January 11, 2002 14:55:55

Sample Type: QC Duplicate

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31676.2058

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL): 5e+001

Diluted To Volume (mL): 1e+002

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.667	24.744	0.000	0.00229	488.65	ug/L
[>	Sc-1	45	621060.912	1.409	621060.912			ug/L
[	Mn	55	1262247.629	3.593	2.031	140.82122	2.21	ug/L
[>	In	115	1282440.991	0.186	1282440.991			ug/L
[	Sb	121	436.015	4.903	0.000	0.09720	5.15	ug/L
[	Sb	123	326.441	6.879	0.000	0.09483	7.66	ug/L
[>	Ho	165	1783438.671	0.551	1783438.671			ug/L
[	Tl	205	285.341	2.280	0.000	0.00509	10.61	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					92.022
[>	Sc-1	45		97.814			
[	Mn	55					✓ 8.528
[>	In	115		97.245			
[	Sb	121					1.739
[	Sb	123					1.316
[>	Ho	165		99.299			
[	Tl	205					82.117

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31676

Report Date/Time: Friday, January 11, 2002 14:56:28

Page 1

Method 200.8 - Summary Report

Sample ID: AD31676

Sample Date/Time: Friday, January 11, 2002 14:57:59

Sample Type: Spike - 2

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31676.2059

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	23806.057	3.342	0.040	50.79536	1.19	ug/L
[>	Sc-1	45	592266.188	2.246	592266.188			ug/L
[	Mn	55	2866950.711	1.880	4.840	167.79140	0.62	ug/L
[>	In	115	1211104.303	1.040	1211104.303			ug/L
[	Sb	121	417741.898	1.128	0.345	53.52308	0.26	ug/L
[	Sb	123	319552.698	0.720	0.264	53.35472	0.56	ug/L
[>	Ho	165	1680248.769	1.042	1680248.769			ug/L
[	Tl	205	1263409.103	1.166	0.752	48.60514	2.19	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9			101.578		
[>	Sc-1	45		93.279			
[	Mn	55			76.977		
[>	In	115		91.836			
[	Sb	121			106.848		
[	Sb	123			106.517		
[>	Ho	165		93.553			
[	Tl	205			97.186		

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31676

Report Date/Time: Friday, January 11, 2002 14:58:32

Page 1

Method 200.8 - Summary Report

Sample ID: AD31676

Sample Date/Time: Friday, January 11, 2002 15:00:02

Sample Type: Spike - 2

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AD31676.2060

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	23396.337	4.224	0.041	51.22714	2.78	ug/L
[>	Sc-1	45	577148.741	2.107	577148.741			ug/L
[	Mn	55	2853120.778	2.134	4.942	171.34488	0.06	ug/L
[>	In	115	1202807.943	0.842	1202807.943			ug/L
[	Sb	121	410766.698	1.478	0.341	52.99223	1.14	ug/L
[	Sb	123	316857.079	0.283	0.263	53.26969	0.57	ug/L
[>	Ho	165	1674997.885	1.830	1674997.885			ug/L
[	Tl	205	1243236.146	3.531	0.742	47.96164	1.71	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9			102.442		
[>	Sc-1	45		90.898			
[	Mn	55			84.083		
[>	In	115		91.206			
[	Sb	121			105.787		
[	Sb	123			106.347		
[>	Ho	165		93.261			
[	Tl	205			95.899		

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AD31676

Report Date/Time: Friday, January 11, 2002 15:00:35

Page 1

Method 200.8 - Summary Report

Sample ID: AE00100

Sample Date/Time: Friday, January 11, 2002 15:05:51

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00100.2061

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	15.667	40.537	0.000	0.01034	123.84	ug/L
[>	Sc-1	45	599658.191	1.766	599658.191			ug/L
[	Mn	55	647576.294	2.742	1.079	37.40574	2.42	ug/L
[>	In	115	1238564.642	1.665	1238564.642			ug/L
[	Sb	121	1029.735	0.716	0.001	0.12488	1.57	ug/L
[	Sb	123	790.443	4.774	0.001	0.12508	6.47	ug/L
[>	Ho	165	1709600.393	2.040	1709600.393			ug/L
[	Tl	205	560.356	11.437	0.000	0.01342	20.49	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		94.443			
[	Mn	55					
[>	In	115		93.918			
[	Sb	121					
[	Sb	123					
[>	Ho	165		95.187			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00100

Report Date/Time: Friday, January 11, 2002 15:06:24

Page 1

Method 200.8 - Summary Report**Sample ID: AE00101**

Sample Date/Time: Friday, January 11, 2002 15:07:55

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00101.2062

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	16.000	22.535	0.000	0.01159	75.17	ug/L
] >	Sc-1	45	595675.240	3.398	595675.240			ug/L
[	Mn	55	700181.216	3.594	1.175	40.73909	4.31	ug/L
] >	In	115	1231160.827	1.124	1231160.827			ug/L
[	Sb	121	728.036	1.584	0.001	0.08761	0.48	ug/L
[	Sb	123	546.738	3.798	0.000	0.08575	4.48	ug/L
] >	Ho	165	1682741.229	2.023	1682741.229			ug/L
[	Tl	205	345.677	11.288	0.000	0.00550	32.12	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
] >	Sc-1	45		93.815			
[	Mn	55					
] >	In	115		93.356			
[	Sb	121					
[	Sb	123					
] >	Ho	165		93.692			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
] >	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
] >	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
] >	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00101

Report Date/Time: Friday, January 11, 2002 15:08:28

Page 1

Method 200.8 - Summary Report

Sample ID: AE00102

Sample Date/Time: Friday, January 11, 2002 15:09:59

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00102.2063

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	14.783	-0.000	-0.00008	4586.25	ug/L
[>	Sc-1	45	581457.225	1.978	581457.225			ug/L
[	Mn	55	2722948.225	1.889	4.682	162.31931	0.76	ug/L
[>	In	115	1175142.559	2.085	1175142.559			ug/L
[	Sb	121	767.373	4.457	0.001	0.09727	6.72	ug/L
[	Sb	123	587.293	5.511	0.000	0.09710	7.76	ug/L
[>	Ho	165	1642230.000	1.295	1642230.000			ug/L
[	Tl	205	305.342	7.500	0.000	0.00422	20.85	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		91.576			
[	Mn	55					
[>	In	115		89.109			
[	Sb	121					
[	Sb	123					
[>	Ho	165		91.436			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00102

Report Date/Time: Friday, January 11, 2002 15:10:32

Page 1

Method 200.8 - Summary Report

Sample ID: AE00104

Sample Date/Time: Friday, January 11, 2002 15:12:03

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00104.2064

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.000	27.273	0.000	0.00113	598.98	ug/L
[>	Sc-1	45	588776.000	2.254	588776.000			ug/L
[	Mn	55	645916.684	1.943	1.096	38.00025	0.60	ug/L
[>	In	115	1231869.636	0.936	1231869.636			ug/L
[	Sb	121	560.690	5.470	0.000	0.06646	5.18	ug/L
[	Sb	123	419.572	3.614	0.000	0.06480	3.21	ug/L
[>	Ho	165	1660239.674	1.908	1660239.674			ug/L
[	Tl	205	269.673	3.733	0.000	0.00271	21.51	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		92.729			
[	Mn	55					
[>	In	115		93.410			
[	Sb	121					
[	Sb	123					
[>	Ho	165		92.439			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00104

Report Date/Time: Friday, January 11, 2002 15:12:36

Page 1

Method 200.8 - Summary Report

Sample ID: AE00105

Sample Date/Time: Friday, January 11, 2002 15:14:06

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00105.2065

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14.000	12.372	0.000	0.00741	54.01	ug/L
] >	Sc-1	45	591051.932	1.849	591051.932			ug/L
[	Mn	55	683409.335	1.471	1.156	40.06032	2.41	ug/L
] >	In	115	1198244.652	1.315	1198244.652			ug/L
[	Sb	121	520.020	6.190	0.000	0.06318	5.94	ug/L
[	Sb	123	381.437	3.855	0.000	0.06033	5.10	ug/L
] >	Ho	165	1639483.414	2.542	1639483.414			ug/L
[	Tl	205	247.339	1.235	0.000	0.00196	18.25	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
] >	Sc-1	45		93.087			
[	Mn	55					
] >	In	115		90.860			
[	Sb	121					
[	Sb	123					
] >	Ho	165		91.283			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
] >	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
] >	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
] >	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00105

Report Date/Time: Friday, January 11, 2002 15:14:39

Page 1

Method 200.8 - Summary Report

Sample ID: AE00106

Sample Date/Time: Friday, January 11, 2002 15:16:10

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00106.2066

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	33.986	0.000	0.00039	2080.65	ug/L
[>	Sc-1	45	571088.123	2.083	571088.123			ug/L
[	Mn	55	2821501.719	2.571	4.941	171.28299	2.97	ug/L
[>	In	115	1185769.445	1.089	1185769.445			ug/L
[	Sb	121	608.693	2.008	0.000	0.07550	1.55	ug/L
[	Sb	123	479.406	2.752	0.000	0.07769	1.98	ug/L
[>	Ho	165	1655715.729	0.797	1655715.729			ug/L
[	Tl	205	237.005	4.975	0.000	0.00145	27.08	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		89.943			
[	Mn	55					
[>	In	115		89.914			
[	Sb	121					
[	Sb	123					
[>	Ho	165		92.187			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00106

Report Date/Time: Friday, January 11, 2002 15:16:43

Page 1

Method 200.8 - Summary Report

Sample ID: AE00108

Sample Date/Time: Friday, January 11, 2002 15:18:14

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00108.2067

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.333	29.565	-0.000	-0.00019	3450.53	ug/L
[>	Sc-1	45	583279.797	1.512	583279.797			ug/L
[	Mn	55	606774.360	0.816	1.039	36.03688	2.00	ug/L
[>	In	115	1196477.276	2.604	1196477.276			ug/L
[	Sb	121	493.352	3.922	0.000	0.05988	6.10	ug/L
[	Sb	123	368.818	1.454	0.000	0.05831	4.17	ug/L
[>	Ho	165	1631352.939	0.587	1631352.939			ug/L
[	Tl	205	226.005	8.637	0.000	0.00116	69.71	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		91.863			
[	Mn	55					
[>	In	115		90.726			
[	Sb	121					
[	Sb	123					
[>	Ho	165		90.831			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00108

Report Date/Time: Friday, January 11, 2002 15:18:47

Page 1

Method 200.8 - Summary Report

Sample ID: AE00109

Sample Date/Time: Friday, January 11, 2002 15:20:18

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00109.2068

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	8.333	24.980	-0.000	-0.00467	107.15	ug/L
[>	Sc-1	45	592019.875	3.912	592019.875			ug/L
[	Mn	55	665581.585	1.432	1.124	38.97595	3.63	ug/L
[>	In	115	1235964.077	1.055	1235964.077			ug/L
[	Sb	121	465.017	1.197	0.000	0.05423	2.02	ug/L
[	Sb	123	346.860	1.719	0.000	0.05268	1.04	ug/L
[>	Ho	165	1684284.666	0.723	1684284.666			ug/L
[	Tl	205	210.338	2.866	0.000	0.00028	95.98	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		93.240			
[	Mn	55					
[>	In	115		93.721			
[	Sb	121					
[	Sb	123					
[>	Ho	165		93.778			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00109

Report Date/Time: Friday, January 11, 2002 15:20:50

Page 1

Method 200.8 - Summary Report

Sample ID: AE00110

Sample Date/Time: Friday, January 11, 2002 15:22:22

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\AE00110.2069

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	11.333	28.364	0.000	0.00172	358.39	ug/L
[>	Sc-1	45	587257.079	3.227	587257.079			ug/L
[	Mn	55	2500758.855	1.409	4.259	147.66269	2.32	ug/L
[>	In	115	1197700.932	1.457	1197700.932			ug/L
[	Sb	121	594.692	2.356	0.000	0.07290	2.71	ug/L
[	Sb	123	447.156	5.258	0.000	0.07145	5.95	ug/L
[>	Ho	165	1658013.959	2.342	1658013.959			ug/L
[	Tl	205	229.005	8.960	0.000	0.00113	72.50	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		92.490			
[	Mn	55					
[>	In	115		90.819			
[	Sb	121					
[	Sb	123					
[>	Ho	165		92.315			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: AE00110

Report Date/Time: Friday, January 11, 2002 15:22:54

Page 1

Method 200.8 - Summary Report

Sample ID: CCV CAL BK

Sample Date/Time: Friday, January 11, 2002 15:40:26

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV CAL BK.2070

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	10.000	10.000	-0.000	-0.00067	292.27	ug/L
[>	Sc-1	45	577203.712	2.777	577203.712			ug/L
[	Mn	55	600.026	4.207	0.000	0.00187	85.41	ug/L
[>	In	115	1200550.282	2.040	1200550.282			ug/L
[	Sb	121	276.007	6.191	0.000	0.03152	6.40	ug/L
[	Sb	123	211.690	7.658	0.000	0.03158	7.36	ug/L
[>	Ho	165	1645249.294	1.113	1645249.294			ug/L
[	Tl	205	252.339	5.049	0.000	0.00212	22.25	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		90.906			
[	Mn	55					
[>	In	115		91.035			
[	Sb	121					
[	Sb	123					
[>	Ho	165		91.605			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV CAL BK

Report Date/Time: Friday, January 11, 2002 15:40:59

Page 1

Method 200.8 - Summary Report**Sample ID: CCV STD1 30 PPB**

Sample Date/Time: Friday, January 11, 2002 15:42:30

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV STD1 30 PPB.2071

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	14230.615	3.694	0.024	30.76786	4.65	ug/L
[>	Sc-1	45	584691.640	2.404	584691.640			ug/L
[	Mn	55	506712.436	1.109	0.866	30.02102	2.54	ug/L
[>	In	115	1199106.846	1.206	1199106.846			ug/L
[	Sb	121	232370.438	0.948	0.194	30.07190	1.66	ug/L
[	Sb	123	178011.444	1.152	0.148	30.01688	0.07	ug/L
[>	Ho	165	1656575.043	1.569	1656575.043			ug/L
[	Tl	205	779707.040	2.253	0.471	30.41574	0.88	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		92.086			
[	Mn	55					
[>	In	115		90.926			
[	Sb	121					
[	Sb	123					
[>	Ho	165		92.235			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV STD1 30 PPB

Report Date/Time: Friday, January 11, 2002 15:43:03

Page 1

Method 200.8 - Summary Report**Sample ID: CCV STD2 60 PPB**

Sample Date/Time: Friday, January 11, 2002 15:44:34

Sample Type: Sample

Sample Description:

Number of Replicates: 3

Batch ID:

Method File: d:\elandata\Method\NEWE200.8.mth

Dataset File: d:\elandata\Dataset\daily performance\CCV STD2 60 PPB.2072

Sample Prep Volume (mL):

Initial Sample Quantity (mg):

Aliquot Volume (mL):

Diluted To Volume (mL):

Concentration Results

	Analyte	Mass	Meas. Intens. Mean	Meas. Intens. RSD	Net Intens. Mean	Conc. Mean	Conc. RSD	Sample Unit
[	Be	9	28098.649	1.929	0.047	60.04454	5.58	ug/L
[>	Sc-1	45	592422.472	4.317	592422.472			ug/L
[	Mn	55	989634.094	1.353	1.672	57.95021	4.53	ug/L
[>	In	115	1213723.414	0.874	1213723.414			ug/L
[	Sb	121	465931.461	1.848	0.384	59.57429	2.27	ug/L
[	Sb	123	361968.977	1.711	0.298	60.30413	1.28	ug/L
[>	Ho	165	1657051.685	1.811	1657051.685			ug/L
[	Tl	205	1524090.706	1.483	0.920	59.44995	0.41	ug/L

QC Calculated Values

	Analyte	Mass	QC Std % Recovery	Int Std % Recovery	Spike % Recovery	Dilution % Diff	Dup. Rel. % Diff
[	Be	9					
[>	Sc-1	45		93.303			
[	Mn	55					
[>	In	115		92.034			
[	Sb	121					
[	Sb	123					
[>	Ho	165		92.262			
[	Tl	205					

Calibration

	Analyte	Mass	Curve Type	Correlation Coefficient
[	Be	9	Linear Thru Zero	0.999992
[>	Sc-1	45	Linear Thru Zero	
[	Mn	55	Linear Thru Zero	0.999983
[>	In	115	Linear Thru Zero	
[	Sb	121	Linear Thru Zero	0.999998
[	Sb	123	Linear Thru Zero	1.000000
[>	Ho	165	Linear Thru Zero	
[	Tl	205	Linear Thru Zero	0.999977

Sample ID: CCV STD2 60 PPB

Report Date/Time: Friday, January 11, 2002 15:45:06

Page 1