

Method 508

Date Extracted 9/15/01

Page 79

GC/MS Unknown

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	count		
	<u>LB</u>	<u>Seal</u>		<u>1 ml</u>			<u>2 ml</u>					
<u>22485</u>	<u>Sample</u>									<u>1</u>		
<u>22486</u>										<u>2</u>		
<u>22487</u>										<u>3</u>		
<u>22488</u>										<u>4</u>		
<u>22489</u>										<u>5</u>		
<u>22490</u>										<u>6</u>		
<u>22491</u>										<u>7</u>		
<u>22492</u>										<u>8</u>		
<u>22493</u>										<u>9</u>		
<u>22494</u>										<u>10</u>		
<u>22495</u>										<u>11</u>		
<u>22496</u>										<u>12</u>		
<u>22485</u>	<u>Sample spike</u>				<u>✓</u>						<u>NEW</u>	
<u>22485</u>	<u>1 dup</u>				<u>✓</u>							
	<u>M&L Chloro</u>				<u>✓</u>							
	<u>Ref. Chloro</u>				<u>✓</u>							
	<u>PRL 1232</u>				<u>✓</u>							
Remarks												
<u>Spike + spike dup : 10000</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors												
Extraction type - <u>sep. funnel</u> Solv. lot# <u>Chloro</u> NaCl lot# <u>T12612</u> Na2SO4 lot# <u>V18612</u>												
STANDARD	INT STD LOG#											
CALIBRATION	<u>40X: 1809</u>	<u>chl: 1808</u>	<u>rest: 1820d</u>									
REF/SPIKE	<u>40X: 1678a</u>	<u>rest: 2062a</u>	<u>chl: 1786</u>									
MDL/MDRF	<u>—</u>	<u>—</u>	<u>—</u>									
Q. CHECK	<u>—</u>											
INT. STD.	<u>—</u>											
SURROGATE	<u>1834a</u>											
LPC	<u>1867d</u>											
ENDRIN/DDT	<u>1866</u>											

Analyst MAH. H.
revision 10/15/99

Verified by My

Supv. approved Dr
9/24/01

Sequence Name: C:\HPCHEM\1\SEQUENCE\SEP2101.S
 Comment: 208
 Operator: GALOS
 Data Path: C:\HPCHEM\1\DATA\sep2101\
 Pre-Seq Cmd:
 Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch
☒ Full Method ☒ Inject Anyway
☐ Reprocessing Only ☐ Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	HP060501	
2	Sample	2	22470	HP060501	GC/MS unknown
3	Sample	3	22471	HP060501	GC/MS unknown
4	Sample	4	22472	HP060501	GC/MS unknown
5	Sample	5	22485	HP060501	GC/MS unknown
6	Sample	6	22486	HP060501	GC/MS unknown
7	Sample	7	22487	HP060501	GC/MS unknown
8	Sample	8	22488	HP060501	GC/MS unknown
9	Sample	9	22489	HP060501	GC/MS unknown

Sequence Name: C:\HPCHEM\1\SEQUENCE\SEP2201.S

Comment: 208

Operator: GALOS

Data Path: C:\HPCHEM\1\DATA\sep2101\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

☒ Full Method☒ Inject Anyway☐ Reprocessing Only☐ Don't Inject

Line Type	Vial	DataFile	Method	Sample Name
1 Sample	1	DFTPP	HP060501	
2 Sample	2	22491	HP060501	GC/MS unknown
3 Sample	3	22490	HP060501	GC/MS unknown

DFTPP

Data File : C:\HPCHEM\1\DATA\SEP2201\DFTPP.D

Acq On : 22 Sep 2001 11:14 am

Sample :

Misc :

MS Integration Params: RTEINT.P

Method : C:\HPCHEM\1\METHODS\NP072501.M (RTE Integrator)


Title : 208 625/TCLP calibration table 7/25/01 C1 IS 5

Vial: 1

Operator: GALOS

Inst : GC/MS Ins

Multiplr: 1.00



Spectrum Information: Average of 26.508 to 26.545 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	30.6	144323	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	41.5	195808	PASS
70	69	0.00	2	0.4	791	PASS
127	198	40	60	41.5	195896	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	471539	PASS
199	198	5	9	6.6	31167	PASS
275	198	10	30	26.3	123798	PASS
365	198	1	99	3.1	14706	PASS
441	443	1	99	82.1	65381	PASS
442	198	40	110	84.8	399718	PASS
443	442	17	23	19.9	79664	PASS

DFTPP.D NP072501.M

Mon Sep 24 09:55:48 2001