

Method 508

Date Extracted 11/8/01

Page

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count
	LB	500mL		7mL			2mL		5.7	
	Ref chl + Pest									
	Ref TOX									
26833	S TOX				1mL					
26833	S D TOX									
26845	S CHL									
26845	S D CHL									
26846	S Pest									
26846	S D Pest									
	REF 1248									
	26833									1
	26845									2
	26846									3
	26847									4
	26882									5
	26883									6
	26884									7
	26885									8
	26886									9
	26887									10
26988	S CHL				1mL					
26988	S D CHL									
26991	S TOX									
26991	S D TOX									
26992	S Pest									
26992	S D Pest									
Remarks										
Reviewed by <u>45</u> 3/13/02										
drip rate must be between 5-15 ml/min for liq-liq extractors										
drip rate -										
Extraction type - <u>Sep funnel</u> Solv. lot# <u>V42917</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V32603</u>										
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	<u>TOX: 1992</u>	<u>CHL: 1993</u>	<u>Pest 20026</u>							
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE	<u>1839b</u>									
LPC										
ENDRIN/DDT										

Analyst Cm/ml/AD/E
revision 4/12/01

Verified by

Reviewed by DW
Date 12/24/11

