

Method 508

Date Extracted 10/9/01

Page

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH = 7	Count
	LB	500ml		1ml			2ml			
	Ref chl + Rest									
	Ref TOX				1ml					
27055	S chl									
27055	SD chl									
27056	S TOX									
27057	SD TOX									
27013	S Peot									
27013	SD Peot									
	PHL 1254									
	27013									1
	27014									2
	27015									3
	27032									4
	27055									5
	27056									6
	27110									7
	27111									8
	27112									9
	27113									10
Remarks										
27032 - missing										
27013 Spk drop - lab accident - lost.										
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type - <u>sp funnel</u> Solv. lot# <u>V42E17</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V32603</u>										
STANDARD	INT STD LOG#									
CALIBRATION										
REF/SPIKE	<u>TOX: 1992</u>	<u>chl: 1993</u>	<u>Peot: 20026</u>							
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE	<u>1839b</u>									
LPC										
ENDRIN/DDT										

Reviewed by AB 3/13/12

Analyst MLAD/CM/E
revision 4/12/01

Verified by

Reviewed by DN
Date 12/28/11