

AUG-12-2002 17:01
Aug 09 02 11:14sDDC
PK Microbiology718 391 1499
856 489 4085
P.08
p.2**P & K Microbiology Services, Inc.**

Tel: 856-489-4455

Fax: 856-489-4085

The Environmental Microbiology Specialists

7 Allison Drive, Cherry Hill, New Jersey 08003

Client: ATC Associates, Inc., New York, NY

Project ID: WTC/DDC

Date sampled: July 25, 2002

Date of inoculation: July 26, 2002

Samples submitted By: Wagdi

Date characterization completed: August 6, 2002

P&K Report No.: 020726-081

Fungal / Bacterial Analysis (Culture Method)**Wipe Samples**

Wipe Samples				Fungal / Bacterial ID	Colony counts	Conc. ** (CFU / in2)	Percentage* (%)
Sample ID	Area (inch2)	Medium used	Dilution factor				
0725-WTC-W1	1	MEA	100X	Fungi			
				Acromonium strictum	37	3,700	76
				Penicillium	11	1,100	22
				Trichoderma koningii	1	100	2
				Total: 4,900			
		TSA	16400X	Bacteria			
				Bacillus	8	131,200	3
				gram negative bacteria and others	236	3,870,400	97
				Total: 4,001,600			
0725-WTC-W2	1	MEA	100X	Fungi			
				Penicillium	7	700	70
				Trichoderma koningii	3	300	30
				Total: 1,000			
				TSA	16400X	Bacteria	
		gram negative bacteria and others	278			4,559,200	100
		Total: 4,559,200					
		0725-WTC-W3	1	MEA	100X	Fungi	
Acromonium strictum	28					2,800	27
Aspergillus versicolor	5					500	5
Cladosporium	9					900	9
Penicillium	56					5,600	55
Rhodotorula	4					400	4
Total: 10,200							
TSA	16400X			Bacteria			
				Bacillus	2	32,800	< 1
				gram negative bacteria and others	297	4,870,800	98
				Staphylococcus	4	65,600	1
				Total: 4,969,200			
0725-WTC-W4	1	MEA	100X	Fungi			
				Penicillium	7	700	64
				Trichoderma koningii	4	400	36
				Total: 1,100			
		TSA	100X	Bacteria			
				Bacillus	4	400	100
				Total: 400			

AUG-12-2002 17:01
Aug 09 02 11:14aDDC
PK Microbiology718 391 1499
855 489 4085P.09
P.3

P & K Microbiology Services, Inc.

P&K Report No.: 020726-081 Page 2

Wipe Samples

Sample ID	Area (inch ²)	Medium used	Dilution factor	Fungal / Bacterial ID	Colony counts	Conc. ** (CFU / in ²)	Percentage* (%)
0725-WTC-W5	1	MEA	100X	Fungi			
				Aspergillus versicolor	24	2,400	77
				Penicillium	3	300	10
				sterile fungi	3	300	10
				Trichoderma koningii	1	100	3
						Total: 3,100	
		TSA	4100X	Bacteria			
				gram negative bacteria and others	312	1,279,200	65
				Pseudomonas sp. non aeruginosa	169	692,900	35
						Total: 1,972,100	
0725-WTC-W6	1	MEA	100X	Fungi			
				Aspergillus versicolor	1	100	< 1
				Cladosporium	24	2,400	17
				Penicillium	3	300	2
				yeasts	116	11,600	81
						Total: 14,400	
		TSA	100X	Bacteria			
				Bacillus	9	900	4
				gram negative bacteria and others	236	23,600	93
				Micrococcus luteus	4	400	2
				Rhodococcus	6	600	2
						Total: 25,500	
0725-WTC-W7	1	MEA	100X	Fungi			
				Basidiomycetes	1	100	25
				Penicillium	3	300	75
						Total: 400	
		TSA	100X	Bacteria			
				gram negative bacteria and others	72	7,200	97
				Micrococcus luteus	2	200	3
						Total: 7,400	

* Percentage of each group of fungi / bacteria in total population.

** Concentration is (CFU/Sample) if sample amount is NA.

Media types: Cellulose agar (CA), Czapek cellulose agar (CCA), cornmeal agar (CMA), 2% malt extract agar (MEA), 2% malt extract agar plus 20% sucrose (MEA+S), inhibitory mold agar (IMA), Pseudomonas isolation agar (PIA), rose bengal agar (RBA), sabouraud dextrose agar (SDA), tryptic soy agar (TSA), nutrient agar (NTA), Blood agar (BA), Staphylococcus Medium 110 (Staphy).

The detection limit of fungal and bacterial analysis using culture methods is one colony. The quantitation limits vary from analysis to analysis and from processing procedure to processing procedure. Contact us to determine your quantitation limits.

Approved by: _____

Chin S. Yang, Ph.D., Microbiologist

Quality control checked by: _____