


Original Extraction

Method 508 mod.

Date Extracted 2/12/02

Page 83

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
LB		500ml		1ml			1ml	7B	7			
REF					1ml			2B				
REF					↓			4A				
3075								1A				
3076								3A				
3077								4B				
3078								6B				
3178								6A				
3179								8A				
3180								3B				
3181								7A				
3182								8B				
3183		↓		↓			↓	2A	↓			
Remarks												
Reviewed by 												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type <u>Sep Fan</u> Solv. lot# <u>V41470</u> NaCl lot# <u>V39618</u> Na2SO4 lot# <u>V39400</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	91759											
MDL/MDRF												
Q. CHECK												
INT. STD.	022702											
SURROGATE	GCMS	02	12	02								
LPC												
ENDRIN/DDT												

Analyst CH

Revision 4140100

Reviewed by

3/2/02
OW

Method 508 mod.

Date Extracted 3/6/02

2nd Extraction unnecessary
Page 134

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Sol add	IS add	Final vol.	T.V. pos.	pH	Count		
LB		500ml		1ml			1ml	10A	7			
REF ₁					1ml			10B				
REF ₂					↓			6A				
3075								9B				
3076								8B				
3077								4B				
3078	FB							2B				
3178	FB							2A				
3179								1A				
3180								9A				
3181								4A				
3182								1B				
3183								8A				
Remarks												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type <u>Sep fun</u>		Solv. lot# <u>V41470</u>		NaCl lot# <u>V39618</u>		Na2SO4 lot# <u>V46623</u>						
STANDARD	INT STD LOG#											
CALIBRATION												
REV SPIKE	<u>2193</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	<u>GCMS 2160 D</u>											
LPC												
ENDRIN/DDT												

Analyst EA
Revision 4/12/01

Verified by _____

Reviewed by DN
Date 3/22/07

DFTPP

•
•
•
Data File : C:\MSDCHEM\1\DATA\031802\DFTPP.D
Acq On : 19 Mar 2002 2:11 am
Sample :
Misc :
MS Integration Params: BENZO.P

Vial: 1
Operator: McLean
Inst : Instrumen
Multiplr: 1.00

Method : C:\MSDCHEM\1\5973N\5080302.M (RTE Integrator)
Title : Quantitation For Method 508gcscan

Spectrum Information: Scan 2288

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	10	80	50.5	373824	PASS
68	69	0.00	2	1.8	5665	PASS
69	198	0.00	100	41.8	309248	PASS
70	69	0.00	2	0.4	1098	PASS
127	198	10	80	45.7	338560	PASS
197	198	0.00	2	0.0	0	PASS
198	198	50	100	100.0	740608	PASS
199	198	5	9	6.1	44816	PASS
275	198	10	60	21.4	158464	PASS
365	198	1	100	2.7	19792	PASS
441	443	0.01	100	81.5	94384	PASS
442	198	50	100	84.9 ✓	628928	PASS
443	442	15	24	18.4	115800	PASS

DFTPP.D 5080302.M Tue Mar 19 21:40:54 2002