

Method 508 mod

Date Extracted 12/20/01

Page 80

ref

\$ 6cms 16369

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count											
LB	✓	500ml		1ml			1ml		77												
Ref	✓				0.5ml																
30605	✓ - ✓																				
27288	P bottle	not included																			
31019	✓ -																				
30903	✓ - ✓																				
31021	✓ - ✓																				
31022	✓ - ✓																				
30902	✓ - ✓																				
31023	✓ - ✓																				
27288 - restructured internal already																					
Reviewed by <i>[Signature]</i> 5/15/02																					
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -																					
Extraction type - <u>sep fun</u>		Solv. lot#	<u>V39466</u>		NaCl lot#	<u>V11599</u>		Na2SO4 lot# <u>V39H00</u>													
STANDARD	INT STD LOG#																				
CALIBRATION																					
REF/SPIKE	<u>91759</u>																				
MDL/MDRF																					
Q. CHECK																					
INT. STD.																					
SURROGATE	<u>121901</u>																				
LPC																					
ENDRIN/DDT																					

Analyst LF
revision 4/12/01

Verified by _____

Reviewed by DN
Date 1/18/02

Sequence Name: C:\HPCHEM\1\SEQUENCE\jan0802.S

Comment: INST 209; 30Mx0.25umX0.50

Operator: 209 GALOS

Data Path: C:\HPCHEM\1\data\jan0802\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0103	
2	Sample	2	CAL40	GCMS0103	calib mix
3	Sample	3	28355	GCMS0103	
4	Sample	5	28327	GCMS0103	508
5	Sample	6	29466	GCMS0103	
6	Sample	4	28355A	GCMS0103	
7	Sample	7	LB	GCMS0103	gcms scan 12/20
8	Sample	8	REF	GCMS0103	gcms scan 12/20
9	Sample	9	31076	GCMS0103	gcms scan 12/20
10	Sample	10	31073	GCMS0103	gcms scan 12/20
11	Sample	11	31075	GCMS0103	gcms scan 12/20
12	Sample	12	31074	GCMS0103	gcms scan 12/20
13	Sample	13	31017	GCMS0103	gcms scan 12/20
14	Sample	14	31018	GCMS0103	gcms scan 12/20
15	Sample	15	31020	GCMS0103	gcms scan 12/20
16	Sample	16	31024	GCMS0103	gcms scan 12/20
17	Sample	17	30605	GCMS0103	gcms scan 12/20
18	Sample	18	LBB	GCMS0103	gcms scan 12/20a
19	Sample	19	REFB	GCMS0103	gcms scan 12/20a
20	Sample	20	30902	GCMS0103	gcms scan 12/20a
21	Sample	1	DFTPP1	GCMS0103	
22	Sample	21	31022	GCMS0103	gcms scan 12/20a
23	Sample	22	31021	GCMS0103	gcms scan 12/20a
24	Sample	23	30903	GCMS0103	gcms scan 12/20a
25	Sample	24	31019	GCMS0103	gcms scan 12/20a
26	Sample	25	31023	GCMS0103	gcms scan 12/20a
27	Sample	2	CAL40A	GCMS0103	calib mix
28	Sample	26	27288	GCMS0103	gcms scan 12/20a
29	Sample	4	LBC	GCMS0103	gcms scan 12/24
30	Sample	5	REFC	GCMS0103	gcms scan 12/24
31	Sample	6	31278	GCMS0103	gcms scan 12/24
32	Sample	7	31129	GCMS0103	gcms scan 12/24
33	Sample	8	31132	GCMS0103	gcms scan 12/24
34	Sample	9	31255	GCMS0103	gcms scan 12/24
35	Sample	10	31253	GCMS0103	gcms scan 12/24
36	Sample	11	31344	GCMS0103	gcms scan 12/24
37	Sample	12	31254	GCMS0103	gcms scan 12/24
38	Sample	13	31248	GCMS0103	gcms scan 12/24
39	Sample	14	31252	GCMS0103	gcms scan 12/24
40	Sample	15	31130	GCMS0103	gcms scan 12/24
41	Sample	16	31251	GCMS0103	gcms scan 12/24
42	Sample	17	LBD	GCMS0103	gcms scan 12/21

Last Modified: Thu Jan 10 08:16:55 2002

Page: 1

DFTPP

Data File : C:\HPCHEM\1\DATA\JAN0802\DFTPP1.D

Acq On : 9 Jan 2002 8:50 am

Sample :

Misc :

MS Integration Params: GOOFY.P

Method : C:\HPCHEM\1\METHODS\GCMS1126.M (RTE Integrator)


Title : 209 625/TCLP calibration table

Vial: 1

Operator: 209 GALOS

Inst : GC/MS 209

Multiplr: 1.00



Spectrum Information: Average of 26.051 to 26.079 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	42.3	122572	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	51.9	150636	PASS
70	69	0.00	2	0.2	247	PASS
127	198	40	60	48.2	139872	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	289984	PASS
199	198	5	9	6.6	19144	PASS
275	198	10	30	20.9	60566	PASS
365	198	1	100	2.4	7018	PASS
441	443	0.01	100	89.2	29128	PASS
442	198	40	110	58.8	170494	PASS
443	442	17	23	19.1	32644	PASS