

Method 508 modified Date Extracted 12/18/01

Page 76

REF GGCMS-16707

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	SPK add	IS add	Final vol.	T.V. pos.	pH	Count		
LB ✓		500M		1M			1M	3B	7			
REF ✓					5M			7A				
29231 ✓	See 12/3							7B				
29232 ✓	See 12/3							3A				
30644 ✓	See 2/1							8A				
30643 ✓								6A				
30646 ✓	See 2/4							7A				
30647 ✓								7B				
30645 ✓								6B				
30649 ✓	See 2/1							8B				
30706 ✓								7B				
29322 ✓								8B				
29324 ✓	See 2/1							3B				
28323								7A				
Remarks												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type <u>Sip Run</u> Solv. lot# <u>V39466</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V39400</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	<u>91759</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	<u>GGCMS</u>	<u>112601</u>										
LPC												
ENDRIN/DDT												

Analyst J. Death
revision 4/12/01

Verified by _____

Reviewed by DW
Date 2/26/02


Sequence Name: C:\HPCHEM\1\SEQUENCE\Dec2701.S

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	6	29228	GCMS1126 GC/MS	SCAN 12/7
45	Sample	7	29286	GCMS1126 GC/MS	SCAN 12/7
46	Sample	8	29292	GCMS1126 GC/MS	SCAN 12/7
47	Sample	9	29284	GCMS1126 GC/MS	SCAN 12/7
48	Sample	10	29465	GCMS1126 GC/MS	SCAN 12/7
49	Sample	11	29344	GCMS1126 GC/MS	SCAN 12/7
50	Sample	12	29343	GCMS1126 GC/MS	SCAN 12/7
51	Sample	13	29519	GCMS1126 GC/MS	SCAN 12/7
52	Sample	14	29562	GCMS1126 GC/MS	SCAN 12/7
53	Sample	15	30485	GCMS1126 GC/MS	SCAN
54	Sample	16	30488	GCMS1126 GC/MS	SCAN
55	Sample	17	30256	GCMS1126 GC/MS	SCAN
56	Sample	18	29735	GCMS1126 GC/MS	SCAN
57	Sample	19	29791	GCMS1126 GC/MS	SCAN
58	Sample	1	DFTPP3	GCMS1126	
59	Sample	2	CAL40	GCMS1126	
60	Sample	38	LBD	GCMS1126 GC/MS	SCAN 12/18
61	Sample	39	REFD	GCMS1126 GC/MS	SCAN 12/18
62	Sample	40	30704	GCMS1126 GC/MS	SCAN 12/18
63	Sample	41	30707	GCMS1126 GC/MS	SCAN 12/18
64	Sample	42	30641	GCMS1126 GC/MS	SCAN 12/18
65	Sample	43	30703	GCMS1126 GC/MS	SCAN 12/18
66	Sample	44	30649	GCMS1126 GC/MS	SCAN 12/18
67	Sample	45	30642	GCMS1126 GC/MS	SCAN 12/18
68	Sample	46	30650	GCMS1126 GC/MS	SCAN 12/18
69	Sample	47	30709	GCMS1126 GC/MS	SCAN 12/18
70	Sample	48	30710	GCMS1126 GC/MS	SCAN 12/18
71	Sample	49	30708	GCMS1126 GC/MS	SCAN 12/18
72	Sample	50	30651	GCMS1126 GC/MS	SCAN 12/18
73	Sample	51	LBE	GCMS1126 GC/MS	SCAN 12/18 A
74	Sample	52	REFE	GCMS1126 GC/MS	SCAN 12/18 A
75	Sample	53	30645	GCMS1126 GC/MS	SCAN 12/18 A
76	Sample	2	CAL40A	GCMS1126	
77	Sample	54	30706	GCMS1126 GC/MS	SCAN 12/18 A
78	Sample	55	28324	GCMS1126 GC/MS	SCAN 12/18 A
79	Sample	56	30648	GCMS1126 GC/MS	SCAN 12/18 A
80	Sample	57	29231	GCMS1126 GC/MS	SCAN 12/18 A
81	Sample	58	30647	GCMS1126 GC/MS	SCAN 12/18 A
82	Sample	59	30643	GCMS1126 GC/MS	SCAN 12/18 A
83	Sample	60	29232	GCMS1126 GC/MS	SCAN 12/18 A
84	Sample	61	30644	GCMS1126 GC/MS	SCAN 12/18 A
85	Sample	62	28322	GCMS1126 GC/MS	SCAN 12/18 A
86	Sample	63	28323	GCMS1126 GC/MS	SCAN 12/18 A
87	Sample	64	30646	GCMS1126 GC/MS	SCAN 12/18 A

DFTPP

Data File : M:\INSTRU~1\209\DATA\DEC2701\DFTPP3.D
 Acq On : 29 Dec 2001 5:42 pm
 Sample :
 Misc :
 MS Integration Params: RTEINT.P
 Method : C:\HPCHEM\1\METHODS\GCMS1126.M (RTE Integrator)
 Title : 209 625/TCLP calibration table

Vial: 1
 Operator: 209 GALOS
 Inst : GC/MS 209
 Multiplr: 1.00



Spectrum Information: Average of 25.816 to 26.927 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	49.6	29676	PASS
68	69	0.00	2	0.0	9	PASS
69	198	0.00	100	60.1	35971	PASS
70	69	0.00	2	0.2	79	PASS
127	198	40	60	54.2	32448	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	59866	PASS
199	198	5	9	6.7	4013	PASS
275	198	10	30	20.6	12323	PASS
365	198	1	100	2.7	1611	PASS
441	443	0.01	100	92.2	5264	PASS
442	198	40	110	49.0	29340	PASS
443	442	17	23	19.5	5710	PASS