

Method 508 MOD.

Date Extracted 4/23/02

Page 72

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
L.B.		500 mL		1 mL			1 mL	1B	7			
Ref 1					1 mL			2A				
Ref 2					1 mL			9A				
9568	F.B.							8A				
9569								6B				
9570								7B				
9571								1A				
9572	F.B.							8B				
9574								2B				
9575								7A				
9576								10A				
9577								9B				
9578								10B				
Remarks: Many low Surrogates will be reextracted												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type - SEP. Fun Solv. lot# <u>X01E19</u> NaCl lot# <u>V39614</u> Na2SO4 lot# <u>V46648</u>												
STANDARD	INT STD LOG#											
CALIBRATION	<u>1105020</u>											
REF SPIKE	<u>2279E</u>											
MDL/MDRF												
Q. CHECK												
INT. STD.	<u>042602</u>											
SURROGATE												
LPC												
ENDRIN/DDT												

Analyst A.D.
revision 4/12/01

Verified by my

Reviewed by _____
Date _____

Comment:
 Operator:
 Data Path: C:\HPCHEM\1\DATA\MAY0302\
 Pre-Seq Cmd:
 Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch
 (X) Full Method (X) Inject Anyway
 () Reprocessing Only () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0501	
2	Sample	2	CAL20	GCMS0501	
3	Sample	3	LB	GCMS0501	GC/MS SCAN 4/24
4	Sample	4	REF1	GCMS0501	GC/MS SCAN 4/24
5	Sample	5	REF2	GCMS0501	GC/MS SCAN 4/24
6	Sample	6	8837	GCMS0501	RE-EXT.
7	Sample	7	8838	GCMS0501	RE-EXT.
8	Sample	8	8934	GCMS0501	RE-EXT.
9	Sample	9	9579	GCMS0501	GC/MS SCAN 4/24
10	Sample	10	9580	GCMS0501	GC/MS SCAN 4/24
11	Sample	11	9581	GCMS0501	GC/MS SCAN 4/24 FB
12	Sample	12	9585	GCMS0501	GC/MS SCAN 4/24
13	Sample	13	9586	GCMS0501	GC/MS SCAN 4/24
14	Sample	14	9587	GCMS0501	GC/MS SCAN 4/24 FB
15	Sample	15	LBA	GCMS0501	GC/MS SCAN 4/22
16	Sample	16	REFA1	GCMS0501	GC/MS SCAN 4/22
17	Sample	17	REFA2	GCMS0501	GC/MS SCAN 4/22
18	Sample	18	9517	GCMS0501	GC/MS SCAN 4/22
19	Sample	19	9518	GCMS0501	GC/MS SCAN 4/22
20	Sample	20	9519	GCMS0501	GC/MS SCAN 4/22
21	Sample	21	9520	GCMS0501	GC/MS SCAN 4/22
22	Sample	22	8271	GCMS0501	RE-EXT. FB
23	Sample	23	8430	GCMS0501	RE-EXT.
24	Sample	24	8485	GCMS0501	RE-EXT.
25	Sample	25	8487	GCMS0501	RE-EXT.
26	Sample	26	9446	GCMS0501	RE-EXT.
27	Sample	1	DFTPP2	GCMS0501	
28	Sample	100	CAL40	GCMS0501	
29	Sample	27	LBB	GCMS0501	GC/MS SCAN 4/26
30	Sample	28	REFB1	GCMS0501	GC/MS SCAN 4/26
31	Sample	29	REFB2	GCMS0501	GC/MS SCAN 4/26
32	Sample	30	9151	GCMS0501	RE-EXT.
33	Sample	31	9154	GCMS0501	RE-EXT.
34	Sample	32	9155	GCMS0501	RE-EXT.
35	Sample	33	9331	GCMS0501	RE-EXT.
36	Sample	34	9332	GCMS0501	RE-EXT.
37	Sample	35	10062	GCMS0501	GC/MS SCAN 4/26
38	Sample	36	10063	GCMS0501	GC/MS SCAN 4/26
39	Sample	37	10064	GCMS0501	GC/MS SCAN 4/26
40	Sample	38	10065	GCMS0501	GC/MS SCAN 4/26
41	Sample	39	LBC	GCMS0501	GC/MS SCAN 4/23
42	Sample	40	REFC1	GCMS0501	GC/MS SCAN 4/23
43	Sample	41	REFC2	GCMS0501	GC/MS SCAN 4/23

Line	Type	Vial	DataFile	Method	Sample Name
44	Sample	42 9568	GCMS0501	GC/MS	SCAN 4/23 FB
45	Sample	43 9569	GCMS0501	GC/MS	SCAN 4/23
46	Sample	44 9570	GCMS0501	GC/MS	SCAN 4/23
47	Sample	45 9571	GCMS0501	GC/MS	SCAN 4/23
48	Sample	46 9572	GCMS0501	GC/MS	SCAN 4/23 FB
49	Sample	47 9574	GCMS0501	GC/MS	SCAN 4/23
50	Sample	48 9575	GCMS0501	GC/MS	SCAN 4/23
51	Sample	49 9576	GCMS0501	GC/MS	SCAN 4/23
52	Sample	50 9577	GCMS0501	GC/MS	SCAN 4/23
53	Sample	51 9578	GCMS0501	GC/MS	SCAN 4/23

DFTPP

Data File : C:\HPCHEM\1\DATA\MAY0302\DFTPP2.D

Acq On : 4 May 2002 3:34 pm

Sample :

Misc :

MS Integration Params: GOOFY.P

Method : C:\HPCHEM\1\METHODS\GCMS0501.M (RTE Integrator)


Title : 209 625/TCLP calibration table

Vial: 1

Operator:

Inst : 209

Multiplr: 1.00



Spectrum Information: Scan 2371

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	43.3	102688	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	51.4	121912	PASS
70	69	0.00	2	0.0	0	PASS
127	198	40	60	41.4	98256	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	237056	PASS
199	198	5	9	6.5	15505	PASS
275	198	10	30	20.7	49064	PASS
365	198	1	100	1.4	3390	PASS
441	443	0.01	100	84.8	24704	PASS
442	198	40	110	61.1	144896	PASS
443	442	17	23	20.1	29120	PASS

DFTPP2.D GCMS0501.M

Tue May 07 14:11:52 2002