

Page 113

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
2R		500ml		1ml			1ml	9A	7			
Rep 1					1ml			9B				
Rep 2					1ml			10A				
10901								10B				
10980								1A				
10981								2A				
10982								6A				
10983								6B				
10984								2B				
10985								7A				
Remarks												
Reviewed by CB 8/22/22												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate - Extraction type - SER Fwd Solv. lot# X06E27 NaCl lot# V39614 Na2SO4 lot# V396648												
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	2279 E											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	TCMX	22849	5/7/22									
LPC												
ENDRIN/DDT												

Analyst Z.D.
revision 4/12/01

Verified by

Reviewed by _____
Date _____

DFTPP

Data File : C:\MSDCHEM\1\DATA\051402\DFTPP.D

Acq On : 14 May 2002 9:31 am

Sample :

Misc :

MS Integration Params: EVENTS.E

Vial: 1

Operator: Mclean

Inst : Instrumen

Multiplr: 1.00

Method : C:\MSDCHEM\1\METHODS\050102.M (Chemstation Integrator)

Title : 508 Modified GC/MS scan

Spectrum Information: Scan 3602

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	47.1	141120	PASS
68	69	0.00	2	1.7	1967	PASS
69	198	0.00	100	39.5	118376	PASS
70	69	0.00	2	0.5	586	PASS
127	198	40	60	49.8	149248	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	299904	PASS
199	198	5	9	6.6	19672	PASS
275	198	10	30	19.7	59024	PASS
365	198	1	99	2.3	6779	PASS
441	443	1	100	75.4	29168	PASS
442	198	40	110	63.8	191360	PASS
443	442	17	23	20.2	38680	PASS

DFTPP.D 050102.M

Thu May 16 14:33:57 2002