

Method 508 MODDate Extracted 5-10-02

IDATC#

18599

Page

45

REF

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
	LAB BLANK	500µl		1ml			1ml	6A	7			
	REF1				1ml			6B				
	REF2				1ml			7A				
	10851							7B				
	10852							8A				
	10890							8B				
	10891							6B				
	10892							6A				
	10893							7A				
	10894							7B				
Remarks												
<u>two</u> USED 1 ML GRADUATED PIPETS TO ADD SURROGATES AND REFERENCE.												
Reviewed by, <u>CB 8/22/02</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type <u>SEP. Funnel</u> Solv. lot# <u>X06E27</u> NaCl lot# Na2SO4 lot# <u>V46648</u>												
STANDARD INT STD LOG#												
CALIBRATION												
REF SPIKE <u>2279E</u>												
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE <u>TCMX 22849 5/2/02</u>												
LPC												
ENDRIN/DDT												

Analyst BN
revision 4/12/01

Verified by _____

Reviewed by DW
Date 5/21/02

Sequence Name: C:\HPCHEM\1\SEQUENCE\MAY1502.S

Comment:

Operator:

Data Path: C:\HPCHEM\1\DATA\MAY1502\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch


☒ Full Method☒ Inject Anyway☐ Reprocessing Only☐ Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0510	
2	Sample	2	CAL40	GCMS0510	
3	Sample	3	LB	GCMS0510	GC/MS SCAN 5/10
4	Sample	4	REF1	GCMS0510	GC/MS SCAN 5/10
5	Sample	5	REF2	GCMS0510	GC/MS SCAN 5/10
6	Sample	6	10851	GCMS0510	GC/MS SCAN 5/10
7	Sample	7	10852	GCMS0510	GC/MS SCAN 5/10
8	Sample	8	10890	GCMS0510	GC/MS SCAN 5/10
9	Sample	9	10891	GCMS0510	GC/MS SCAN 5/10
10	Sample	10	10892	GCMS0510	GC/MS SCAN 5/10
11	Sample	11	10893	GCMS0510	GC/MS SCAN 5/10
12	Sample	12	10894	GCMS0510	GC/MS SCAN 5/10

DFTPP

Data File : C:\HPCHEM\1\DATA\MAY1502\DFTPP.D
 Acq On : 15 May 2002 8:34 am
 Sample :
 Misc :
 MS Integration Params: GOOFY.P
 Method : C:\HPCHEM\1\METHODS\GCMS0510.M (RTE Integrator)
 Title : 209 625/TCLP calibration table

Vial: 1
 Operator:
 Inst : 209
 Multiplr: 1.00



Spectrum Information: Average of 26.734 to 26.752 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	58.6	90981	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	63.1	98072	PASS
70	69	0.00	2	0.2	162	PASS
127	198	40	60	49.9	77592	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	155371	PASS
199	198	5	9	6.8	10519	PASS
275	198	10	30	21.1	32821	PASS
365	198	1	100	1.7	2654	PASS
441	443	0.01	100	92.5	25709	PASS
442	198	40	110	91.7	142520	PASS
443	442	17	23	19.5	27792	PASS