

Method 508 MOD.

Date Extracted 5/13/02

Page 119

BGCMS/18526

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
L.B.		500 mL		1 mL			1 mL	3B	7			
Ref 1					1 mL			4B				
Ref 2					1 mL			1A				
11157								3A				
11158								2A				
11159	FB							10A				
11278								10B				
11279								9B				
11280								2B				
11281	F.B.							4A				
Remarks												
Reviewed by, <u>C/D 8/23/02</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type - Sep. Fun.		Solv. lot# X05E10		NaCl lot# V39614		Na2SO4 lot# V46648						
STANDARD	INT STD LOG#											
CALIBRATION												
REF SPIKE	2279E											
MDL/MDRF												
Q. CHECK												
INT. STD.	2286F											
PROBATE	TCHX 22849 (5/7/02)											
LPC												
DRIN/DDT												

yst A.D.
ion 4/12/01

Verified by _____

Reviewed by

Dy 5/22/02

Sequence Name: C:\MSDCHEM\1\SEQUENCE\020515.S

Comment: GCMS scan

Operator: Nefliu

Data Path: C:\MSDchem\1\DATA\020515\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run On A Barcode Mismatch

☒ Full Method ☒ Inject Anyway☐ Reprocessing Only ☐ Don't Inject


Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	508FAST	
2	Sample	2	0515C5	508SCAN1	
3	Sample	3	0513LB1	508SCAN1	05/13 lab blk
4	Sample	4	0513LB2	508SCAN1	05/13 lab blk
5	Sample	5	0513REF1	508SCAN1	05/13 ref
6	Sample	6	0513REF2	508SCAN1	05/13 ref
7	Sample	7	0513REF3	508SCAN1	05/13 ref
8	Sample	8	0513REF4	508SCAN1	05/13 ref
9	Sample	9	11157	508SCAN1	sample
10	Sample	10	11158	508SCAN1	sample
11	Sample	11	11159	508SCAN1	sample
12	Sample	12	11278	508SCAN1	sample
13	Sample	13	11279	508SCAN1	sample
14	Sample	14	11280	508SCAN1	sample
15	Sample	15	11281	508SCAN1	sample
16	Sample	16	10986	508SCAN1	sample
17	Sample	17	10987	508SCAN1	sample
18	Sample	18	11152	508SCAN1	sample
19	Sample	19	11153	508SCAN1	sample
20	Sample	20	11154	508SCAN1	sample
21	Sample	21	11155	508SCAN1	sample
22	Sample	22	11156	508SCAN1	sample

Lab #	Filename	Description	Vial #		Inj. ul	Standards used - Lot #			Remarks	Initials
			front	rear		Calibration	Reference	Surrogate		
	SETPP	dftpp	1		1					MN
	075C5	dev	2			22816/22844/22860				
	073L61	eds blk	3					22846	22864	
	073L62	↓	4							
	073R41	ref	5				22790			
	↓ 2	↓	6							
	↓ 3	↓	7							
	↓ 4	↓	8							
	11157	sample	9							
	11158		10							
	11159		11							
	11278		12							
	11279		13							
	11280		14							
	11281		15							
	10986		16							
	10987		17							
	11152		18							
	11153		19							
	11154		20							
	11155		21							
	11156	↓	22							

DFTPP

Data File : C:\MSDCHEM\1\DATA\020515\DFTPP.D
 Acq On : 15 May 2002 11:01 am
 Sample :
 Misc :
 MS Integration Params: rteint.e
 Method : C:\MSDCHEM\1\METHODS\SCAN020507.M (RTE Integrator)
 Title : Method 508 GC/MS Scan

Vial: 1
 Operator: Nefliu
 Inst : Instrume
 Multiplr: 1.00



Spectrum Information: Average of 17.813 to 17.836 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	53.3	229226	PASS
68	69	0.00	2	1.5	3100	PASS
69	198	0.00	100	46.8	201304	PASS
70	69	0.00	2	0.5	1033	PASS
127	198	40	60	51.6	221866	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	430208	PASS
199	198	5	9	6.5	27912	PASS
275	198	10	30	24.8	106570	PASS
365	198	1	99	3.3	14184	PASS
441	443	1	100	75.1	52906	PASS
442	198	40	110	83.9	360938	PASS
443	442	17	23	19.5	70493	PASS

DFTPP.D SCAN020507.M

Thu May 16 10:01:21 2002

DFTPP

Data File : C:\MSDCHEM\1\DATA\020515\DFTPP.D

Vial: 1

Acq On : 15 May 2002 11:01 am

Operator: Nefliu

Sample :

Inst : Instrume


Misc :

Multiplr: 1.00

MS Integration Params: rteint.e

Method : C:\MSDCHEM\1\METHODS\SCAN020507.M (RTE Integrator)

Title : Method 508 GC/MS Scan



Spectrum Information: Average of 17.813 to 17.836 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	53.3	229226	PASS
68	69	0.00	2	1.5	3100	PASS
69	198	0.00	100	46.8	201304	PASS
70	69	0.00	2	0.5	1033	PASS
127	198	40	60	51.6	221866	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	430208	PASS
199	198	5	9	6.5	27912	PASS
275	198	10	30	24.8	106570	PASS
365	198	1	99	3.3	14184	PASS
441	443	1	100	75.1	52906	PASS
442	198	40	110	83.9	360938	PASS
443	442	17	23	19.5	70493	PASS

DFTPP.D SCAN020507.M

Thu May 16 10:01:26 2002