

Method 508 modifiedDate Extracted 11/21/01Page 12

ref

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
✓ LB ✓		500 mL		1.0 mL		✓	1 mL		7			
✓ Ref ✓					0.5 mL							
✓ 28263 ✓												
✓ 28269 ✓												
✓ 28267 ✓												
✓ 28271 ✓												
✓ 28275 ✓												
✓ 28273 ✓												
✓ 28277 ✓												
✓ 28413 ✓												
✓ 28433 ✓												
✓ 28431 ✓												
✓ 28432 ✓												
✓ 28434 ✓												
✓ 28430 ✓												
✓ 28427 ✓												
✓ 28471 ✓	spilled ✓	0										
✓ 28468 ✓												
✓ 28469 ✓												
✓ 28429 ✓												
✓ 28470 ✓												
✓ 28428 ✓												
Reviewed by, <u>JS 5/11/02</u> 28263 re extracted 11/27 28273 re extracted												
Remarks 28471 was spilled after being C.D. 28263 had hexane added by co-worker (accident). 15 mL CH₂Cl₂ was added, sample put back on C.D. 28430 vial was incorrectly filled by RW, must re-extract 28267 must be re-extracted												
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -												
Extraction type <u>sep fun</u> Solv. lot# <u>V39466</u> NaCl lot# <u>V1159</u> Na₂SO₄ lot# <u>V32603</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE	91759											
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE	24000											
LPC												
ENDRIN/DDT												

Analyst

LF/RW

revision 4/12/01

Verified by _____

Reviewed by

Date

DN

12/6/11