

Method 508 MOD.

Date Extracted 4/29/02

Page 87

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
L.B.		500ml		1ml			1ml	9A	7			
Ref 1					1ml			3B				
Ref 2					1ml			3A				
10066		↓						4B				
10067		475ml						1A				
10068	F.B.	500ml						2A				
10130		↓						10A				
10131								4A				
10132								9B				
10133								1B				
✓ 10134								5A				
10135								5B				
10136		↓		↓			↓	2B	↓			
Remarks												
All re-est except 10134.												
Reviewed by, <u>5/8/21/n</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors												
drip rate -												
Extraction type - <u>SEP. Fun.</u> Solv. lot# <u>X05E10</u> NaCl lot# <u>V 39614</u> Na2SO4 lot# <u>V 46648</u>												
STANDARD	INT STD LOG#											
CALIBRATION												
REMARKS	2279E											
MOI. DRF												
CHECK												
INT. STD.												
SURROGATE	GCMS -											
LPC												
ENDRIN/DDT												

Analyst A.D
revision 4/12/01

Verified by _____

Reviewed by 5/8/21/n
Date _____

Data File : C:\MSDCHEM\1\DATA\050202\DFTPA.D
Acq On : 2 May 2002 10:27 am
Sample :
Misc :
MS Integration Params: EVENTS.E

Vial: 1
Operator: Mclean
Inst : Instrumen
Multiplr: 1.00

Method : C:\MSDCHEM\1\METHODS\050102.M (Chemstation Integrator)
Title : 508 Modified GC/MS scan

Spectrum Information: Scan 3603

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	47.5	91000	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	41.1	78656	PASS
70	69	0.00	2	0.6	460	PASS
127	198	40	60	51.0	97680	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	191424	PASS
199	198	5	9	6.7	12756	PASS
275	198	10	30	19.2	36808	PASS
365	198	1	99	2.3	4381	PASS
441	443	1	100	74.6	18160	PASS
442	198	40	110	65.7	125688	PASS
443	442	17	23	19.4	24328	PASS

DFTPA.D 050102.M

Mon May 06 14:09:40 2002