

Page 175

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count
L.B.		500ml		1ml			1ml	1A	7	
Ref 1					1ml			1B		
Ref 2					1ml			2B		
10869								3A		
10810								9A		
10811								9B		
10812								10A		
10813								10B		
10814	FB							8A		
10815								8B		
Remarks										
Reviewed by, <i>[Signature]</i> 5/23/02										
drip rate must be between 5-15 ml/min for liq-liq extractors drip rate -										
Extraction type - SEP. Fw.		Solv.lot# X05E10		NaCl lot# V37614		Na2SO4 lot# V46648				
STANDARD INT STD LOG#										
CALIBRATION										
(REF) SPIKE 2279 E										
MDL/MDRF										
Q. CHECK										
INT. STD.										
SURROGATE TCMX (IST 440)		2284 5/2/02								
LPC										
ENDRIN/DDT										

Analyst KB.
revision 4/12/01

Verified by _____

Reviewed by _____
Date _____

Sequence Name: C:\MSDCHEM\1\SEQUENCE\051402.S

Comment:

Operator: Mclean

Data Path: C:\MSDCHEM\1\DATA\051402\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

☒ (X) Full Method☒ (X) Inject Anyway☐ () Reprocessing Only☐ () Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	508SCAN	
2	Sample	2	CAL50	508SCAN	
3	Sample	3	0510LB	508SCAN	5/10 kb
4	Sample	4	510REF1	508SCAN	5/10 kb
5	Sample	5	510REF2	508SCAN	5/10kb
6	Sample	6	TCMX	508SCAN	40
7	Sample	7	MC	508SCAN	
8	Sample	8	10901	508SCAN	5/10 kb
9	Sample	9	10980	508SCAN	5/10 kb
10	Sample	10	10981	508SCAN	5/10 kb
11	Sample	11	10982	508SCAN	5/10 kb
12	Sample	12	10983	508SCAN	5/10 kb
13	Sample	13	10984	508SCAN	5/10 kb
14	Sample	14	10985	508SCAN	5/10 kb
15	Sample	1	DFTPPB	508SCAN	
16	Sample	15	CAL20	508SCAN	
17	Sample	16	0507LB	508SCAN	5/7 kb
18	Sample	17	507REF1	508SCAN	5/7 kb
19	Sample	18	507REF2	508SCAN	5/7 kb
20	Sample	19	10809	508SCAN	5/7 kb
21	Sample	20	10810	508SCAN	5/7 kb
22	Sample	21	10811	508SCAN	5/7 kb
23	Sample	22	10812	508SCAN	5/7 kb
24	Sample	23	10813	508SCAN	5/7 kb
25	Sample	24	10814	508SCAN	5/7 kb
26	Sample	25	10815	508SCAN	5/7 kb
27	Sample	2	CAL50B	508SCAN	5/7 kb

DFTPP

Data File : C:\MSDCHEM\1\DATA\051402\DFTPPB.D
Acq On : 15 May 2002 12:55 am
Sample :
Misc :
MS Integration Params: EVENTS.E

Vial: 1
Operator: Mclean
Inst : Instrumen
Multiplr: 1.00

Method : C:\MSDCHEM\1\METHODS\050102.M (Chemstation Integrator)
Title : 508 Modified GC/MS scan

Spectrum Information: Scan 3601

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	48.0	120872	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	40.5	102032	PASS
70	69	0.00	2	0.6	662	PASS
127	198	40	60	50.3	126680	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	251776	PASS
199	198	5	9	7.0	17504	PASS
275	198	10	30	19.2	48240	PASS
365	198	1	99	2.1	5394	PASS
441	443	1	100	73.0	20896	PASS
442	198	40	110	55.9	140864	PASS
443	442	17	23	20.3	28632	PASS

DFTPPB.D 050102.M

Fri May 17 10:27:04 2002