

508 MOD.

Date Extracted 1-3-02

Page 7

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Ref add	IS add	Final vol.	T.V. pos.	pH	Count
L.B.		500 ml		1 ml			1 ml	5B	7	
Ref.					0.5 ml			5A		
31663	1 sec 2/11							4A		
7								4B		
31662	1 sec 2/11							2B		
6	1 sec 2/11							7B		
102								3B		
8								3A		
9								6A		
5								6B		
10	F.B.							1B		
107	F.B.							1A		
104								7A		
Remarks										
Reviewed by, <u>C/D 5/17/2</u>										
drip rate must be between 5-15 ml/min for liq-liq extractors										
drip rate -										
Extraction type - <u>Sep. F.</u> Solv. lot# <u>V39466</u> NaCl lot# <u>V4589</u> Na2SO4 lot# <u>V391100</u>										
STANDARD	INT STD LOG#									
CALIBRATION	<u>HPO1020</u>									
REF/SPIKE	<u>91759</u>									
MDL/MDRF										
Q. CHECK										
INT. STD.	<u>1591020</u>									
SURROGATE	<u>129101</u>									
LPC										
ENDRIN/DDT										

Analyst A.D.
revision 4/12/01

Verified by WJ

Reviewed by DN
Date 2/13/2

File Name: C:\HPCHEM\1\SEQUENCE\JAN2402.S
 Comment: INST 209; 30Mx0.25umX0.50
 Injector: 209 GALOS
 Path: C:\HPCHEM\1\data\jan2402\
 Cmd:
 Seq Cmd:


Method Sections To Run	On A Barcode Mismatch
(X) Full Method	(X) Inject Anyway
() Reprocessing Only	() Don't Inject

Line	Type	Vial	DataFile	Method	Sample Name
1	Sample	1	DFTPP	GCMS0103	
2	Sample	2	CAL25	GCMS0103	
3	Sample	16	31134	GCMS0103	REINJECT
4	Sample	17	31218	GCMS0103	REINJECT
5	Sample	3	LB	GCMS0103	GC/MS SCAN 1/3
6	Sample	31	SURR	GCMS0103	GC/MS SCAN 1/3
7	Sample	4	REF	GCMS0103	GC/MS SCAN 1/3
8	Sample	5	31663	GCMS0103	GC/MS SCAN 1/3
9	Sample	6	7	GCMS0103	GC/MS SCAN 1/3
10	Sample	7	31662	GCMS0103	GC/MS SCAN 1/3
11	Sample	8	6	GCMS0103	GC/MS SCAN 1/3
12	Sample	9	102	GCMS0103	GC/MS SCAN 1/3
13	Sample	10	8	GCMS0103	GC/MS SCAN 1/3
14	Sample	11	9	GCMS0103	GC/MS SCAN 1/3
15	Sample	12	5	GCMS0103	GC/MS SCAN 1/3
16	Sample	13	10	GCMS0103	GC/MS SCAN 1/3
17	Sample	14	107	GCMS0103	GC/MS SCAN 1/3
18	Sample	15	104	GCMS0103	GC/MS SCAN 1/3
19	Sample	18	LBA	GCMS0103	GC/MS SCAN 1/2
20	Sample	19	REFA	GCMS0103	GC/MS SCAN 1/2
21	Sample	20	31549	GCMS0103	GC/MS SCAN 1/2
22	Sample	21	31587	GCMS0103	GC/MS SCAN 1/2
23	Sample	22	31499	GCMS0103	GC/MS SCAN 1/2
24	Sample	23	31498	GCMS0103	GC/MS SCAN 1/2
25	Sample	24	31551	GCMS0103	GC/MS SCAN 1/2
26	Sample	25	31591	GCMS0103	GC/MS SCAN 1/2
27	Sample	26	31589	GCMS0103	GC/MS SCAN 1/2
28	Sample	27	31590	GCMS0103	GC/MS SCAN 1/2
29	Sample	28	31588	GCMS0103	GC/MS SCAN 1/2
30	Sample	29	31592	GCMS0103	GC/MS SCAN 1/2
31	Sample	30	31630	GCMS0103	GC/MS SCAN 1/2

DFTPP

: C:\HPCHEM\1\DATA\JAN2402\DFTPP.D
 : 24 Jan 2002 10:06 am
 :
 :
 MS Integration Params: GOOFY.P
 Method : C:\HPCHEM\1\METHODS\GCMS0103.M (RTE Integrator)
 Title : 209 625/TCLP calibration table

Vial: 1
 Operator: 209 GALOS
 Inst : GC/MS 209
 Multiplr: 1.00



Spectrum Information: Average of 25.997 to 26.043 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	45.4	111292	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	55.1	134905	PASS
70	69	0.00	2	0.1	176	PASS
127	198	40	60	49.5	121262	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	244956	PASS
199	198	5	9	6.7	16322	PASS
275	198	10	30	20.8	50959	PASS
365	198	1	100	2.4	5999	PASS
441	443	0.01	100	91.9	24587	PASS
442	198	40	110	58.2	142443	PASS
443	442	17	23	18.8	26768	PASS