

Method 508 MOD.

Date Extracted 1-10-02

Page 20

| Lab #                                                                                                        | Description  | Sample<br>vol/wt | Drip rate<br>ml/min | Surr<br>add | <del>Spl</del><br>add | IS<br>add | Final<br>vol. | T.V.<br>pos. | pH | Count |  |  |
|--------------------------------------------------------------------------------------------------------------|--------------|------------------|---------------------|-------------|-----------------------|-----------|---------------|--------------|----|-------|--|--|
| L.B.                                                                                                         |              | 500 ml           |                     | 1 ml        |                       |           | 1 ml          | 1B           | 7  |       |  |  |
| Ref.                                                                                                         |              |                  |                     |             | 0.5 ml                |           |               | 2A           |    |       |  |  |
| 7030.                                                                                                        |              |                  |                     |             |                       |           |               | 2B           |    |       |  |  |
| 700.                                                                                                         |              |                  |                     |             |                       |           |               | 3A           |    |       |  |  |
| 704.                                                                                                         |              |                  |                     |             |                       |           |               | 4A           |    |       |  |  |
| 701.                                                                                                         |              |                  |                     |             |                       |           |               | 3B           |    |       |  |  |
| 702.                                                                                                         |              |                  |                     |             |                       |           |               | 4B           |    |       |  |  |
| 697 FB F.B.                                                                                                  |              |                  |                     |             |                       |           |               | 1A           |    |       |  |  |
| 740 FB F.B.                                                                                                  |              |                  |                     |             |                       |           |               | 5A           |    |       |  |  |
| 727.                                                                                                         |              |                  |                     |             |                       |           |               | 6B           |    |       |  |  |
| 738.                                                                                                         |              |                  |                     |             |                       |           |               | 5B           |    |       |  |  |
| 739.                                                                                                         |              |                  |                     |             |                       |           |               | 6A           |    |       |  |  |
| 699                                                                                                          |              |                  |                     |             |                       |           |               | 7A           |    |       |  |  |
| Remarks                                                                                                      |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| BATCH #                                                                                                      |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| \$GCMS-16473                                                                                                 |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| drip rate must be between 5-15 ml/min for liq-liq extractors                                                 |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| drip rate -                                                                                                  |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| Extraction type - <u>SEP FUN.</u> Solv. lot# <u>V39466</u> NaCl lot# <u>V11599</u> Na2SO4 lot# <u>V39H00</u> |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| STANDARD                                                                                                     | INT STD LOG# |                  |                     |             |                       |           |               |              |    |       |  |  |
| CALIBRATION                                                                                                  |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| (REP)/SPIKE                                                                                                  | 91759        |                  |                     |             |                       |           |               |              |    |       |  |  |
| MDL/MDRF                                                                                                     |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| O. CHECK                                                                                                     |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| INT. STD.                                                                                                    |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| SURROGATE                                                                                                    | GC HS 010302 |                  |                     |             |                       |           |               |              |    |       |  |  |
| LPC                                                                                                          |              |                  |                     |             |                       |           |               |              |    |       |  |  |
| ENDRIN/DDT                                                                                                   |              |                  |                     |             |                       |           |               |              |    |       |  |  |


Analyst A.D.  
Vision 4/12/01

Verified by \_\_\_\_\_

Reviewed by DN  
Date 1/29/02

## DFTPP

Data File : C:\MSDCHEM\1\5973N\02JAN16\0116DFT.D Vial: 1  
Acq On : 16 Jan 2002 8:21 am Operator:  
Sample : dftpp Inst : Instrume  
Misc : January 16, 2002 Multiplr: 1.00  
MS Integration Params: SPLIT.P  
Method : C:\MSDCHEM\1\METHODS\METHODS\SCAN508.M (RTE Integrator)  
Title : 508 scan method



Spectrum Information: Average of 24.118 to 24.141 min.

| Target Mass | Rel. to Mass | Lower Limit% | Upper Limit% | Rel. Abn% | Raw Abn | Result Pass/Fail |
|-------------|--------------|--------------|--------------|-----------|---------|------------------|
| 51          | 198          | 30           | 60           | 43.5      | 235946  | PASS             |
| 68          | 69           | 0.00         | 2            | 0.5       | 1072    | PASS             |
| 69          | 198          | 0.00         | 100          | 39.8      | 216192  | PASS             |
| 70          | 69           | 0.00         | 2            | 0.5       | 1092    | PASS             |
| 127         | 198          | 40           | 60           | 52.4      | 284565  | PASS             |
| 197         | 198          | 0.00         | 1            | 0.5       | 2858    | PASS             |
| 198         | 198          | 100          | 100          | 100.0     | 543018  | PASS             |
| 199         | 198          | 5            | 9            | 6.7       | 36136   | PASS             |
| 275         | 198          | 10           | 30           | 22.4      | 121384  | PASS             |
| 365         | 198          | 1            | 99           | 2.5       | 13484   | PASS             |
| 441         | 443          | 1            | 99           | 74.5      | 74488   | PASS             |
| 442         | 198          | 40           | 110          | 94.7      | 514282  | PASS             |
| 443         | 442          | 17           | 23           | 19.4      | 99946   | PASS             |

0116DFT.D SCAN508.M

Wed Jan 16 10:14:25 2002