

Method 504.1 508 508A 515.1
Instrument 113

File Directory C:\NCHem\5 Data\ 508

Storage Directory M:\Instruments\113\Data\ 508

MS Chem
011226

011226

Lab #	Filename	Description	Vial #		Inj.	Standards Used - Lot #			Remarks	Initials
			front	rear		ul	Calibration	Reference		
	400	Chex	1							AN
	1226GND	quadrupole	2				2472			
	70x C1	cho cdu1	3			1926/2450	21		} cov	
	chl C1	che cdu1	4			2472/2450				
	Pent C1	pent cdu1	5			19680/2450				
	1127 R1 T	ref for	22				2109	2064		
	1128 R1 T	↓	23				↓			
	1127 R1 C	ref chl	24				2108			
	1128 R1 C	↓	25				↓			
	27828 SP	spike pent	26				2062			
	27828 SP	↓	27				↓			
	27848 SP	↓	28				↓			
	27848 SP	↓	29				↓			
	27817 ST	spike for	30				2109			
	27817 ST	↓	31				↓			
	28320 ST	↓	32				↓			
	28320 ST	↓	33				↓			
	27273 SC	spike chl	34				2108			
	27273 SC	↓	35				↓			
	28323 SC	↓	36				↓			
	28323 SC	↓	37				↓			
	27770	sample	38							
	27773	↓	39							
	27774	↓	40							
	27775	↓	41							
	27778	↓	42							
	27779	↓	43							
	27817	↓	44							

Date 12/26/01

Extraction date 11/27, 28/01; 12/26/01

Page 20

cont'd on next page

Instrument 113

File Directory C:\IPChem\4\Data\

Storage Directory M:\Instruments\113\Data\

Date 12/26/07

Extraction date 11/27, 28/09 ; 12/26/01

Page 21