

Method 308

Date Extracted 11/5/01

Page 75

Lab #	Description	Sample	Drip rate	Surr	Spk	IS	Final	T.V.	pH	Count		
		vol/wt	ml/min	add	add	add	vol.	pos.	=7			
	Lab Bldg 500ml			1ml		-	2ml					
	Ref TOXaphen											
	Ref Pest + Chlor											
	PEL 122											
	26203									1		
	26335									2		
	26336									3		
	26337									4		
	26386									5		
	26387									6		
	26388									7		
	26389									8		
	26390									9		
	26391									10		
26335	spike chlor + dup	Pest		1ml								
26336	spike TOX + dup	Chlor		↓								
26337	spike Pest + dup	TOX										
	26392									1		
	26393									2		
	26394									3		
	26433									4		
	26434									5		
	26470									6		
	26471									7		
	26472									8		
	26473									9		
	26474	✓								10		
Remarks												
26203 - neutral												
drip rate must be between 5-15 ml/min for liq-liq extractors												
Extraction type - Solv. lot# NaCl lot# Na2SO4 lot#												
STANDARD	INT STD LOG#											
CALIBRATION												
REF/SPIKE												
MDL/MDRF												
Q. CHECK												
INT. STD.												
SURROGATE												
LPC												
ENDRIN/DDT												

Reviewed by: B 3/12/02

Analyst CM/BD/ML
revision 4/12/01

Verified by _____

Reviewed by DN
Date 12/21/11

