



AIR MONITORING DATA AND CHAIN OF CUSTODY FORM

Page _____ of _____

TURNAROUND TIME

☐ 2 hrs ☐ 6 hrs ☐ 48 hrs ☐ Standard
☐ 4 hrs ☐ 24 hrs ☐ 72 hrs ☐ Other _____

PROJECT INFORMATION

1 Client		2 Project Name <u>QUALITY AIR</u>		3 Project #		4 Project Monitor <u>Keen</u>		4c Rotameter Number <u>SD087</u>	
		2a Project Address <u>100 CHAMBER ST NYP</u>		3a Task # <u>00001</u>		4a Air Sampler <u>ET PATRICK</u>			
5 Date <u>09-19-01</u>		6 Abatement Location <u>27TH FLOOR</u>		7 ☒ PCM (0.8 micron MCE) Cassette/Filter Manufacture <u>URGAL</u> Lot # _____		8 ☐ TEM (0.45 micron MCE) Cassette/Filter Manufacture _____ Lot # _____		9 Type a ☐ Background c ☐ Abatement b ☐ Pre-Abate d ☐ Post-Abate e ☐ OSHA f ☐ Environmental g ☒ Ambient h ☐ Other	

DAILY AIR SAMPLE RECORD

10 Sample I.D. Number	11 Pump Number	12 Sample Location		12c Sample Coordinates	13. Time (24 hrs clock)			14 Flow Rate (L/min)			15 Total Air Vol (liters)	16 # Fibers/Fields Minus Blank	17 Fiber Conc F/CC
		12a IWA	12b OWA		13a Start	13b End	13c Total	14a Start	14b End	14c Avg			
01				Field Blank								0/100	
02				Field Blank								0/100	
03				INSIDE APT R 27 TH FLOOR	9:00 pm	11:00 pm	120	15	15	15	1800	1/100	40.003
04				INSIDE APT D 27 TH FLOOR	9:02 pm	11:02 pm	120	15	15	15	1800	2/100	40.003
05				INSIDE CENTER OF 27 TH FLOOR	9:04 pm	11:04 pm	120	15	15	15	1800	1/100	40.003
06				SOUTHSIDE OF 27 TH FLOOR	9:07 pm	11:06 pm	120	15	15	15	1800	1/100	40.003
07				EASTSIDE OF 27 TH FLOOR	9:10 pm	11:10 pm	120	15	15	15	1800	1.5/100	40.003

CHAIN OF CUSTODY

18 Relinquished By		19 Received By		20 Date	21 Time
I <u>Kevin F. [Signature]</u>		II <u>A. Elbert</u>		09.20.01	7:00
III					
26 ATC Project Manager <u>TREV B</u>		27 Result To _____ Beeper/Phone #'s _____			

LAB INFORMATION (ATC ELAP # 10879, ATC NVLAP # 1187)

22. Lab Name		23 Date	24 Time	25 Microscope Make/Serial # <u>238913</u> ☐ Nikon Labophot 238913 ☐ Nikon Labophot 238917 ☐ Nikon Labophot 239079
a Analyzed By <u>A. Elbert</u>		09/20/01	8:00	
b QC By				
c Lab Batch # <u>97016</u>				
28 Comments				Analysis Methodology ☒ PCM - NIOSH 7400A ☐ TEM - AHERA

WHITE - File

YELLOW - Client/Report

PINK - Lab

GOLDENROD - Field

NY-1 (1/01)

CITY_21MC100_002848730